

Interpret

NGS Analysis Software



A Sysmex Group Company

OGT Handbook

Interpret User Guide




Table of contents

Introduction	5
Notices	5
Overview	5
Resources Required to Operate Interpret	5
Preparing Data For Analysis	6
Accessing Interpret	7
Login to Interpret	7
The Dashboard	7
Logging Out of Interpret	8
The Dashboard View	9
Loading FASTQ Files	10
Renaming a Sample ID	13
Pairing 'unpaired' FASTQ Files	14
Proceeding with File Upload	16
Remote Copy of FASTQ files	18
Viewing Samples	20
Adding User-defined Variables	25
Running an Analysis	29
Viewing Analysis Batches	39
Deleting Batches	41
Individual Batches	42
Batch QC	43
MultiQC	44
Sample QC	45
Batch Functions	46
Viewing Analysis QC	52
Sample QC	57



Viewing Analysis Results By Sample	59
Viewing SNV and Indel Events	62
SNV Options	66
Viewing CNV and LOH Events	79
Manual Creation of CNVs	90
Merging CNV calls	94
Separating Merged CNV calls	96
Aneuploidy Plots	97
Viewing Translocation Events	99
Viewing a Sample in IGV	106
Viewing Analysis Results By Variant	117
Administration Controls	128
Administration Overview	129
User Controls	129
Current User	129
Add User	131
Analysis	135
Manage Samples	135
Current Analyses	139
Protocols	141
Panels	155
Region Lists	158
Variant Lists	165
Classifications	168
Metric Sets	170
Manage Links	174
Filters	176
Preferred Transcripts	183
Reports	186
Guidelines	188
Software	190
Advanced Settings	190
Annotation	192



Plug-ins	192
Software Overview	193
Reporting	194
Product-specific Guidance	194
Measurable Residual Disease	194
Overview	194
Discovery Mode	194
Hotspot Monitoring	195
Monitoring Mode	198
Selecting Hotspots	198
Appendix	199
Attribute Definitions	199
Legal Information	207
Interpret Software	207
Customer's obligations	207



Introduction

User Guide for Software Version 3.8

Released April 2025



This guide is a manual for using OGT's Interpret NGS analysis software and is designed to be used in conjunction with OGT's range of NGS panel products.

Notices

Limitations of Use:

This software product is classified as **For Research Use Only**. It is not intended for use in diagnosis or treatment of human or animal diseases.

Symbols Used In This Guide:

Symbol	Meaning
	Attention - Denotes critical information that users need to be aware of.
	Error Icons - Such as the one shown below, are highlighted in red and to continue they must be removed. This is easily accomplished by simply clicking on the icon.

Overview

Interpret is a powerful bioinformatic tool designed to allow easy and comprehensive analysis of NGS data generated from OGT's NGS Panel products. The software functionality is limited solely to OGT's NGS panels – other panels CANNOT be loaded. The input is raw data in the form of FASTQ files.

Resources Required to Operate Interpret

Access to the software is via a web browser, OGT recommend Google Chrome, however the data is processed and stored by OGT's NGS analysis pipeline. OGT provides an installer to manage the process of deploying the software. Please consult the separate relevant installation guide for further information on this.

Users must ensure that browsers have pop-ups enabled.

The computational resources required for processing data with Interpret will depend on the number of samples being processed and the depth of sequencing. As a guideline we would recommend:

	Standard Panels		Myeloid MRD	Myeloid Fusion
Requirement	Minimum	Recommended	Minimum	Minimum
Memory (RAM)	16 GB	24 GB	64 GB	48 GB
CPU Cores	8	16	36	16
Storage (Disk space)	500 GB	2 TB	500 GB	500 GB
Operating System	Windows: Any version supporting Hyper-V Unix: Any version supporting Docker			

Table 1: Minimum and recommended hardware requirements for running Interpret

Preparing Data For Analysis

• Demultiplexed FASTQ files – REQUIRED

For each sample to be processed a pair of corresponding paired end FASTQ files are required. **The software expects the FASTQ files to be compressed with gzip.**

Software to demultiplex FASTQ files is not part of the functionality provided by Interpret and must be implemented prior to loading of the FASTQ files into Interpret. It is assumed that this will be provided by the sequencing instrument vendor.

• Target Regions File – REQUIRED

The Target Regions File is supplied by OGT with each of our NGS Panel products. The Target Regions File defines the regions covered by its associated panel. Interpret will use the regions within the panel file to define range of the analysis.

The Target Regions Files supplied by OGT are in a proprietary format; **ONLY** files supplied by OGT can be used by Interpret.

• Protocol File – REQUIRED

The Protocol File defines settings used by the NGS analysis pipeline. A default protocol is supplied with the panel but additional protocols can be created and stored within the software.

Accessing Interpret

Interpret is accessed through a web browser, all up to date browsers should work, however we recommend using Google Chrome.

To discover the correct web address of the software please consult the Installation Guide.

Login to Interpret

The login screen for Interpret will appear as shown below. To login submit the user name and password for your account.

If you do not have an account please contact your administrator to have one created.




Figure 1: The login window of Interpret

The Dashboard

After successful login the default dashboard page will be displayed as shown below:




Figure 2: The dashboard view of Interpret

Logging Out of Interpret

To logout of the software move the mouse to the user icon on the top right of the dashboard page and in the drop-down select 'Logout'




Figure 3: Accessing the logout option of Interpret

The Dashboard View

The dashboard view displayed below comprises 3 sections:

- Menu Bar
- Dashboard buttons to provide function shortcuts
- User account options




Figure 4: Annotated view of the dashboard

Menu Bar

The menu bar provides access to the functionality:

- **Home** - Link back to the Dashboard View
- **Batches** - Setting up and reviewing analysis batches
- **Samples** - Sample related functions
- **Variants** - Provides a means to view all data from a variant centric view
- **Help & Support** - A means to provide feedback as well request support
- **Tools** - Access to any additional tools




Figure 5: The menu bar from the dashboard

Dashboard Buttons



These provide shortcuts to the common actions required by users.

- **Upload FASTQs** – Select and upload FASTQ files.
- **Run Batch** – Run an analysis of a batch of loaded sample files
- **View Batches** – View the results of the batch analyses




Figure 6: Shortcut icons on the dashboard view

User Options

The User Options drop down menu gives the user access to a range of administration tools. Additionally this section of the dashboard displays the build of the genome being used as well as the version of the software. In this case it is GRCh37 and v3.3.61.

The drop down options are as follows:

- **My Account** – Your account details
- **Admin Controls** – Additional options described in detail in the Admin Options section of the guide
- **Logout** – Return to the Login page




Figure 7: User account options

Loading FASTQ Files

To load FASTQ files, on the dashboard either select 'Upload FASTQs' in the drop down from the 'Batches' menu item.




Figure 8: Accessing FASTQ uploads from the menu bar

Or, click on the 'Upload FASTQs' icon on the dashboard page




Figure 9: Accessing FASTQ uploads via the dashboard short-cut

Either choice opens the Upload FASTQs window.

Paired FASTQs		
Read 1	Read 2	Sample ID
Drop files here		

Figure 10: Upload FASTQs form

Initially the window shows an empty table with the heading "Paired FASTQs" with column headers for Read1, Read2 and Sample ID.

Paired FASTQs		
Read 1	Read 2	Sample ID

Figure 11: Initial view of the FASTQ upload table

To load FASTQ files users can either click on the 'Select FASTQ files' button to open a file browser or they can drag and drop files directly.




Figure 12: Methods of uploading files to Interpret

! The software requires all FASTQ files to be compressed by gzip. Any file without this file extension will not be loaded.

When FASTQ files are loaded the software will try to automatically pair them into Read 1 and Read 2. This is based on the filenames automatically generated by Illumina sequencers.

The file matching protocol assumes that a file pair shares the same name up to the annotation for whether the file is for Read1 data or Read2 data which for Illumina would be either `_R1` or `_R2`.

For example, with the following pair of FASTQ files:

Read1: Sample-400_R1.fastq.gz

Read2: Sample-400_R2.fastq.gz

The software would be automatically select the Sample ID as the portion of the file name highlighted in magenta.

Sample-400_R1.fastq.gz

Sample-400_R2.fastq.gz

As there is a Read1 file and Read2 file with the matching file names the software will automatically pair them. In the first instance the Sample ID is set to Sample-400, though, this can be easily changed.

These can be seen now populating the "Paired FASTQs" table.





Figure 13: Automatic pairing of FASTQ files

Renaming a Sample ID


As you can see in the previous screenshot a Sample ID is initially displayed with a red background as this is the ID that has been automatically generated.

If the Sample ID is not correct then it can be easily modified by clicking on the Sample ID that needs to be modified.


In the popup, the current Sample ID will be shown.



In order to change this, enter the required name in the text box and press enter.



The updated name is then displayed.




Selecting Done updates the display with the new Sample ID
In the example below the background colour has changed to green to represent the Sample ID for that pair of files has been modified by the user.

<div>Select FASTQ Files</div>			
Paired FASTQs		Sample ID	
Read 1			
Sample-400_R1.fastq.gz	0%	Sample-400_R2.fastq.gz	0%
		Updated Sample Name	Different_Name-400

Figure 14: Table with updated Sample ID

Pairing 'unpaired' FASTQ Files

Sometimes it may be the case that files have been paired incorrectly.
If this is the case and there are files that are not a pair then clicking on the bin icon  will remove the samples from the Paired FASTQs table and move them to a new table called "Unpaired FASTQs".

Upload FASTQs

Select FASTQ Files


Unpaired FASTQs

File Name	Size (bytes)	
Sample-400_R1.fastq.gz	89,503,947	<div></div>
Sample-400_R2.fastq.gz	94,390,832	<div></div>

Paired FASTQs

Read 1	Read 2	Sample ID

Figure 15: Table with unmatched FASTQ files

From here they can easily combine as a pair by highlighting one file and then clicking on the second. When this is done, they are considered to be a pair again and moved back to the Paired FASTQs table.
Alternatively clicking on the bin icon  for files in the Unpaired FASTQs table removes them from the Upload FASTQs page.

When files without matching file names are loaded, they are initially displayed in a panel for Unpaired FASTQs.

However, an alternative pair of read files without matching names such as:

This_Fastq_R1.fastq.gz

and

That_Fastq_R2.fastq.gz

would not be automatically paired. This needs to be completed by the user.




Figure 16: Table with unmatched FASTQ files highlighted

Selecting two of these files by shift clicking automatically denotes them to be a pair and they are moved to the Paired FASTQ panel.

The software will select the Read1 file in the pair to be the first FASTQ file that is selected, so in the example below ONE_Fastq_R1.fastq.gz is selected so it becomes the Read1 file.

However, initial selection the other file ANOTHER_Fastq_R2.fastq.gz would have resulted in the alternative situation




Figure 17: Selecting an unpaired FASTQ file

The software creates a name for the Sample ID based on the file name for the file denoted as the Read1 file, but this can be changed by following the protocol for changing a Sample ID.




Figure 18: Table showing the pairing of unmatched FASTQ files and the automatically assigned name

Proceeding with File Upload

Now that the files have been paired and correct Sample IDs have been assigned, they are ready to be uploaded.




Figure 19: Manually paired FASTQ files ready for upload

At the bottom of the Paired FASTQs table are 3 parameters that allow tracking of the upload.

These are:

- An estimate of the amount of time remaining.
- The amount of data uploaded as an amount in Mb out of the total as well as a percentage completed statistic.
- The proportion of files uploaded.




Figure 20: Three different metrics for upload of FASTQ files can be monitored

Selection of the "Upload Paired FASTQ Files" button starts the upload process, though there is an additional check run by the software when there is a Sample ID that has been automatically assigned.

In this case the user must confirm the Sample ID is correct before the upload is initiated.




Figure 21: Confirmation of automatic Sample ID assignment

Once it has started users are able to track progress of the upload by looking at the blue progress bars for each of the files. In order for the download to progress it is important to keep the browser tab open (or in a separate window) while uploading FASTQs, otherwise the upload will pause.




Figure 22: Progress bars tracking file uploading

Alternatively there are the tracking metrics displayed which can give a more precise estimation of how the upload is progressing.




Figure 23: Metrics for tracking file uploading

Once the upload is complete, the software asks whether the user would like to create a new analysis run or to wait. Setting up an analysis will be covered in the next section.




Figure 24: Popup menu displayed once FASTQ file upload is complete

Remote Copy of FASTQ files

In addition to previous method of loading FASTQ files it is also possible to run a remote copy of files accessible on the network mounted to Interpret.

This is accessed by Batches drop-down in the menu bar.




Figure 25: Accessing the 'Remote Copy FASTQs' function from the menu bar




Figure 26: The remote FASTQ file copy interface

If input FASTQ files are located in a folder which is accessible by the web application, it is possible to navigate the file system and select FASTQ files for upload to the system directly from this folder. Upload using this method has the advantage of avoiding the need to maintain an active web browser during the upload process, as files are copied in the background. Whether such folders are accessible to the web application will depend on many factors, including how the software has been installed. It may be the case that it is not possible to access input data folders, but if the required folders are not accessible, contact OGT for assistance.

To identify FASTQ files, expand the file tree to the required location. FASTQ files will be displayed in green, and can be added to the list by double-clicking. Alternatively, double-click on a folder to add all FASTQ files in that folder. In some browsers, double-clicking on a folder will instead set the "root" location of the file tree to the selected folder, and it may be necessary to double-click a second time in order to add the files to the list.

Files can be removed from the list of Paired FASTQ files by selecting the associated checkbox and clicking the red delete button. Once the required set of FASTQ files have been selected, click the **Copy Files** button to begin the upload process. As with the **Upload FASTQ Files** section, users then have option of creating a batch from the selected FASTQ files.



Sample IDs

Please note that the **Remote Copy** interface does not provide a means to modify the sample ID associated with each FASTQ file - these are instead automatically determined by the name of the FASTQ file based on the standard nomenclature.



Viewing Samples

Once the FASTQ files have been uploaded the samples will be available to view in the Samples page

Accessing this is via the Samples button on the dashboard menu bar shown in the figure below.



Figure 27: Selection of Samples from the Dashboard menu bar

Samples and status are displayed on the left hand side of the window and when a sample is selected further information is displayed on the right hand side.

When a sample is first loaded there is no additional information present.





Figure 28: The initial view of samples in the Sample page

Samples can be searched by entering a part of the sample name in the search box. In the example below all samples containing "10" are displayed.

 **Samples**

⏮

◀

Page 1 of 1 (1 - 3 of 3)

▶

⏭

Page Size: 25 ▾

Sample	Status
<input type="text" value="10"/>	<input type="text"/>
10384	To Do
8210	To Do
10847	To Do

Figure 29: Searching for samples containing 10

The status of a sample can be updated. When first loaded, the status will be set to "To Do" and will be updated to "Running Pipeline" once processing has begun. Once sample processing is complete, the status will be further updated to "In Review". Users can assess the results if the analysis and manually update the sample status to "Completed" as required.

Sample

10384

Status

To Do

▼

Mother

To Do

Running Pipeline

In Review

Completed

Father

Runs

▼

Figure 30: Modifying the status of a sample in the Sample view

It is also possible to specify the mother and father of a sample if they are also loaded in Interpret.

Sample 10384

Status To Do ▼

Mother ▼

Father ▼

Runs ▼

- 10847
- 11516
- 12878
- 14130
- 4315
- 5881
- 6937
- 7408
- 8210

Figure 31: Specifying the mother of a sample

Initially, before any analysis, the run drop-down list will be empty.

Sample 10384

Status To Do ▼

Mother ▼

Father ▼

Runs ▼

Figure 32: A sample that has not been processed yet having no run data listed

When a sample has been analysed, each run can be accessed from the drop-down menu. Each run will have a set of data which is displayed by 3 tabs. These are for general run information, QC metrics and results of the analysis.

The General tab displays basic information about the analysis and provides a link to batch view.

Sample 10384

Status In Review ▼

Mother ▼

Father ▼

Runs 14 Jan 2021 12:14:27 ▼

General QC Metrics Results

Date	14 Jan 2021 12:14:27
Batch	CytoSure NGS Batch 00001
Status	Completed
Panel	CytoSure NGS Panel
Protocol	Default Protocol
ISCN	seq[GRCh37] (1-22,X)x2
FASTQ read1	10384_7_L001_R1_001.fastq.gz
FASTQ read2	10384_7_L001_R2_001.fastq.gz

Figure 33: Viewing the General tab for a sample run

The QC Metrics tab gives an overview of the metrics of the sample. The data will be colour-coded according to the metric set that was defined for the analysis protocol.

There is further information on metric sets in the section of the guide that covers the admin options.




Figure 34: Viewing the QC Metrics tab for a sample run

Finally, the Results tab provides links (in green) to download files from the analysis as well to view (in blue) the different variants that have been detected.




Figure 35: Viewing the Results tab for a sample run

Adding User-defined Variables

In order to enable the user to capture and report custom information related to samples processed in Interpret, the admin controls section provides a means to create variables of different data types via Admin Controls > Analysis > Manage Samples > Variables.




Figure 36: The manage samples page in the Admin Controls

Selecting **Add New Variable Category** provides a text box to name the new variable and clicking **Add** adds a new sub-tab to the **Variables** tab.




Figure 37: Creating a new variable category named Additional

With the new category called "Additional" generated, users can create variables associated with the category by selecting **Add New Variable**.




Figure 38: An empty custom category named "Additional"

To delete an empty category, click the **Remove Empty Category** button. To create new variables in the category, click the **Add New Variable** button, assign a **Name** to the variable, confirm the **Category** with which it should be associated, and select the appropriate data type from the **DataType** drop-down box, then click the **Add** button.





Figure 39: Creating a new variable named "Phenotype" in a category named "Additional".



The image shows a dialog box titled "Add New Variable...". It has three main fields: "Name", "Category", and "DataType". The "Name" field contains the text "Phenotype". The "Category" field is a dropdown menu currently showing "Additional". The "DataType" field is also a dropdown menu, which is currently open, displaying a list of options. The options in the "DataType" dropdown are: "Integer", "Floating Point Number", "Percentage", "Text - Single Line", "Text - Multiple Lines" (which is highlighted with a blue background), "True or False", "Date", "Date and Time", and "List - Select One".

Figure 40: Selecting the appropriate data type for the new variable - in this case, "Text - Multiple Lines"

Once a variable has been created, it will be listed, along with its data type, in the appropriate category in the **Manage Sample Variables** section. To delete a variable, click on the cross next to the variable name.

Manage Samples

Overview [Variables](#)

Manage Sample Variables

[Additional](#)

[+ Add New Variable](#)

Phenotype (Text - Multiple Lines) [×](#)

[+ Add New Variable Category](#)

[Save Changes](#) [Discard Changes](#) [← return](#)

Figure 41: A custom field named "Phenotype" listed under the "Additional" category

Having been created in the system, custom fields may be populated for each sample in the **Samples** view, and will also be displayed in the sample run page whenever the sample has been processed in a batch (accessible by clicking on the sample row in the **Completed Samples** table in the **Batch Overview** page).

Samples

Page 1 of 1 (1 - 20 of 20) Page Size: 25

Sample	Status
CR007-012	In Review
CR007-006	In Review
CR007-001	In Review
CR007-008	In Review
CR007-011	In Review
CR007-009	In Review
CR007-007	In Review
CR007-010	In Review
CR007-002	In Review
CR007-005	In Review
Sample-4	In Review
Sample-3	In Review
Sample-2	In Review
Sample-1	In Review

Sample-1

Status:

Mother:

Father:

Runs:

General QC Metrics Results [Additional](#)

Phenotype: The phenotype has been defined for this patient:
 HP:0000488 Retinopathy Noninflammatory retina disease
 HP:0000556 Retinal dystrophy

Reviewed: ☒ Yes ☐ No

[Save](#)

Figure 42: Editing the content of the "Phenotype" field in the Samples page

Interpret also provides a framework enabling the development of plug-ins to import sample data in bulk from other sources, such as spreadsheets, text files or a LIMS. If you are interested in importing data in bulk, contact OGT - a suitable plug-in may be available, or it may be possible to develop a plug-in to satisfy your requirements.

Running an Analysis

On the dashboard either select "Run Batch" in the drop down from the 'Batches' menu item.




Figure 43: Selection of Run Batch from the Dashboard menu bar Batches drop down menu

Or, click on the 'Run Batch' icon on the dashboard page





Figure 44: Selection of Run Batch from the dashboard short-cut buttons

Either choice leads to the initial Run Batch page is as follows:



Run Batch

Batch Name * 2020-01-07 17:44:08 Panel * Protocol *

Search for Samples by Sample ID

ID	Sample ID	FASTQ Read 1	FASTQ Read 2
1	5881	5881_7_L001_R1_001.fastq.gz	5881_7_L001_R2_001.fastq.gz
2	6937	6937_7_L001_R1_001.fastq.gz	6937_7_L001_R2_001.fastq.gz
3	7408	7408_7_L001_R1_001.fastq.gz	7408_7_L001_R2_001.fastq.gz
4	8210	8210_7_L001_R1_001.fastq.gz	8210_7_L001_R2_001.fastq.gz
5	10384	10384_7_L001_R1_001.fastq.gz	10384_7_L001_R2_001.fastq.gz
6	10847	10847_7_L001_R1_001.fastq.gz	10847_7_L001_R2_001.fastq.gz
7	11516	11516_7_L001_R1_001.fastq.gz	11516_7_L001_R2_001.fastq.gz

Selected Samples (0)

ID	Sample ID	FASTQ Read 1	FASTQ Read 2
----	-----------	--------------	--------------


→ Run Analysis

Figure 45: Initial view of the Run Batch window

Besides showing the list of available samples there are additional text fields and drop-down menus.

In order to run an analysis the user needs to

1. Select samples for the analysis
2. Select the correct panel for the samples
3. Select the analysis protocol



Run Batch

Batch Name * 2020-01-06 15:33:39 Panel * Protocol *

Search for Samples by Sample ID

Annotations: Set Batch Name (points to Batch Name field), Select Panel (points to Panel dropdown), Select Protocol (points to Protocol dropdown), Search for Samples (points to search bar).

Figure 46: Input fields for the Run Batch window

Optionally users can specify a name for the batch analysis. A default batch name is provided with the date followed by the time in the format YYYY-MM-DD HH:MM:SS.

In the example below the user has created the batch name CytoSure NGS Batch 00001



Run Batch

Batch Name * CytoSure NGS Batch 00001 Panel * Protocol *

Search for Samples by Sample ID

Figure 47: Entering a batch name

The samples have been processed with OGT's CytoSure NGS panel so that is the selection to make from the Panel dropdown menu.

Run Batch

Batch Name * CytoSure NGS Batch 00001

Panel * CytoSure NGS Panel

Protocol * Default Protocol

Search for Samples by Sample ID

Figure 48: Selecting the appropriate Panel

The user now specifies the protocol that will be used, in this case the Default Protocol.

Run Batch

Batch Name * CytoSure NGS Batch 00001

Panel * CytoSure NGS Panel

Protocol * Default Protocol

Search for Samples by Sample ID

Figure 49: Selecting the protocol to use for processing the batch



Panel-Protocol Compatibility

Only protocols whose **Pipeline Type** are included in the list of pipeline types supported by the selected **Panel** will be listed in the **Protocol** drop-down list. Additionally, if any **Pipeline Capabilities** supported by the protocol are not supported by the selected panel, a warning will be displayed indicating which processes will not be run.

Lastly, the user specifies the samples to be analysed.

There may be a large number of samples loaded into the system, so to enable easier sample selection it is possible to add a search term. In this case the user is looking for all samples containing the number 5. Additionally, search terms are independent of the case used.

Run Batch

Batch Name * CytoSure NGS Batch 00001

Panel * CytoSure NGS Panel

Protocol * Default Protocol

Search for Samples by Sample ID

ID	Sample ID	FASTQ Read 1	FASTQ Read 2
1	5881	5881_7_L001_R1_001.fastq.gz	5881_7_L001_R2_001.fastq.gz
7	11516	11516_7_L001_R1_001.fastq.gz	11516_7_L001_R2_001.fastq.gz
10	4315	4315_7_L001_R1_001.fastq.gz	4315_7_L001_R2_001.fastq.gz

Figure 50: Filtering loaded samples with a search term

Selecting the checkbox next to a sample moves a loaded sample into the Selected Samples table.

Run Batch

Batch Name * CytoSure NGS Batch 00001 Panel * CytoSure NGS Panel Protocol * Default Protocol

Search for Samples by Sample ID

ID	Sample ID	FASTQ Read 1	FASTQ Read 2
1	5881	5881_7_L001_R1_001.fastq.gz	5881_7_L001_R2_001.fastq.gz
2	6937	6937_7_L001_R1_001.fastq.gz	6937_7_L001_R2_001.fastq.gz
3	7408	7408_7_L001_R1_001.fastq.gz	7408_7_L001_R2_001.fastq.gz
4	8210	8210_7_L001_R1_001.fastq.gz	8210_7_L001_R2_001.fastq.gz
5	10384	10384_7_L001_R1_001.fastq.gz	10384_7_L001_R2_001.fastq.gz
6	10847	10847_7_L001_R1_001.fastq.gz	10847_7_L001_R2_001.fastq.gz
7	11516	11516_7_L001_R1_001.fastq.gz	11516_7_L001_R2_001.fastq.gz

Selected Samples (4) *

ID	Sample ID	FASTQ Read 1	FASTQ Read 2
1	5881	5881_7_L001_R1_001.fastq.gz	5881_7_L001_R2_001.fastq.gz
2	6937	6937_7_L001_R1_001.fastq.gz	6937_7_L001_R2_001.fastq.gz
3	7408	7408_7_L001_R1_001.fastq.gz	7408_7_L001_R2_001.fastq.gz
4	8210	8210_7_L001_R1_001.fastq.gz	8210_7_L001_R2_001.fastq.gz

→ Run Analysis

Figure 51: Adding a sample to an analysis batch

Clicking on the minus icon  will remove the sample from the Selected Samples tables.

Run Batch

Batch Name * CytoSure NGS Batch 00001 Panel * CytoSure NGS Panel Protocol * Default Protocol

Search for Samples by Sample ID

ID	Sample ID	FASTQ Read 1	FASTQ Read 2
1	5881	5881_7_L001_R1_001.fastq.gz	5881_7_L001_R2_001.fastq.gz
2	6937	6937_7_L001_R1_001.fastq.gz	6937_7_L001_R2_001.fastq.gz
3	7408	7408_7_L001_R1_001.fastq.gz	7408_7_L001_R2_001.fastq.gz
4	8210	8210_7_L001_R1_001.fastq.gz	8210_7_L001_R2_001.fastq.gz
5	10384	10384_7_L001_R1_001.fastq.gz	10384_7_L001_R2_001.fastq.gz
6	10847	10847_7_L001_R1_001.fastq.gz	10847_7_L001_R2_001.fastq.gz
7	11516	11516_7_L001_R1_001.fastq.gz	11516_7_L001_R2_001.fastq.gz

Selected Samples (10) *

ID	Sample ID	FASTQ Read 1	FASTQ Read 2
1	5881	5881_7_L001_R1_001.fastq.gz	5881_7_L001_R2_001.fastq.gz
2	6937	6937_7_L001_R1_001.fastq.gz	6937_7_L001_R2_001.fastq.gz
3	7408	7408_7_L001_R1_001.fastq.gz	7408_7_L001_R2_001.fastq.gz
4	8210	8210_7_L001_R1_001.fastq.gz	8210_7_L001_R2_001.fastq.gz
5	10384	10384_7_L001_R1_001.fastq.gz	10384_7_L001_R2_001.fastq.gz
6	10847	10847_7_L001_R1_001.fastq.gz	10847_7_L001_R2_001.fastq.gz
7	11516	11516_7_L001_R1_001.fastq.gz	11516_7_L001_R2_001.fastq.gz
8	12878	12878_7_L001_R1_001.fastq.gz	12878_7_L001_R2_001.fastq.gz

Remove Sample from Analysis

→ Run Analysis

Figure 52: Removing a sample from an analysis batch

When all selections have been made the run can be started by selecting 

Run Batch

Batch Name * CytoSure NGS Batch 00001 Panel * CytoSure NGS Panel Protocol * Default Protocol

Search for Samples by Sample ID

ID	Sample ID	FASTQ Read 1	FASTQ Read 2
1	5881	5881_7_L001_R1_001.fastq.gz	5881_7_L001_R2_001.fastq.gz
2	6937	6937_7_L001_R1_001.fastq.gz	6937_7_L001_R2_001.fastq.gz
3	7408	7408_7_L001_R1_001.fastq.gz	7408_7_L001_R2_001.fastq.gz
4	8210	8210_7_L001_R1_001.fastq.gz	8210_7_L001_R2_001.fastq.gz
5	10384	10384_7_L001_R1_001.fastq.gz	10384_7_L001_R2_001.fastq.gz
6	10847	10847_7_L001_R1_001.fastq.gz	10847_7_L001_R2_001.fastq.gz
7	11516	11516_7_L001_R1_001.fastq.gz	11516_7_L001_R2_001.fastq.gz

Selected Samples (10) *

ID	Sample ID	FASTQ Read 1	FASTQ Read 2
1	5881	5881_7_L001_R1_001.fastq.gz	5881_7_L001_R2_001.fastq.gz
2	6937	6937_7_L001_R1_001.fastq.gz	6937_7_L001_R2_001.fastq.gz
3	7408	7408_7_L001_R1_001.fastq.gz	7408_7_L001_R2_001.fastq.gz
4	8210	8210_7_L001_R1_001.fastq.gz	8210_7_L001_R2_001.fastq.gz
5	10384	10384_7_L001_R1_001.fastq.gz	10384_7_L001_R2_001.fastq.gz
6	10847	10847_7_L001_R1_001.fastq.gz	10847_7_L001_R2_001.fastq.gz
7	11516	11516_7_L001_R1_001.fastq.gz	11516_7_L001_R2_001.fastq.gz
8	12878	12878_7_L001_R1_001.fastq.gz	12878_7_L001_R2_001.fastq.gz

[Run Analysis](#)

Figure 53: Starting an analysis

! If the selected protocol has "Enable CNV and LOH Calling" set to "Yes", CNVs will be detected by comparison with a set of reference samples which need to be defined in the protocol as either "All Batch Samples" or a specific set of reference samples whose FASTQ files have already been uploaded and designated to the system by the user.

In the latter case, OGT may provide a set of data files that can be used as a reference set for CNV analysis. As more samples are processed users may extend the reference pool by adding any samples they believe are suitable as controls for CNV calling. A user can modify samples designated as reference pool in the protocol in Admin Controls-Manage Samples-Protocols.

If CNV calling is enabled without a reference data set being defined, then, on selecting [Run Analysis](#), the following error will be displayed.




Figure 54: Error message for using an invalid protocol

Click on the message to remove the warning and select **Admin Controls > Analysis > Protocols** to set reference samples. More details are in the Protocols section of this User Guide.

Otherwise, a popup presents the chosen files and selected parameters. Following this there is a request for confirmation and upon confirmation the analysis run will be initiated.

Are you sure you would like to start the following Analysis?

Name CytoSure NGS Batch 00001

SureSeq Panel CytoSure NGS Panel

Protocol Default Protocol

Samples (10)

ID	Sample ID	FASTQ Read 1	FASTQ Read 2
1	5881	5881_7_L001_R1_001.fastq.gz	5881_7_L001_R2_001.fastq.gz
2	6937	6937_7_L001_R1_001.fastq.gz	6937_7_L001_R2_001.fastq.gz
3	7408	7408_7_L001_R1_001.fastq.gz	7408_7_L001_R2_001.fastq.gz
4	8210	8210_7_L001_R1_001.fastq.gz	8210_7_L001_R2_001.fastq.gz
5	10384	10384_7_L001_R1_001.fastq.gz	10384_7_L001_R2_001.fastq.gz
6	10847	10847_7_L001_R1_001.fastq.gz	10847_7_L001_R2_001.fastq.gz
7	11516	11516_7_L001_R1_001.fastq.gz	11516_7_L001_R2_001.fastq.gz
8	12878	12878_7_L001_R1_001.fastq.gz	12878_7_L001_R2_001.fastq.gz
9	14130	14130_7_L001_R1_001.fastq.gz	14130_7_L001_R2_001.fastq.gz
10	4315	4315_7_L001_R1_001.fastq.gz	4315_7_L001_R2_001.fastq.gz

Cancel **OK**

Figure 55: Window requesting confirmation to run an analysis

Selecting **OK** will start the analysis and the display will change to show information about batch being analysed.

Within this there is an overview window providing an overview of the analysis and a sample window giving information about the status of each sample.




Figure 56: The batch processing view with the overview window and sample window highlighted

Initially the status of the samples will be listed in the overview window as "Waiting to start"




Figure 57: Initial batch status before analysis starts

Waiting for a reference to be generated

If a reference pool needs to be generated the status shown in the batch overview will report this and provide a means to track progress of the reference pool creation.



Figure 58: The analysis status in the overview window, highlighted, showing reporting the status of pre-processing of the reference samples

The status of reference building can also be tracked in the View Batches window which is discussed in the View Batches section of this User Guide.




Figure 59: Reference building status being shown in the View Batches window

If the protocol performs CNV analysis and samples in the analysis are to be used to generate the reference pool against which to make CNV calls then the overview will report the combining of the reference samples.

Once the reference samples have been aligned and counted, they are combined into a pool for the CNV analysis




Figure 60: The analysis in the over window, highlighted, reporting the combining of the reference samples into a pool

Samples will be queued until there is capacity available in the pipeline. Once this is available the software will start processing the samples sequentially. The stage of the process is updated and the overall progress can be monitored in the progress bar.




Figure 61: The batch view showing progress of analysis

Once analysis started the stage of each sample is displayed and can be followed




Figure 62: Tracking progress of a sample processing

Once a sample has been analysed the overview updates the count and a summary of the analysis is displayed in a Completed Samples table.




Figure 63: The first completed sample is displayed below the samples to be processed

When all samples are completed




Figure 64: An analysis with all samples analysed

There is no need to wait until all samples have been processed to view the results for a completed sample.

This will be discussed in the Viewing Analysis Results section of the manual.

Viewing Analysis Batches

On the dashboard either select "View Batches" in the drop down from the 'Batches' menu item.




Figure 65: Selecting View Batches from the menu bar drop down menu

Or, click on the 'View Batches' icon on the dashboard page





Figure 66: Selecting View Batches from the dashboard shortcut buttons

The Batches are presented in a table as below.




The Batches window displays a table with the following data:

Name	Panel	Protocol	Number of Samples	Status	Date	User
CytoSure NGS Batch 0002	CytoSure NGS Panel	Default Protocol	3	Waiting to start	10 Jan 2020 11:38:27	admin
CytoSure NGS Batch 0001	CytoSure NGS Panel	Default Protocol	10	Completed	8 Jan 2020 14:51:53	admin

Figure 67: Initial view of the Batches window


As with other tables in Interpret where there is a column selector icon  a user can add or remove columns from the display



The Batches window shows the column selector icon circled in red. A red box labeled "Column Selector" points to the icon.

Figure 68: Column selection options for the Batches window

Column names annotated with a tick are in the current display and changes can easily be made to add or remove columns



The "Select Displayed Columns" dialog shows the following columns with checkmarks indicating they are currently displayed:

- ✓ Name
- ✓ Panel
- ✓ Protocol
- ✓ Number of Samples
- ✓ Status
- ✓ Date
- ✓ User
- Samples

Figure 69: Selection of columns to display in the Batches window

By default, all batches are presented in the first instance but these can easily be filtered.

Where the column header has a text field, users can type in a search term and all batches with that text contained somewhere in the name, will be retained. The text search is independent of lower- or upper-case letters, "Demo" will return the same samples as "demo".

Alternatively, where there is a drop-down menu selecting one of the values in the menu will lead only to the batches matching the selection being displayed, for example, below only batches that have completed will be displayed.




Figure 70: Filtering batches on status

Lastly, there are date fields, allowing selection of batches run within a set time frame.




Figure 71: Filtering batches on date of processing

Deleting Batches

In the batch view it is possible to delete batches. When first opened there is a greyed out Delete button in the display.

If a batch is selected it is highlighted in blue and Delete button is now active, Clicking the delete button will delete the batch from the software.




Figure 72: Selection of a batch to delete highlights the delete button

If the delete button is selected there will be a popup box requesting confirmation of the deletion.

Selecting  will lead to the batch being deleted.



Once a batch is deleted it **CANNOT** be recovered.




Figure 73: Popup box requesting confirmation of batch deletion

Individual Batches

Clicking on a row in the View Batches page will open a new page showing the selected batch in more detail.

There are 3 parts to the information provided,

1. Overview

The overview provides information about the analysis

2. Batch Functions

The batch functions allow users to download files from the analysis, repeat the analysis or generate a report

3. Sample Details

In this part there is the headline information from the run about each sample such as the number of SNVs and CNVs called.

Batch - CytoSure NGS Batch 00001

Overview **Batch Functions** **Sample results**

▼ Overview

Operator admin Date 14 Jan 2021 12:14:27

Panel CytoSure NGS Panel Status Completed

Protocol [Default Protocol](#) MultiQC Report [1 MultiQC](#)

File Status Copy Batch Stop Batch Generate Report

For selected, view: Completed Samples

Sample	View	# SNVs	# CNVs	# LOH	# Translocations	Report	QC
10384	SNVs CNVs/LOH Translocations Logs	2,774	2	31	0	View	View
10847	SNVs CNVs/LOH Translocations Logs	2,786	2	33	0	View	View
11516	SNVs CNVs/LOH Translocations Logs	2,685	5	38	0	View	View
12878	SNVs CNVs/LOH Translocations Logs	2,731	6	34	0	View	View
14130	SNVs CNVs/LOH Translocations Logs	2,734	4	31	0	View	View
4315	SNVs CNVs/LOH Translocations Logs	3,514	4	16	0	View	View
5881	SNVs CNVs/LOH Translocations Logs	2,883	5	33	0	View	View
6937	SNVs CNVs/LOH Translocations Logs	2,815	8	34	0	View	View
7408	SNVs CNVs/LOH Translocations Logs	2,846	6	25	0	View	View
8210	SNVs CNVs/LOH Translocations Logs	2,780	5	39	0	View	View

Figure 74: The sections of the batch analysis window

Batch QC

Included in the Batch page are two QC reports.

In the batch overview there is a link to a MultiQC report which gives an overview of all the samples that were in the batch.

Additionally, each sample in the completed table has a FastQC report for each read file.

Examples of both of these QC reports are shown below.

Batch - CytoSure NGS Batch 00001

Overview

Operator admin Date 14 Jan 2021 12:14:27

Panel CytoSure NGS Panel Status Completed

Protocol [Default Protocol](#) MultiQC Report [1 MultiQC](#)

File Status Copy Batch Stop Batch Generate Report

For selected, view:

Completed Samples

Sample	View	# SNVs	# CNVs	# LOH	# Translocations	Report	QC
10384	SNVs CNVs/LOH Translocations Logs	2,774	2	31	0	View	View
10847	SNVs CNVs/LOH Translocations Logs	2,786	2	33	0	View	View
11516	SNVs CNVs/LOH Translocations Logs	2,685	5	38	0	View	View
12878	SNVs CNVs/LOH Translocations Logs	2,731	6	34	0	View	View
14130	SNVs CNVs/LOH Translocations Logs	2,734	4	31	0	View	View
4315	SNVs CNVs/LOH Translocations Logs	3,514	4	16	0	View	View
5881	SNVs CNVs/LOH Translocations Logs	2,883	5	33	0	View	View
6937	SNVs CNVs/LOH Translocations Logs	2,815	8	34	0	View	View
7408	SNVs CNVs/LOH Translocations Logs	2,846	6	25	0	View	View
8210	SNVs CNVs/LOH Translocations Logs	2,780	5	39	0	View	View

Figure 75: Links to QC reports for a batch and a sample

MultiQC

MultiQC is a reporting tool for the whole batch of samples. It parses summary statistics from results and log files generated by other bioinformatics tools.

When you launch MultiQC, it recursively searches through any provided file paths for specific files. These files are parsed for relevant information and used to generate a single stand-alone HTML report file. It also saves a directory of data files with all parsed data for further use downstream. To save MultiQC report to user's computer, right click on the page, and choose "Save as..."

Additional information about MultiQC can be found in the next section of this guide.




Figure 76: Example of a MultiQC report

Sample QC

FastP is used for sample QC data generation. Clicking on the  button on the sample view will open up a new tab in the web browser with the sample QC details.

1

Summary

General

fastp version:	0.20.0 (https://github.com/OpenGene/fastp)
sequencing:	paired end (151 cycles + 151 cycles)
mean length before filtering:	147bp, 147bp
mean length after filtering:	147bp, 147bp
duplication rate:	8.103519%
Insert size peak:	173

Before filtering

total reads:	32.386588 M
total bases:	4.776467 G
Q20 bases:	4.399816 G (92.114459%)
Q30 bases:	4.130451 G (86.475024%)
GC content:	47.035643%

After filtering

total reads:	30.184820 M
total bases:	4.440014 G
Q20 bases:	4.185280 G (94.262757%)
Q30 bases:	3.945594 G (88.864444%)
GC content:	46.795308%

Filtering result

reads passed filters:	30.184820 M (93.201606%)
reads with low quality:	2.094360 M (6.466751%)
reads with too many N:	52.744000 K (0.162858%)
reads too short:	54.664000 K (0.168786%)

Figure 77: Example of a FastP report

Batch Functions

Below the overview section there are a set of buttons providing a set of option - when the batch has finished processing the Stop Batch button is disabled.




Figure 78: Batch options




File Status

File status provides shows the files that have been generated for each sample during the analysis.

Files provided are:

1. Alignment files
2. QC files
3. VCF files
4. CGH files for loading into CytoSure Interpret
5. Log files

Where a green tick is displayed, that file is available for download and this can be achieved by clicking on the  button.

File Status of Batch: CytoSure NGS Batch 0001

← Return to Batch

Bulk Download






























































Id	Sample	BAM	BAI	QC	VCF	CGH	Log	Actions
21	5881	✓ 	✓ 	✓ 	✓ 	✓ 	✓ 	
22	6937	✓ 	✓ 	✓ 	✓ 	✓ 	✓ 	
23	7408	✓ 	✓ 	✓ 	✓ 	✓ 	✓ 	
24	8210	✓ 	✓ 	✓ 	✓ 	✓ 	✓ 	
25	10384	✓ 	✓ 	✓ 	✓ 	✓ 	✓ 	
26	10847	✓ 	✓ 	✓ 	✓ 	✓ 	✓ 	
27	11516	✓ 	✓ 	✓ 	✓ 	✓ 	✓ 	
28	12878	✓ 	✓ 	✓ 	✓ 	✓ 	✓ 	
29	14130	✓ 	✓ 	✓ 	✓ 	✓ 	✓ 	
30	4315	✓ 	✓ 	✓ 	✓ 	✓ 	✓ 	

Figure 79: Status of files generated by the pipeline for each sample

It is possible to download all files, or selected files, simultaneously via the bulk download button 

Id	<input type="checkbox"/> Sample	BAM <input type="checkbox"/>	BAI <input type="checkbox"/>	QC <input type="checkbox"/>	VCF <input type="checkbox"/>	CGH <input type="checkbox"/>	Log <input type="checkbox"/>
21	<input type="checkbox"/> 5881	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
22	<input type="checkbox"/> 6937	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
23	<input type="checkbox"/> 7408	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24	<input type="checkbox"/> 8210	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
25	<input type="checkbox"/> 10384	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
26	<input type="checkbox"/> 10847	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
27	<input type="checkbox"/> 11516	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
28	<input type="checkbox"/> 12878	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
29	<input type="checkbox"/> 14130	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
30	<input type="checkbox"/> 4315	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Download Selected (0) Grouped By Sample ▾

Figure 80: Bulk download file selector

Specific files can be selected as below

Id	<input type="checkbox"/> Sample	BAM <input type="checkbox"/>	BAI <input type="checkbox"/>	QC <input type="checkbox"/>	VCF <input type="checkbox"/>	CGH <input type="checkbox"/>	Log <input type="checkbox"/>
21	<input checked="" type="checkbox"/> 5881	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
22	<input type="checkbox"/> 6937	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
23	<input type="checkbox"/> 7408	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24	<input type="checkbox"/> 8210	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
25	<input type="checkbox"/> 10384	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
26	<input type="checkbox"/> 10847	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
27	<input type="checkbox"/> 11516	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
28	<input type="checkbox"/> 12878	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
29	<input type="checkbox"/> 14130	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
30	<input type="checkbox"/> 4315	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Download Selected (6) Grouped By Sample ▾

Figure 81: Bulk download with single sample selected

Alternatively, all files can be selected for download




Figure 82: Bulk download with all files selected

Files downloaded in bulk can be grouped by sample or file type




Figure 83: Bulk download selecting grouping by Sample of downloaded files

Copy Batch

This function allows the user to repeat the batch analysis with the same settings. When selected a Run Batch window opens and if the user selects to Run Analysis the processing will be repeated. The software will automatically update the Batch Name but otherwise nothing is changed including the time stamp.

Run Batch

Batch Name * CytoSure NGS Batch 0001(2) SureSeq Panel * CytoSure NGS Panel Protocol * Default Protocol

Search for Samples by Sample ID

ID	Sample ID	FASTQ Read 1	FASTQ Read 2
1	5881	5881_7_L001_R1_001.fastq.gz	5881_7_L001_R2_001.fastq.gz
2	6937	6937_7_L001_R1_001.fastq.gz	6937_7_L001_R2_001.fastq.gz
3	7408	7408_7_L001_R1_001.fastq.gz	7408_7_L001_R2_001.fastq.gz
4	8210	8210_7_L001_R1_001.fastq.gz	8210_7_L001_R2_001.fastq.gz
5	10384	10384_7_L001_R1_001.fastq.gz	10384_7_L001_R2_001.fastq.gz
6	10847	10847_7_L001_R1_001.fastq.gz	10847_7_L001_R2_001.fastq.gz
7	11516	11516_7_L001_R1_001.fastq.gz	11516_7_L001_R2_001.fastq.gz

Selected Samples (10) *

ID	Sample ID	FASTQ Read 1	FASTQ Read 2
1	5881	5881_7_L001_R1_001.fastq.gz	5881_7_L001_R2_001.fastq.gz
2	6937	6937_7_L001_R1_001.fastq.gz	6937_7_L001_R2_001.fastq.gz
3	7408	7408_7_L001_R1_001.fastq.gz	7408_7_L001_R2_001.fastq.gz
4	8210	8210_7_L001_R1_001.fastq.gz	8210_7_L001_R2_001.fastq.gz
5	10384	10384_7_L001_R1_001.fastq.gz	10384_7_L001_R2_001.fastq.gz
6	10847	10847_7_L001_R1_001.fastq.gz	10847_7_L001_R2_001.fastq.gz
7	11516	11516_7_L001_R1_001.fastq.gz	11516_7_L001_R2_001.fastq.gz
8	12878	12878_7_L001_R1_001.fastq.gz	12878_7_L001_R2_001.fastq.gz

Updated Batch Name

Run Analysis

Figure 84: A batch analysis being repeated using the Copy Batch option showing the updated batch name

Report Generation

Report Generation shows a drop down in which the user can select the report to be generated

Report Type

Generate Report

Figure 85: Initial view of the report options

Currently, the only template loaded is the Batch Report

Report Type

Batch Report with template

Figure 86: Selecting a report type for the QC of the run




Figure 87: Selection of a template for the report

When the report is generated the output is a table with a set of metrics for each sample in the batch.

Sam ple	Percent Reads Aligned	Percent Duplicat ion	Mean Target Coverage	Targets Not Covered	Aligned Reads GC	Aligned Reads Per Base Quality	Usable On Target Reads	Usable On Target Bases
5881	99.4	43.3	536	0	40	37.4	53.9219	35.5563
693 7	99.4	43.1	627	0	40	37.4	55.4457	36.6627
740 8	99.1	50.2	440	0	40	37.3	39.6108	25.9758
821 0	99.1	50.7	418	0	41	37.3	40.9108	26.9532

Table 2: Example output of the QC report for a batch

Selecting a completed sample or samples allows viewing of the variant information and this is described in Viewing Analysis Results section.

Batch - CytoSure NGS Batch 00001

▼ Overview

Operator admin Date 14 Jan 2021 12:14:27

Panel CytoSure NGS Panel Status Completed

Protocol [Default Protocol](#) MultiQC Report [1 MultiQC](#)

File Status
 Copy Batch
 Stop Batch
 Generate Report

For selected, view:

- SNVs/Indels
- Translocations
- CNVs/LOH Calls

Completed Samples

Sample	View	# SNVs	# CNVs	# LOH	# Translocations	Report	QC
10384	SNVs CNVs/LOH Translocations Logs	2,774	2	31	0		View
10847	SNVs CNVs/LOH Translocations Logs	2,786	2	33	0		View
11516	SNVs CNVs/LOH Translocations Logs	2,685	5	38	0		View
12878	SNVs CNVs/LOH Translocations Logs	2,731	6	34	0		View
14130	SNVs CNVs/LOH Translocations Logs	2,734	4	31	0		View
4315	SNVs CNVs/LOH Translocations Logs	3,514	4	16	0		View
5881	SNVs CNVs/LOH Translocations Logs	2,883	5	33	0		View
6937	SNVs CNVs/LOH Translocations Logs	2,815	8	34	0		View
7408	SNVs CNVs/LOH Translocations Logs	2,846	6	25	0		View
8210	SNVs CNVs/LOH Translocations Logs	2,780	5	39	0		View

Figure 88: Selecting a batch to view

Viewing Analysis QC

When a batch of samples is processed, besides individual sample metrics that were discussed in the previous section, there is a batch QC report generated. This uses MultiQC and fastp to collate a set of metrics for each sample and merge into a set of graphs and tables.

The report can be accessed from in the batch overview displayed once a batch has completed analysis.

Batch - CytoSure NGS Batch 00001

▼ Overview

Operator admin Date 14 Jan 2021 12:14:27

Panel CytoSure NGS Panel Status Completed

Protocol [Default Protocol](#) MultiQC Report [1 MultiQC](#)

File Status
 Copy Batch
 Stop Batch
 Generate Report

MultiQC Report

Figure 89: Accessing the Batch QC report

When the user clicks on the MultiQC Report link a new tab opens up in the browser displaying the QC report. The view is divided into 3 parts – the quality control report for the batch, which comprises the bulk of the display, and 2 tabs that come into the

view from the left and the right of the page. These tabs can be viewed and hidden by clicking on their respective buttons. The second tab provides the MultiQC toolbox for:

1. The Quality Control report for the batch
2. The report short cut tab
3. The tool box tab

At the head is the quality control report; this provides general information about the analysis such as the date of the analysis and which user performed it.



Figure 90: Example batch overview details

OGT's Summary Table

Each sample has a row in the table with some key metrics.

Sample Name	Mapped reads	Duplicated reads	Avg. quality	Avg. insert size	Insert size std.	MPQ=0	On-target	±250bp	Off-target	Mean cov.	Not covered	Evenness score	Fold.80	Sex
10_bam	99.83%	8.93%	33.10	206.20	66.30	2.71%	73.08%	82.55%	17.45%	314.04	136	87.56%	1.346	Female
1_bam	99.57%	8.51%	33.10	200.80	61.50	2.99%	72.81%	81.56%	16.44%	315.89	133	87.56%	1.344	Female
2_bam	99.62%	9.27%	33.00	216.80	66.60	2.75%	72.06%	82.59%	17.41%	408.41	126	87.31%	1.357	Female
3_bam	99.60%	8.45%	33.00	213.30	65.30	2.61%	72.91%	83.09%	16.91%	390.11	121	87.43%	1.350	Female
4_bam	99.60%	8.40%	32.60	216.00	66.90	2.67%	72.51%	82.91%	17.09%	397.52	125	87.34%	1.357	Female
5_bam	99.61%	8.05%	33.10	203.10	61.90	2.89%	73.23%	82.24%	17.76%	349.09	132	87.68%	1.336	Female
6_bam	99.55%	9.07%	33.10	210.50	70.30	2.53%	71.96%	82.62%	17.18%	345.20	129	87.76%	1.338	Female
7_bam	99.63%	10.72%	33.10	206.20	64.20	2.64%	72.19%	81.59%	16.41%	558.37	117	87.87%	1.329	Female
8_bam	99.63%	8.58%	33.10	210.80	66.10	2.70%	72.40%	82.21%	17.79%	382.89	128	87.26%	1.357	Female
9_bam	99.58%	8.46%	33.10	204.40	65.80	2.67%	73.48%	82.81%	17.19%	357.73	125	87.60%	1.345	Female

Figure 91: Example sample QC summary table

The column names from the summary table are listed in the table below with some additional detail as to their meaning.

Column Name	Description
Mapped reads	The percentage of reads that mapped to the reference genome
Duplicate reads	The percentage of reads that were duplicated
Avg. quality	The average read quality reported by Samtools stats
Avg. insert size	The average insert size reported by Samtools stats
Insert size std	The standard deviation of the insert size reported by Samtools stats
MPQ = 0	The percentage of reads that were mapped that have a mapping quality of 0
On-target	The percentage of reads that map on target that are not duplicate reads
± 250bp	The percentage of reads that overlap target regions extended by 250bp



Column Name	Description
Off-target	The percentage of reads that are neither on target nor within the specified flanking region
Mean cov.	The mean target coverage
Not covered	The number of targets with a coverage of less than 1
Evenness score	The fraction of the whole sequencing output that is correctly distributed
Fold-80	The fold of additional sequencing that would be required to ensure that 80% of targeted bases achieve the mean target coverage.
Sex	The chromosomal sex of the sample predicted from the distribution of reads that map to the sex chromosomes

Table 3: Column names and their description from the QC summary table

Targets Not Covered

Any targets not covered are detailed, providing that they are not within a segmental duplication.



Figure 92: Example targets not covered summary

Coverage Efficiency

The efficiency of coverage as a measure of depth are displayed.




Figure 93: Example QC report summary

Insert Sizes Samtools

The distribution of insert sizes for each sample is displayed.




Figure 94: Example QC report summary

Percent Mapped

The percentage of base calls at each position for which an N was called.




Figure 95: Example QC report summary

Alignment Metrics

The alignment metrics for all the samples in the batch are plotted.




Figure 96: Example of the alignment metrics

Filtered Reads

The filtered reads graph shows the number or percentage of reads that have been removed by the filter.




Figure 97: Example QC report summary

Duplication Rates

The relative level of duplication found for each sample as a percentage.




Figure 98: Example QC report summary

Sequence Quality

The mean sequence quality or Phred score of each base in the insert for each sample.




Figure 99: Example QC report summary

GC Content




Figure 100: Example QC report summary

N Content




Figure 101: Example QC report summary

Sample QC

Sample QC information can also be access via the Sample page. For a particular run selecting the QC metrics tab will provide the relevant information. Colours are defined by the metric set used in the analysis protocol.




Figure 102: Viewing the metrics for an individual sample

The screenshot shows the 'Samples' page in the OpenGene interface. On the left, a table lists samples with their IDs and status. Sample 10384 is selected. On the right, there are tabs for 'General', 'QC Metrics', and 'Results'. The 'Results' tab is highlighted with a red box and labeled 'Results Tab'. Below the 'Results' tab, there are links for 'View SNVs', 'View CNVs/LOH Calls', and 'View Translocations'. Each link has a corresponding 'Download' button. A red box labeled 'Sample QC Report' points to the 'Download' button for the 'QC Translocations' link.

Figure 103: Accessing the Sample QC report from the sample results tab

10

Summary

General

fastp version:	0.20.0 (https://github.com/OpenGene/fastp)
sequencing:	paired end (151 cycles + 151 cycles)
mean length before filtering:	146bp, 146bp
mean length after filtering:	145bp, 145bp
duplication rate:	7.220084%
Insert size peak:	171

Before filtering

total reads:	28.905316 M
total bases:	4.227335 G
Q20 bases:	3.894346 G (92.122952%)
Q30 bases:	3.658105 G (86.534532%)
GC content:	47.412412%

After filtering

total reads:	26.944780 M
total bases:	3.928181 G
Q20 bases:	3.703276 G (94.274565%)
Q30 bases:	3.492954 G (88.920399%)
GC content:	47.151145%

Filtering result

reads passed filters:	26.944780 M (93.217386%)
reads with low quality:	1.862392 M (6.443078%)
reads with too many N:	46.288000 K (0.160137%)
reads too short:	51.856000 K (0.179400%)

Figure 104: Start of a FastP report for an individual sample



Viewing Analysis Results By Sample

Viewing a Sample

Access to the results from running the pipeline are described in the previous section "View Analysis Batches".

Within each batch are the samples processed in that batch comprising analysed variants and QC metrics.

For selected, view:


SNVs/Indels

Translocations

CNVs/LOH Calls

Completed Samples					# SNVs	# CNVs	# LOH	# Translocations	Report	QC	% Reads Aligned	
<input type="checkbox"/>	View											
<input type="checkbox"/>	<div><div>SNVs</div><div>CNVs/LOH</div><div>Translocations</div></div>	<div><div>Logs</div></div>	7	2	0	1	<div><div></div></div>	<div><div>View</div></div>	99.52			
<input type="checkbox"/>	<div><div>SNVs</div><div>CNVs/LOH</div><div>Translocations</div></div>	<div><div>Logs</div></div>	5	2	0	1	<div><div></div></div>	<div><div>View</div></div>	99.47			
<input type="checkbox"/>	<div><div>SNVs</div><div>CNVs/LOH</div><div>Translocations</div></div>	<div><div>Logs</div></div>	5	2	0	1	<div><div></div></div>	<div><div>View</div></div>	99.51			
<input type="checkbox"/>	<div><div>SNVs</div><div>CNVs/LOH</div><div>Translocations</div></div>	<div><div>Logs</div></div>	7	1	0	1	<div><div></div></div>	<div><div>View</div></div>	99.49			
<input type="checkbox"/>	<div><div>SNVs</div><div>CNVs/LOH</div><div>Translocations</div></div>	<div><div>Logs</div></div>	7	1	0	1	<div><div></div></div>	<div><div>View</div></div>	99.6			
<input type="checkbox"/>	<div><div>SNVs</div><div>CNVs/LOH</div><div>Translocations</div></div>	<div><div>Logs</div></div>	4	2	0	1	<div><div></div></div>	<div><div>View</div></div>	99.31			
<input type="checkbox"/>	<div><div>SNVs</div><div>CNVs/LOH</div><div>Translocations</div></div>	<div><div>Logs</div></div>	4	3	0	1	<div><div></div></div>	<div><div>View</div></div>	99.51			

Figure 105: View of a set of processed samples in the batch view

As with other tables in Interpret, where there is a column selection icon  users can use it to configure which columns are being displayed.

Completed Samples

<input type="checkbox"/>	Sample	View	# SNVs	# CNVs	# LOH	Report	QC
<input type="checkbox"/>	5881	<div>SNVSCNVs/LOHVCF</div> <div>Logs</div>	2,754	8	16	<div></div>	<div>View</div>

Column Selector




Figure 106: Column selector button for configuring columns to view in display

The column options for this view are shown in the Figure.



Select Displayed Columns

Sample
SNVs
CNVs
LOH
Translocations
% Reads Aligned
Average Quality
% Duplication
Mean Target Coverage
Targets Not Covered
% Usable On Target Reads
% Usable On And Near Target Reads
Off Target Reads
% Reads Mapping Quality 0
Average Insert Size
Insert Size Std
Eveness
Uniformity
Sample Sex
Exon Targets Not Covered
Seg Dup Exon Targets Not Covered
I SCN

Figure 107: Columns available for selection to display

There is one completed sample per row and for each sample there is a range of information available to view.

SNVs

CNVs/LOH

Translocations

Logs

Report

QC

For selected, view:	Completed Samples										
SNVs/Indels	View										
Translocations	SNVs	CNVs/LOH	Translocations	Logs	# SNVs	# CNVs	# LOH	# Translocations	Report	QC	% Reads Aligned
CNVs/LOH Calls	SNVs	CNVs/LOH	Translocations	Logs	7	2	0	1	View	View	99.52
	SNVs	CNVs/LOH	Translocations	Logs	5	2	0	1	View	View	99.47
	SNVs	CNVs/LOH	Translocations	Logs	5	2	0	1	View	View	99.51

Figure 108: Information available for each sample

Variants for a sample can be viewed by selecting the SNVs or CNVs/LOH links present in each row.

SNVs Translocations CNVs/LOH

Figure 109: Selecting multiple samples to view simultaneously

Variants

[← Return to Batch](#) [CNVs / LOH Calls](#)

▼ Protocol Filter: Default SNV Filter ⚙ Hide

2946 >
 Alt Strand Bias < 80% +
 Reads Placed Left > 4 +
 Reads Placed Right > 4 +
 Quality Score > 1.0 >
 2754

⏪ ⏩ Page 1 of 276 (1 - 10 of 2754) Page Size: 10 ▼ ⚙ Actions ▼

<input type="checkbox"/>	HGVSc (Gene Symbol)	Chromosome ▼	Start ▼	End ▼	Ref ▼	Alt ▼	Allele Frequency ▼	Type ▼	Ref Depth ▼	Alt Depth ▼	Total Depth ▼	Quality Score ▼	Ref Quality ▼	Alt Quality ▼	Sample ▼
<input type="checkbox"/>	DVL1:c.*347T>C	1	1335795	1335795	A	G	100%	SNV	0	138	138	4122.59	0.00	4771.00	5881
<input type="checkbox"/>	DVL1:c.366A>G	1	1342153	1342153	T	C	99.57%	SNV	1	232	233	6750.47	32.00	7624.00	5881
<input type="checkbox"/>	ATAD3A:c.671G>A	1	1520153	1520153	G	A	48.9%	SNV	185	177	362	4259.04	6194.00	6050.00	5881
<input type="checkbox"/>	ATAD3A:c.*154C>T	1	1534226	1534226	C	T	43.24%	SNV	21	16	37	275.18	738.00	515.00	5881
<input type="checkbox"/>	ATAD3A:c.*330C>C	1	1534402	1534402	G	C	51.07%	SNV	183	191	374	4665.67	6338.00	6565.00	5881
<input type="checkbox"/>	SKI:c.*1802T>C	1	2308567	2308567	T	C	100%	SNV	0	470	470	14005.10	0.00	15743.00	5881
<input type="checkbox"/>	SKI:c.*3182dup	1	2309937	2309937	C	CT	86.16%	Insertion	22	137	166	6276.51	752.00	8035.00	5881
<input type="checkbox"/>	CAMTA1:c.288C>T	1	6887707	6887707	C	T	40.99%	SNV	190	132	322	3142.37	6523.00	4644.00	5881
<input type="checkbox"/>	CAMTA1:c.*419T>C	1	6888144	6888144	T	C	47.55%	SNV	203	184	387	4498.84	7148.00	6375.00	5881
<input type="checkbox"/>	CAMTA1:c.1350G>A	1	7663897	7663897	G	A	51.5%	SNV	194	206	400	5107.83	6530.00	7083.00	5881

IGV hg38 chr6 chr6:1-170,805,979 170 mb

Cursor Guide Center Line Track Labels

mb 10 mb 20 mb 30 mb 40 mb 50 mb 60 mb 70 mb 80 mb 90 mb 100 mb 110 mb 120 mb 130 mb 140 mb 150 mb 160 mb 170 mb

At the top of the page are buttons that allow the user to toggle between the CNV and SNV views.




Figure 111: Toggling between SNV, CNV and Translocation reports

The Variants viewer is divided into three parts:

1. The Protocol Filter, initially this will be showing the Protocol Filter used in the analysis. The filter is modifiable in the Admin Controls (Admin Controls > Analysis > Protocols)
2. The Variant Table, showing the variants, one on each row. In the header there is a drop down "Actions" menu options; these are discussed below.
3. The Integrated Genome Viewer (IGV) that has been embedded in the software. Further details on using IGV are below.

Figure 112: The sections of a sample results view page

Viewing SNV and Indel Events

The variant table has a column selector icon  allowing user to configure which columns are displayed.

There are different columns available depending on whether you are viewing the SNV variants page or the CNV/LOH variant page. The options for SNVs are displayed below.




Figure 113: Columns available to display in the SNVs variant page

Selection of a variant will load the alignment file in IGV allowing review of the alignment.

A range of variants can be displayed and examples of each of these are:

1. SNV




Figure 114: Example of a SNV being displayed in the IGV browser

2. Deletion




Figure 115: Example of a deletion being displayed in the IGV browser

3. Insertion




Figure 116: Example of an insertion being displayed in the IGV browser

4. Complex




Figure 117: Example of a complex event being displayed in the IGV browser

5. Multi Nucleotide Polymorphism (MNP)




Figure 118: Example of a MNP variant being displayed in the IGV browser

6. Partial Tandem Duplication (PTD)




Figure 119: Example of a PTD being displayed in the IGV browser; the duplication event is highlighted by the transparent red box

SNV Options

Right clicking on a row will generate a popup menu with a range of options.




Figure 120: Options available for each SNV or Indel variant

Add to a Shortlist




Figure 121: Adding a variant to shortlist

Once a variant has been added to the shortlist it will be annotated with a tick.





Figure 122: Annotation of a selected variant


A variant can be removed from the shortlist using the Remove from Shortlist command.



Sample	HGVSc (Gene Symbol)	Chromosome	Start	End	Ref	Alt	Allele Frequency	Type	Ref Depth	Alt Depth	Total Depth	Ref Quality	Alt Quality	Quality Score
5881	ABCC9:c.*728_*729del	12	21923083	21923085	GAA	G	100%	Deletion	0	218	257	0.00	12940.00	11022.50
5881	ABCC9:c.2158G>A	12	21872665	21872665	C	T	53.48%	SNV	147	169	316	4930.00	5870.00	4241.90
5881	ABCC9:c.*1200T>C	12	21924253	21924253	A	G	50.84%	SNV	117	121	238	4131.00	4208.00	3009.09
5881	ABCC9:c.*787A>C	12	21923026	21923026	T	G			0	309	309	0.00	10312.00	9218.83
5881	ABCC9:c.*142C>G	12	21925311	21925311	G	C			0	276	276	0.00	9474.00	8450.52
5881	ABCD1:c.*894G>A	X	153744629	153744629	G	A			184	176	360	6225.00	6198.00	4430.94
5881	ABCD1:c.1548G>A	X	153744343	153744343	G	A			144	105	249	4906.00	3543.00	2359.41
5881	ABCD1:c.*59C>T	X	153740151	153740151	G	A			177	207	384	6092.00	7301.00	5364.42
5881	ABCD1:c.*877C>T	X	153725208	153725208	C	T			206	166	372	7170.00	5719.00	3971.18
5881	ABCD1:c.*877C>T	X	153744612	153744612	C	T			0	385	385	0.00	13436.00	12030.80
5881	ABCD1:c.*80G>C	X	153743743	153743743	G	C			1	342	343	36.00	11668.00	10131.30
5881	ABCD1:c.*194G>A	X	153729764	153729764	G	A			256	278	534	8879.00	9750.00	7106.52
5881	ABL1:c.-15_-14insC	9	130835432	130835432	T	TC	100%	Insertion	0	192	194	0.00	11486.00	9599.84

Figure 123: Removing a variant from a shortlist

Subsequently the tick annotation will be removed.




Sample	HGVSc (Gene Symbol)	Chromosome	Start	End	Ref	Alt	Allele Frequency	Type	Ref Depth	Alt Depth	Total Depth	Ref Quality	Alt Quality	Quality Score
5881	ABCC9:c.*728_*729del	12	21923083	21923085	GAA	G	100%	Deletion	0	218	257	0.00	12940.00	11022.50
5881	ABCC9:c.2158G>A	12	21872665	21872665	C	T	53.48%	SNV	147	169	316	4930.00	5870.00	4241.90
5881	ABCC9:c.*1200T>C	12	21924253	21924253	A	G	50.84%	SNV	117	121	238	4131.00	4208.00	3009.09

Figure 124: Annotation denoting shortlisting has been removed

Add to New List

Variants can be added to lists that can be used in software; for example a list of variants can be used in a filter as a means to specifically search for a data set.



Sample	HGVSc (Gene Symbol)	Chromosome	Start	End	Ref	Alt	Allele Frequency	Type	Ref Depth	Alt Depth	Total Depth	Ref Quality	Alt Quality	Quality Score
5881	ABCC9:c.*728_*729del	12	21923083	21923085	GAA	G	100%	Deletion	0	218	257	0.00	12940.00	11022.50
5881	ABCC9:c.2158G>A	12	21872665	21872665	C	T	53.48%	SNV	147	169	316	4930.00	5870.00	4241.90
5881	ABCC9:c.*1200T>C	12	21924253	21924253	A	G	50.84%	SNV	117	121	238	4131.00	4208.00	3009.09
5881	ABCC9:c.*787A>C	12	21923026	21923026	T	G			0	309	309	0.00	10312.00	9218.83
5881	ABCC9:c.*142C>G	12	21925311	21925311	G	C			0	276	276	0.00	9474.00	8450.52
5881	ABCD1:c.*894G>A	X	153744629	153744629	G	A			184	176	360	6225.00	6198.00	4430.94
5881	ABCD1:c.1548G>A	X	153744343	153744343	G	A			144	105	249	4906.00	3543.00	2359.41
5881	ABCD1:c.*59C>T	X	153740151	153740151	G	A			177	207	384	6092.00	7301.00	5364.42
5881	ABCD1:c.*877C>T	X	153725208	153725208	C	T			206	166	372	7170.00	5719.00	3971.18
5881	ABCD1:c.*877C>T	X	153744612	153744612	C	T			0	385	385	0.00	13436.00	12030.80
5881	ABCD1:c.*80G>C	X	153743743	153743743	G	C			1	342	343	36.00	11668.00	10131.30
5881	ABCD1:c.*194G>A	X	153729764	153729764	G	A			256	278	534	8879.00	9750.00	7106.52
5881	ABL1:c.-15_-14insC	9	130835432	130835432	T	TC	100%	Insertion	0	192	194	0.00	11486.00	9599.84

Figure 125: Adding a variant to New List

Initially users will be prompted to create a new variant list by setting the name of the list.




Figure 126: Creating a new variant list

In the example below a list called New List has been created.




Figure 127: Setting the name of a new variant list

The list New List is now available and variants can be added to it.

	Sample	HGVSc (Gene Symbol)	Chromosome	Start	End	Ref	Alt	Allele Frequency	Type	Ref Depth	Alt Depth	Total Depth	Ref Quality	Alt Quality	Quality Score
	5881	ABCC9:c.*728_*729del	12	21923083	21923085	GAA	G	100%	Deletion	0	218	257	0.00	12940.00	11022.50
	5881	ABCC9:c.2158G>A	12	21872665	21872665	C	T	53.48%	SNV	147	169	316	4930.00	5870.00	4241.90
	5881	ABCC9:c.*1200T>C	12	21924253	21924253	-	-	50.84%	SNV	117	121	238	4131.00	4208.00	3099.09
	5881	ABCC9:c.*787A>C	12	21923026	21923026	-	-	100%	SNV	0	309	309	0.00	10312.00	9218.83
	5881	ABCC9:c.*142C>G	12	21925311	21925311	-	-	100%	SNV	0	276	276	0.00	9474.00	8450.52
	5881	ABCD1:c.*608G>A	X	153744343	153744343	-	-	100%	SNV	144	105	249	4906.00	3543.00	2359.41
	5881	ABCD1:c.1548G>A	X	153740151	153740151	-	-	100%	SNV	177	207	384	6092.00	7301.00	5364.42
	5881	ABCD1:c.*59C>T	X	153725208	153725208	-	-	100%	SNV	206	166	372	7170.00	5719.00	3971.18
	5881	ABCD1:c.*877C>T	X	153744612	153744612	-	-	100%	SNV	0	385	385	0.00	13436.00	12030.80
	5881	ABCD1:c.*87G>C	X	153743743	153743743	-	-	99.71%	SNV	1	342	343	36.00	11668.00	10131.30
	5881	ABCD1:c.*194G>A	X	153729764	153729764	-	-	52.06%	SNV	256	278	534	8879.00	9750.00	7106.52
	5881	ABCD1:c.*894G>A	X	153744629	153744629	-	-	48.89%	SNV	184	176	360	6225.00	6198.00	4430.94
	5881	ABL1:c.*1756C>T	9	130887439	130887439	C	T	44.55%	SNV	178	143	321	6208.00	5058.00	3520.45

Figure 128: Adding a variant to the newly created list

Select Transcript

By default the a gene will have the largest canonical transcript set as the preferred transcript.

To change this, users can use the Set Preferred Transcript option to select an alternative transcript from the list available in the database.

	Sample ▼	▲ HGVSc (Gene Symbol)	Chromosome ▼	Start ▼	End ▼	Ref ▼	Alt ▼	Allele Frequency ▼	Type ▼	Ref Depth ▼	Alt Depth ▼	Total Depth ▼	Ref Quality ▼	Alt Quality ▼	Quality Score ▼
	5881	ARHGAP31:c.*1156C>G	3	119417420	119417420	C	G	47.31%	SNV	196	176	372	6676.00	6082.00	4278.64
✓	5881	ARHGAP31:c.3615C>T	3	119415544	119415544	C	T		ARHGAP31:c.3615C>T Add to Shortlist Add to...> Set Preferred Transcript Classify Variant Links... Notes	189	200	389	6474.00	7004.00	5069.57
	5881	ARHGAP31:c.*3968C>T	3	119420232	119420232	C	T			0	366	366	0.00	12771.00	11444.40
	5881	ARHGAP31:c.435C>G	3	119382295	119382295	C	G			157	183	340	5440.00	6257.00	4453.61
	5881	ARHGAP31:c.*1391_135	3	119417648	119417652	GCAC	G			211	192	406	6972.00	11450.00	8522.90
	5881	ARHGAP31:c.*194G>A	3	119294711	119294711	G	A			223	240	463	7445.00	8365.00	6059.01
	5881	ARHGAP31:c.*1176T>G	3	119417440	119417440	T	G								
	5881	ARHGAP31:c.*985_986G	3	119417249	119417249	A	AC								
	5881	ARHGAP31:c.2407G>A	3	119414336	119414336	G	A								
	5881	ARHGAP31:c.*2281C>G	3	119418545	119418545	C	G								
	5881	ARHGAP31:c.-38C>A	3	119294867	119294867	C	A								
	5881	ARHGAP31:c.*1238T>C	3	119417502	119417502	T	C	51.75%	SNV	152	163	315	5360.00	5568.00	3995.19
✓	5881	ARID1A:c.3869C>T	1	26773582	26773582	C	T	48.15%	SNV	252	234	486	8573.00	8071.00	5726.71

Figure 129: Setting a preferred transcript for a gene

Variant Classification

Users are able to classify variants in two ways; firstly, a variant can be directly assigned one of the defined classifications.

Additional classifications can be added via the Admin Controls (Admin Controls > Analysis > Classifications).

	Sample	HGVSc (Gene Symbol)	Chromosome	Start	End	Ref	Alt	Allele Frequency	Type	Ref Depth	Alt Depth	Total Depth	Ref Quality	Alt Quality	Quality Score
	5881	DVL1:c.*347T>C	1	1335795	1335795	A	G	100%	SNV	0	138	138	0.00	4771.00	4122.5
	5881	DVL1:c.366A>G	1	1342153	1342153	T	C	99.44%	SNV	1	232	233	32.00	7624.00	6750.6
	5881	ATAD3A:c.671G>A	1	1520153	1520153	G	A			185	177	362	6194.00	6050.00	4259.0
	5881	ATAD3A:c.*154C>T	1	1534226	1534226	C	T			21	16	37	738.00	515.00	275.1
	5881	ATAD3A:c.*330G>C	1	1534402	1534402	G	C			183	191	374	6338.00	6565.00	4665.6
	5881	SKI:c.*1802T>C	1	2308567	2308567	T	C			0	470	470	0.00	15743.00	14005.0
	5881	SKI:c.*3182dup	1	2309937	2309937	C	CT			22	137	166	752.00	8035.00	6276.5
	5881	CAMTA1:c.288C>T	1	6887707	6887707	C	T			0	123	123	6523.00	4644.00	3142.3
	5881	CAMTA1:c.*419T>C	1	6888144	6888144	T	C			7148	6375	13523	7148.00	6375.00	4498.8
	5881	CAMTA1:c.1350G>A	1	7663897	7663897	G	A			6530	7083	13613	6530.00	7083.00	5107.8
	5881	CAMTA1:c.*890T>A	1	7767381	7767381	T	A			5559	5895	11454	5559.00	5895.00	4270.7
	5881	CAMTA1:c.*1364_1365del	1	7767843	7767843	C	CAA	94.44%	Insertion	291	8624	8915	291.00	8624.00	7109.4
	5881	CAMTA1:c.*1640del	1	7768113	7768114	TA	T	70%	Deletion	856	3439	4295	856.00	3439.00	2613.2
	5881	CAMTA1:c.*1738del	1	7768218	7768219	CA	C	92.83%	Deletion	537	12239	12776	537.00	12239.00	9947.6
	5881	CAMTA1:c.*2522del	1	7769003	7769004	AT	A	55.33%	Deletion	5160	11149	16309	5160.00	11149.00	8450.7
	5881	CAMTA1:c.*2556G>A	1	7769047	7769047	G	A	50.38%	SNV	6783	7084	13867	6783.00	7084.00	5113.6
	5881	CAMTA1:c.*2725A>G	1	7769226	7769226	A	G	45.71%	SNV	6045	4987	11032	6045.00	4987.00	3481.8
	5881	REFSeq: *2212G>A	1	8352875	8352875	C	T	51.22%	SNV	180	189	369	6259.00	6475.00	4658.6

DVL1:c.366A>G

Add to Shortlist

Add to...

Set Preferred Transcript

Classify

Variant Links...

Notes

Using Guidelines...

Clear Classification

Benign

Uncertain significance: likely benign

Uncertain significance

Uncertain significance: likely pathogenic

Pathogenic

Figure 130: Variant classification options

A variant classification is selected from the list that is included by default. These are:

- Benign
- Uncertain significance, likely benign
- Uncertain significance
- Uncertain significance, likely pathogenic
- Pathogenic

	Sample	HGVSc (Gene Symbol)	Chromosome	Start	End	Ref	Alt	Allele Frequency	Type	Ref Depth	Alt Depth	Total Depth	Ref Quality	Alt Quality	Quality Score
	5881	DVL1:c.*347T>C	1	1335795	1335795	A	G	100%	SNV	0	138	138	0.00	4771.00	4122.5
	5881	DVL1:c.366A>G	1	1342153	1342153	T	C	99.37%	SNV	1	232	233	32.00	7624.00	6750.6
	5881	ATAD3A:c.671G>A	1	1520153	1520153	G	A			185	177	362	6194.00	6050.00	4259.0
	5881	ATAD3A:c.*154C>T	1	1534226	1534226	C	T			21	16	37	738.00	515.00	275.1
	5881	ATAD3A:c.*330G>C	1	1534402	1534402	G	C			183	191	374	6338.00	6565.00	4665.6
	5881	SKI:c.*1802T>C	1	2308567	2308567	T	C			0	470	470	0.00	15743.00	14005.
	5881	SKI:c.*3182dup	1	2309937	2309937	C	CT			22	137	166	752.00	8035.00	6276.5
	5881	CAMTA1:c.288C>T	1	6887707	6887707	C	T			0	123	123	6523.00	4644.00	3142.3
	5881	CAMTA1:c.*419T>C	1	6888144	6888144	T	C			7148.00	6375.00	4498.8			
	5881	CAMTA1:c.1350G>A	1	7663897	7663897	G	A			6530.00	7083.00	5107.8			
	5881	CAMTA1:c.*890T>A	1	7767381	7767381	T	A			5559.00	5895.00	4270.7			
	5881	CAMTA1:c.*1364_1365del	1	7767843	7767843	C	CAA			291.00	8624.00	7109.4			
	5881	CAMTA1:c.*1640del	1	7768113	7768114	TA	T	70%	Deletion	856.00	3439.00	2613.2			
	5881	CAMTA1:c.*1738del	1	7768218	7768219	CA	C	92.83%	Deletion	537.00	12239.00	9947.6			
	5881	CAMTA1:c.*2522del	1	7769003	7769004	AT	A	55.33%	Deletion	5160.00	11149.00	8450.7			
	5881	CAMTA1:c.*2556G>A	1	7769047	7769047	G	A	50.38%	SNV	6783.00	7084.00	5113.6			
	5881	CAMTA1:c.*2735A>G	1	7769226	7769226	A	G	45.71%	SNV	6045.00	4987.00	3481.8			
	5881	REFSeq: *2212G>A	1	8352875	8352875	C	T	51.22%	SNV	180	189	369	6259.00	6475.00	4658.6

DVL1:c.366A>G

Add to Shortlist

Add to...

Set Preferred Transcript

Classify

Variant Links...

Notes

Benign

Uncertain significance: likely benign

Uncertain significance

Uncertain significance: likely pathogenic

Pathogenic

Figure 131: Annotating a variant as benign

Once this classification has been made the variant will be annotated with the corresponding colour classification. This colour can be changed in the Admin Controls section of the software (Admin Controls > Analysis > Classifications)

This update will be applied to the variant annotation. As a result where the same variant appears in other samples it will have the same colour coding in the table.

	Sample	HGVSc (Gene Symbol)	Chromosome	Start	End	Ref	Alt	Allele Frequency	Type	Ref Depth	Alt Depth	Total Depth	Ref Quality	Alt Quality	Quality Score
	5881	DVL1:c.*347T>C	1	1335795	1335795	A	G	100%	SNV	0	138	138	0.00	4771.00	4122.59
	5881	DVL1:c.366A>G	1	1342153	1342153	T	C	99.57%	SNV	1	232	233	32.00	7624.00	6750.47
	5881	ATAD3A:c.671G>A	1	1520153	1520153	G	A	48.9%	SNV	185	177	362	6194.00	6050.00	4259.04

Figure 132: Update of the annotation to show a variant as benign

A variant classification can be removed using the clear classification method

	Sample	HGVSc (Gene Symbol)	Chromosome	Start	End	Ref	Alt	Allele Frequency	Type	Ref Depth	Alt Depth	Total Depth	Ref Quality	Alt Quality	Quality Score
	5881	DVL1:c.*347T>C	1	1335795	1335795	A	G	100%	SNV	0	138	138	0.00	4771.00	4122.59
	5881	DVL1:c.366A>G	1	1342153	1342153	T	C	99.57%	SNV	1	232	233	32.00	7624.00	6750.47
	5881	ATAD3A:c.671G>A	1	1520153	1520153	G	A	48.9%	SNV	185	177	362	6194.00	6050.00	4259.04
	5881	ATAD3A:c.*154C>T	1	1534226	1534226	C	T	21	16	37	738.00	515.00	275.1
	5881	ATAD3A:c.*390G>C	1	1534402	1534402	G	C	183	191	374	6338.00	6565.00	4665.6
	5881	SKI:c.*1802T>C	1	2308567	2308567	T	C	0	470	470	0.00	15743.00	14005.0
	5881	SKI:c.*3182dup	1	2309937	2309937	C	CT	22	137	166	752.00	8035.00	6276.5
	5881	CAMTA1:c.288C>T	1	6887707	6887707	C	T	6523.00	4644.00	3142.3
	5881	CAMTA1:c.*419T>C	1	6888144	6888144	T	C	7148.00	6375.00	4498.8
	5881	CAMTA1:c.1350G>A	1	7663897	7663897	G	A	6530.00	7083.00	5107.8
	5881	CAMTA1:c.*890T>A	1	7767381	7767381	T	A	5559.00	5895.00	4270.7
	5881	CAMTA1:c.*1364_*1365del	1	7767843	7767843	C	CAA	94.44%	Insertion	291.00	8624.00	7109.4
	5881	CAMTA1:c.*1640del	1	7768113	7768114	TA	T	70%	Deletion	856.00	3439.00	2613.2
	5881	CAMTA1:c.*1738del	1	7768218	7768219	CA	C	92.83%	Deletion	537.00	12239.00	9947.6
	5881	CAMTA1:c.*2522del	1	7769003	7769004	AT	A	55.33%	Deletion	5160.00	11149.00	8450.7
	5881	CAMTA1:c.*2556G>A	1	7769047	7769047	G	A	50.38%	SNV	6783.00	7084.00	5113.6
	5881	CAMTA1:c.*2735A>G	1	7769226	7769226	A	G	45.71%	SNV	6045.00	4987.00	3481.8
	5881	RRF1:c.*2212G>A	1	8352875	8352875	C	T	51.22%	SNV	180	189	369	6759.00	6475.00	4658.6

Figure 133: Removing a variant classification annotation

The classification will be removed for the variant in the table and all other samples with the same variant will be similarly updated.

	Sample	HGVSc (Gene Symbol)	Chromosome	Start	End	Ref	Alt	Allele Frequency	Type	Ref Depth	Alt Depth	Total Depth	Ref Quality	Alt Quality	Quality Score
	5881	DVL1:c.*347T>C	1	1335795	1335795	A	G	100%	SNV	0	138	138	0.00	4771.00	4122.59
	5881	DVL1:c.366A>G	1	1342153	1342153	T	C	99.57%	SNV	1	232	233	32.00	7624.00	6750.47
	5881	ATAD3A:c.671G>A	1	1520153	1520153	G	A	48.9%	SNV	185	177	362	6194.00	6050.00	4259.04


Figure 134: Update of the annotation to show variant without any classification

Using the American College of Medical Genetics and Genomics (ACMG) Guidelines

An alternative means to derive a classification for a variant is via guidelines described by the ACMG. These guidelines are included with Interpret.

To follow the ACMG guidelines the user provides answers to a specific set of questions. Each answer will navigate the user through the conditions of the guidelines until a classification of the variant can be made.

Selecting to use the guidelines option leads to a new window opening.




Sample	HGVSc (Gene Symbol)	Chromosome	Start	End	Ref	Alt	Allele Frequency	Type	Ref Depth	Alt Depth	Total Depth	Ref Quality	Alt Quality	Quality
5881	DVL1:c.366A>G	1	1335795	1335795	A	G	100%	SNV	0	138	138	0.00	4771.00	4122.5
5881	ATAD3A:c.671G>A	1	1520153	1520153	G	A			185	177	362	6194.00	6050.00	4259.0
5881	ATAD3A:c.*154C>T	1	1534226	1534226	C	T			21	16	37	738.00	515.00	275.1
5881	ATAD3A:c.*330G>C	1	1534402	1534402	G	C			183	191	374	6338.00	6565.00	4665.6
5881	SKI:c.*1802T>C	1	2308567	2308567	T	C			0	470	470	0.00	15743.00	14005.
5881	SKI:c.*3182dup	1	2309937	2309937	C	CT			22	137	166	752.00	8035.00	6276.5
5881	CAMTA1:c.288C>T	1	6887707	6887707	C	T			400	133	533	6523.00	4644.00	3142.3
5881	CAMTA1:c.*419T>C	1	6888144	6888144	T	C						5559.00	5895.00	4270.7
5881	CAMTA1:c.1350G>A	1	7663897	7663897	G	A						291.00	8624.00	7109.4
5881	CAMTA1:c.*890T>A	1	7767381	7767381	T	A						856.00	3439.00	2613.2
5881	CAMTA1:c.*1364,*1365d	1	7767843	7767843	C	CAA	94.44%	Insertion				537.00	12239.00	9947.6
5881	CAMTA1:c.*1740del	1	7768113	7768114	TA	T	70%	Deletion				5160.00	11149.00	8450.7
5881	CAMTA1:c.*1738del	1	7768218	7768219	CA	C	92.83%	Deletion				6783.00	7084.00	5113.6
5881	CAMTA1:c.*2522del	1	7769003	7769004	AT	A	55.33%	Deletion				6045.00	4987.00	3481.6
5881	CAMTA1:c.*2556G>A	1	7769047	7769047	G	A	50.38%	SNV				6259.00	6475.00	4658.6
5881	CAMTA1:c.*2735A>G	1	7769226	7769226	A	G	45.71%	SNV	180	189	369			
5881	REFSeq:*2212G>A	1	8357875	8357875	C	T	51.72%	SNV						

Figure 135: Selection of the inbuilt ACMG classification guidelines

The initial ACMG window, shown below, consists of a progress bar that will report how close to a classification

- Progress Bar – Showing the progress of the classification.
- Questions – These are the questions to be answered.
- Toggle – Allowing the display to be a graph view or a table view.



Variant Classification - SKI:c.*3182dup in 5881

Back In Progress 0%

Is null variant?

Yes No

Evidence History

Classification Details

Decision Tree ACMG Guidelines

Variant [FNST00000378536.5:c.*3182dup](#)

Date Started 30-Jan-2020 16:05:32

Initiated By Admin User

Percentage Complete 0.00%

Provisional Classification Uncertain significance

Sub-Classifications

Graph View Table View Classify Variant

Figure 136: The initial ACMG classification window

As the user answers questions the progress bar will update.

Variant Classification - SKI:c.*3182dup in 5881

Back In Progress 4%

Is null variant? (Yes)

Is LOF a known mechanism of disease?

Question: Is LOF a known mechanism of disease?

Yes No

Evidence History

Classification Details

Decision Tree ACMG Guidelines

Variant [ENST00000378536.5:c.*3182dup](#)

Date Started 30-Jan-2020 16:05:32

Initiated By Admin User

Percentage Complete 3.85%

Provisional Classification Uncertain significance

Sub-Classifications

Graph View Table View Classify Variant

Figure 137: A classification at the start of the guidelines, with 4% progress

Variant Classification - SKI:c.*3182dup in 5881

Back In Progress 67%

Is null variant? (Yes)

Is LOF a known mechanism of disease? (Yes)

PVS1 De novo in a patient with the disease and no family history?

Question: De novo in a patient with the disease and no family history?

Yes No

Evidence History

Classification Details

Decision Tree ACMG Guidelines

Variant [ENST00000378536.5:c.*3182dup](#)

Date Started 30-Jan-2020 16:05:32

Initiated By Admin User

Percentage Complete 66.67%

Provisional Classification Uncertain significance

Sub-Classifications

PVS1

Graph View Table View Classify Variant

Figure 138: A classification with 67% progress

When sufficient questions have been answered to allow a classification the progress bar will update to show 100% and say Ready to be classified.

Variant Classification - SKI:c.*3182dup in 5881

Back Ready to be classified 100%

Question: De novo in a patient with the disease and no family history?

Yes No

Evidence History

Classification Details

Decision Tree ACMG Guidelines

Variant [ENST00000378536.5:c.*3182dup](#)

Date Started 30-Jan-2020 16:05:32

Initiated By Admin User

Percentage Complete 100%

Provisional Classification Uncertain significance

Sub-Classifications

Figure 139: A completed ACMG classification



A window will appear showing the classification and give the user the option of making the classification or cancelling.




Figure 140: A completed classification showing a Pathogenic all has been made

Selecting the classification will update the variant's annotation accordingly.

<input type="checkbox"/>	5881	SKI:c.*1802T>C	1	2308567	2308567	T	C	100%	SNV	0	470	470	0.00	15743.00	14005.00
<input checked="" type="checkbox"/>	5881	SKI:c.*3182dup	1	2309937	2309937	C	CT	86.16%	Insertion	22	137	166	752.00	8035.00	6276.50
<input type="checkbox"/>	5881	CAMTA1:c.288C>T	1	6887707	6887707	C	T	40.99%	SNV	190	132	322	6523.00	4644.00	3142.30

Figure 141: Updated annotation for the variant to show its status as pathogenic

It is possible to review the choices made in the guidelines; using the table view, users can see which questions were asked and how they were answered by whom and when.

Variant Classification - SKI:c.*3182dup in 5881				
<div> Back Completed </div>		<div> <div>100%</div> <div></div> </div>		Classified As: Pathogenic
Unanswered Questions			Answer	Question:
				Click on a Question
Answered Questions		Value	Answered By	Date Answered
1	De novo in a patient with the disease and no family history?	Yes	Admin User	30-Jan-2020 16:07:38
2	Is LOF a known mechanism of disease?	Yes	Admin User	30-Jan-2020 16:07:17
3	Is null variant?	Yes	Admin User	30-Jan-2020 16:07:00

Figure 142: Table view of a completed classification

Selecting a row from the table view allows a result to be modified if that is required. Alternatively, evidence can be added to support the answer to the question

Variant Classification - SKI:c.*3182dup in 5881				
<div> Back Completed </div>		<div> <div>100%</div> <div></div> </div>		Classified As: Pathogenic
Unanswered Questions			Answer	Question:
				Is LOF a known mechanism of disease?
				<div> Yes Edit Answer </div>
				<div> Evidence History </div>
Answered Questions		Value	Answered By	Date Answered
1	De novo in a patient with the disease and no family history?	Yes	Admin User	30-Jan-2020 16:07:38
2	Is LOF a known mechanism of disease?	Yes	Admin User	30-Jan-2020 16:07:17
3	Is null variant?	Yes	Admin User	30-Jan-2020 16:07:00

Figure 143: Reviewing an answer in the table view

Variant Links

The software allows users to link out to external sources of documentation. Currently included are:

- Ensembl
- ClinView
- ExAC

Additional resources can be added in the Admin Controls (Admin Controls > Analysis > Manage Links).

Sample	HGVSc (Gene Symbol)	Chromosome	Start	End	Ref	Alt	Allele Frequency	Type	Ref Depth	Alt Depth	Total Depth	Ref Quality	Alt Quality	Quality Score
5881	DVL1:c.347T>C	1	1335795	1335795	A	G	100%	SNV	0	138	138	0.00	4771.00	4122.59
5881	DVL1:c.366A>G	1	1342153	1342153	T	C	99.57%	SNV	1	233	233	32.00	7624.00	6750.47
5881	ATAD3A:c.671G>A	1	1520153	1520153	G	A	100%	SNV	185	177	362	6194.00	6050.00	4259.04
5881	ATAD3A:c.*154C>T	1	1534226	1534226	C	T	100%	SNV	21	16	37	738.00	515.00	275.18
5881	ATAD3A:c.*330G>C	1	1534402	1534402	G	C	100%	SNV	183	191	374	6338.00	6565.00	4665.67
5881	SKI:c.*1802T>C	1	2308567	2308567	T	C	100%	SNV	0	470	470	0.00	15743.00	14005.10
5881	SKI:c.*3182dup	1	2309937	2309937	C	CT	100%	SNV	22	137	166	752.00	8035.00	6276.51
5881	CAMTA1:c.288C>T	1	6887707	6887707	C	T	100%	SNV	190	132	322	6523.00	4644.00	3142.37
5881	CAMTA1:c.*419T>C	1	6887707	6887707	C	T	100%	SNV	203	184	387	7148.00	6375.00	4498.84
5881	CAMTA1:c.1350G>A	1	7663897	7663897	G	A	100%	SNV	36	400	436	6530.00	7083.00	5107.83
5881	CAMTA1:c.*890T>A	1	7767381	7767381	T	A	100%	SNV	88	325	413	5559.00	5895.00	4270.76
5881	CAMTA1:c.*1364_*1365del	1	7767843	7767843	C	CAA	94.44%	Insertion	53	224	277	291.00	8624.00	7109.48
5881	CAMTA1:c.*1640del	1	7768113	7768114	TA	T	70%	Deletion	3	124	127	856.00	3439.00	2613.25
5881	CAMTA1:c.*1738del	1	7768218	7768219	CA	C	92.83%	Deletion	37	239	276	537.00	12239.00	9947.63
5881	CAMTA1:c.*2522del	1	7769003	7769004	AT	A	55.33%	Deletion	37	346	383	5160.00	11149.00	8450.75
5881	CAMTA1:c.*2556G>A	1	7769047	7769047	G	A	50.38%	SNV	31	399	430	6783.00	7084.00	5113.03
5881	CAMTA1:c.*2735A>G	1	7769226	7769226	A	G	45.71%	SNV	44	315	359	6045.00	4987.00	3481.89
5881	RERE:c.*2212G>A	1	8352875	8352875	C	T	51.22%	SNV	89	369	458	6259.00	6475.00	4658.64
5881	RERE:c.*2158T>C	1	8352929	8352929	A	G	100%	SNV	0	331	331	0.00	11478.00	10231.50

Figure 144: Variant links available in the software

If a source is selected Interpret will show the information in a separate tab in the web browser.

Sample	HGVSc (Gene Symbol)	Chromosome	Start	End	Ref	Alt	Allele Frequency	Type	Ref Depth	Alt Depth	Total Depth	Ref Quality	Alt Quality	Quality Score
5881	DVL1:c.347T>C	1	1335795	1335795	A	G	100%	SNV	0	138	138	0.00	4771.00	4122.59
5881	DVL1:c.366A>G	1	1342153	1342153	T	C	99.57%	SNV	1	233	233	32.00	7624.00	6750.47
5881	ATAD3A:c.671G>A	1	1520153	1520153	G	A	100%	SNV	185	177	362	6194.00	6050.00	4259.04
5881	ATAD3A:c.*154C>T	1	1534226	1534226	C	T	100%	SNV	21	16	37	738.00	515.00	275.18
5881	ATAD3A:c.*330G>C	1	1534402	1534402	G	C	100%	SNV	183	191	374	6338.00	6565.00	4665.67
5881	SKI:c.*1802T>C	1	2308567	2308567	T	C	100%	SNV	0	470	470	0.00	15743.00	14005.10
5881	SKI:c.*3182dup	1	2309937	2309937	C	CT	100%	SNV	22	137	166	752.00	8035.00	6276.51
5881	CAMTA1:c.288C>T	1	6887707	6887707	C	T	100%	SNV	190	132	322	6523.00	4644.00	3142.37
5881	CAMTA1:c.*419T>C	1	6887707	6887707	C	T	100%	SNV	203	184	387	7148.00	6375.00	4498.84
5881	CAMTA1:c.1350G>A	1	7663897	7663897	G	A	100%	SNV	36	400	436	6530.00	7083.00	5107.83
5881	CAMTA1:c.*890T>A	1	7767381	7767381	T	A	100%	SNV	88	325	413	5559.00	5895.00	4270.76
5881	CAMTA1:c.*1364_*1365del	1	7767843	7767843	C	CAA	94.44%	Insertion	53	224	277	291.00	8624.00	7109.48
5881	CAMTA1:c.*1640del	1	7768113	7768114	TA	T	70%	Deletion	3	124	127	856.00	3439.00	2613.25
5881	CAMTA1:c.*1738del	1	7768218	7768219	CA	C	92.83%	Deletion	37	239	276	537.00	12239.00	9947.63
5881	CAMTA1:c.*2522del	1	7769003	7769004	AT	A	55.33%	Deletion	37	346	383	5160.00	11149.00	8450.75
5881	CAMTA1:c.*2556G>A	1	7769047	7769047	G	A	50.38%	SNV	31	399	430	6783.00	7084.00	5113.03
5881	CAMTA1:c.*2735A>G	1	7769226	7769226	A	G	45.71%	SNV	44	315	359	6045.00	4987.00	3481.89
5881	RERE:c.*2212G>A	1	8352875	8352875	C	T	51.22%	SNV	89	369	458	6259.00	6475.00	4658.64
5881	RERE:c.*2158T>C	1	8352929	8352929	A	G	100%	SNV	0	331	331	0.00	11478.00	10231.50

Figure 145: Selection of ExAC as an external resource for the gene selected



Figure 146: Example of the software linking out to an external data source, in this case the GnomAD for the gene containing the variant in Interpret

Add Notes

Interpret allows users to add notes for a variant and to also edit notes on the system. This is accessed through the Notes menu item.

	Sample ▾	HGVSc (Gene Symbol)	Chromosome ▾	Start ▾	End ▾	Ref ▾	Alt ▾	Allele Frequency ▾	Type ▾	Ref Depth ▾	Alt Depth ▾	Total Depth ▾	Ref Quality ▾	Alt Quality ▾	Quality Score ▾
<input type="checkbox"/>	5881	DVL1:c.347T>C	1	1335795	1335795	A	G	100%	SNV	0	138	138	0.00	4771.00	4122.59
<input checked="" type="checkbox"/>	5881	DVL1:c.366A>G	1	1342153	1342153	T	C		<div>DVL1:c.366A>G</div> <div>Add to Shortlist</div> <div>Add to.. ></div> <div>Set Preferred Transcript ></div> <div>Classify ></div> <div>Variant Links.. ></div> <div>Notes ></div>	1	232	233	32.00	7624.00	6750.47
<input type="checkbox"/>	5881	ATAD3A:c.671G>A	1	1520153	1520153	G	A	185		177	362	6194.00	6050.00	4259.04	
<input type="checkbox"/>	5881	ATAD3A:c.*154C>T	1	1534226	1534226	C	T	21		16	37	738.00	515.00	275.18	
<input type="checkbox"/>	5881	ATAD3A:c.*330G>C	1	1534402	1534402	G	C	183		191	374	6338.00	6565.00	4665.67	
<input type="checkbox"/>	5881	SKI:c.*1802T>C	1	2308567	2308567	T	C	0		470	470	0.00	15743.00	14005.10	
<input type="checkbox"/>	5881	SKI:c.*3182dup	1	2309937	2309937	C	CT	22		137	166	752.00	8035.00	6276.51	
<input type="checkbox"/>	5881	CAMTA1:c.288C>T	1	6887707	6887707	C	T	190		132	322	6523.00	4644.00	3142.37	
<input type="checkbox"/>	5881	CAMTA1:c.*419T>C	1	6888144	6888144	T	C	203		184	387	7148.00	6375.00	4498.84	
<input type="checkbox"/>	5881	CAMTA1:c.1350G>A	1	7663897	7663897	G	A	194		206	400	6530.00	7083.00	5107.83	
<input type="checkbox"/>	5881	CAMTA1:c.*890T>A	1	7767381	7767381	T	A	168		325	5559.00	5895.00	4270.76		
<input type="checkbox"/>	5881	CAMTA1:c.*1364_*1365del	1	7767843	7767843	C	CAA	94.44%	Insertion	153	224	291.00	8624.00	7109.48	
<input type="checkbox"/>	5881	CAMTA1:c.*1640del	1	7768113	7768114	TA	T	70%	Deletion	27	63	124	856.00	3439.00	2613.25

Figure 147: Adding a note to a variant

Selecting the Add Note will generate a popup window




Figure 148: The note template window for the selected variant

Users can add the required text, up to 250 characters, in the text box




Figure 149: Addition of a note to a variant

Selecting Add Now will append the note to the variant.




Figure 150: An example of a note on the system

The additional text will now be displayed




Figure 151: Appending text to an existing notation

Notes can be modified by clicking on the pen icon. This makes the text box editable




Figure 152: Adding an update to a note

Once any update has been made, selecting Apply will incorporate the changes.




Figure 153: A note showing the updated annotation

Similarly, a note can be deleted through the red bin icon.

Users are asked to confirm the delete request after which the note will be removed.




Figure 154: Deleting a variant note

Where there is a note for a variant the note can be viewed through the Notes options seen when right clicking on the variant.

Sample	HGVS (Gene Symbol)	Chromosome	Start	End	Ref	Alt	Allele Frequency	Type	Ref Depth	Alt Depth	Total Depth	Ref Quality	Alt Quality	Quality Score
5881	DVL1:c.366A>G	1	1335795	1335795	A	G	100%	SNV	0	138	138	0.00	4771.00	4122.59
5881	ATAD3A:c.671G>A	1	1520153	1520153	G	A			185	177	362	6194.00	6050.00	4259.04
5881	ATAD3A:c.*154C>T	1	1534226	1534226	C	T			21	16	37	738.00	515.00	275.18
5881	ATAD3A:c.*330G>C	1	1534402	1534402	G	C			183	191	374	6338.00	6565.00	4665.67
5881	SKI:c.*1802T>C	1	2308567	2308567	T	C			0	470	470	0.00	15743.00	14005.10
5881	SKI:c.*3182dup	1	2309937	2309937	C	CT			22	137	166	752.00	8035.00	6276.51
5881	CAMTA1:c.288C>T	1	6887707	6887707	C	T			190	132	322	6523.00	4644.00	3142.37
5881	CAMTA1:c.*419T>C	1	6888144	6888144	T	C			203	184	387	7148.00	6375.00	4498.84
5881	CAMTA1:c.1350G>A	1	7663897	7663897	G	A			194	206	400	6530.00	7083.00	5107.83
5881	CAMTA1:c.*890T>A	1	7767381	7767381	T	A			157	168	325	5559.00	5895.00	4270.76
5881	CAMTA1:c.*1364_*1365d	1	7767843	7767843	C	CAA			9	153	224	291.00	8624.00	7109.48
5881	CAMTA1:c.*1640del	1	7768113	7768114	TA	T			27	63	124	856.00	3439.00	2613.25
5881	CAMTA1:c.*1738del	1	7768218	7768219	CA	C			16	207	239	537.00	12239.00	9947.63
5881	CAMTA1:c.*2522del	1	7769003	7769004	AT	A			151	187	346	5160.00	11149.00	8450.75

Figure 155: Selecting a Note to view

The note is displayed on the screen.

Sample	HGVS (Gene Symbol)	Chromosome	Start	End	Ref	Alt	Allele Frequency	Type	Ref Depth	Alt Depth	Total Depth	Ref Quality	Alt Quality	Quality Score
5881	DVL1:c.366A>G	1	1335795	1335795	A	G	100%	SNV	0	138	138	0.00	4771.00	4122.59
5881	ATAD3A:c.671G>A	1	1520153	1520153	G	A			185	177	362	6194.00	6050.00	4259.04
5881	ATAD3A:c.*154C>T	1	1534226	1534226	C	T			21	16	37	738.00	515.00	275.18
5881	ATAD3A:c.*330G>C	1	1534402	1534402	G	C			183	191	374	6338.00	6565.00	4665.67
5881	SKI:c.*1802T>C	1	2308567	2308567	T	C			0	470	470	0.00	15743.00	14005.10
5881	SKI:c.*3182dup	1	2309937	2309937	C	CT			22	137	166	752.00	8035.00	6276.51
5881	CAMTA1:c.288C>T	1	6887707	6887707	C	T			190	132	322	6523.00	4644.00	3142.37
5881	CAMTA1:c.*419T>C	1	6888144	6888144	T	C			203	184	387	7148.00	6375.00	4498.84
5881	CAMTA1:c.1350G>A	1	7663897	7663897	G	A			194	206	400	6530.00	7083.00	5107.83
5881	CAMTA1:c.*890T>A	1	7767381	7767381	T	A			157	168	325	5559.00	5895.00	4270.76
5881	CAMTA1:c.*1364_*1365d	1	7767843	7767843	C	CAA			9	153	224	291.00	8624.00	7109.48
5881	CAMTA1:c.*1640del	1	7768113	7768114	TA	T			27	63	124	856.00	3439.00	2613.25
5881	CAMTA1:c.*1738del	1	7768218	7768219	CA	C			16	207	239	537.00	12239.00	9947.63
5881	CAMTA1:c.*2522del	1	7769003	7769004	AT	A			151	187	346	5160.00	11149.00	8450.75
5881	CAMTA1:c.*2756G>A	1	7769047	7769047	G	A			198	201	399	6783.00	7084.00	5113.03
5881	REER:c.*2212G>A	1	8352875	8352875	C	T			171	144	315	6045.00	4987.00	3481.89
5881	REER:c.*2158T>C	1	8352929	8352929	A	G			180	189	369	6259.00	6475.00	4658.64
5881	REER:c.*1830T>G	1	8353257	8353257	A	C			0	331	331	0.00	11478.00	10231.50

Viewing CNV and LOH Events

The variant table has a column selector icon  allowing user to configure which columns are displayed. The figure below shows the columns available for display.




Figure 156: Columns available to select for display in the CNV/LOH variants page

Selecting a variant will show it in IGV, a track for both CNV and LOH will displayed. Sometimes there will only be a CNV call as in the example below.




Figure 157: Example of a CNV call only

Sometimes there will only be a LOH call




Figure 158: Example of a LOH call only

Sometimes there will be CNV and LOH calls




Figure 159: Example of a sample with a CNV call and a LOH call in the same genomic location

CNV and LOH Options

As with the page displaying SNV and Indel calls there are options available for each variant called by the CNV/LOH pipeline,
Right clicking on a variant will provide a menu of the possible options.

	Chromosome	Start	End	Type	Length	Copy Number	# Markers	Mean	Confidence	Overlap	Genome Build	Classification	Depth	Frequency	Sample
<input checked="" type="checkbox"/>	1	152305335	152305635	Deletion	300b	0	6	-4.20985	High	CDS (target)	GRCh38	Unclassified	31		8210
<input type="checkbox"/>	21	14110138	17533897	Deletion	3.42Mb	1	49	-1.01204	High	CDS (other)	GRCh38	Unclassified	123		8210
<input type="checkbox"/>	7	5973240	5973540	Deletion	300b	1	6	-1.10507	High	CDS (target)	GRCh38	Unclassified	90		8210
<input type="checkbox"/>	16	2102058	2102258	Deletion	200b	1	6	-0.7722	High	CDS (target)	GRCh38	Unclassified	206		8210
<input type="checkbox"/>	7	94597891	94598041	Deletion	150b	1	6	-1.26924	Low	CDS (target)	GRCh38	Unclassified	13		8210
<input type="checkbox"/>	4	44967	3843401	Deletion	3.8Mb	1	43	-1.01016	High	CDS (target)	GRCh38	Unclassified	182	0.62	8210
<input type="checkbox"/>	17	46170855	46225432	Deletion	54.58kb	1	8	-0.714658	High	CDS (target)	GRCh38	Unclassified	579	0.47	8210
<input type="checkbox"/>	6	257160	335276	Deletion	78.12kb	1	8	-0.80957	High	CDS (other)	GRCh38	Unclassified	229		8210
<input type="checkbox"/>	1	152308235	152309085	Duplication	850b	3	17	0.651716	High	CDS (target)	GRCh38	Unclassified	1386	0.46	8210
<input type="checkbox"/>	22	50720285	50720935	Duplication	650b	3	12	0.401733	High	CDS (target)	GRCh38	Unclassified	302		8210

Figure 160: Available options for a selected CNV or LOH call

Adding to a shortlist
Variants added to a shortlist are annotated with a tick

	Chromosome	Start	End	Type	Length	Copy Number	# Markers	Mean	Confidence	Overlap	Genome Build	Classification	Depth	Frequency	Sample
<input checked="" type="checkbox"/>	1	152305335	152305635	Deletion	300b	0	6	-4.20985	High	CDS (target)	GRCh38	Unclassified	31		8210
<input checked="" type="checkbox"/>	21	14110138	17533897	Deletion	3.42Mb	1	49	-1.01204	High	CDS (other)	GRCh38	Unclassified	123		8210
<input type="checkbox"/>	7	5973240	5973540	Deletion	300b	1	6	-1.10507	High	CDS (target)	GRCh38	Unclassified	90		8210

Figure 161: Variants added to the shortlist displayed

Shortlisted variants can be viewed.

Chromosome	Start	End	Type	Length	Copy Number	# Markers	Mean	Confidence	Overlap	Genotype	Sample
1	152305335	152305635	Deletion	300b	0	6	-4.20985	High	CDS (target)	G	8210
21	14110138	17533897	Deletion	3.42Mb	1	49	-1.01204	High	CDS (other)	G	8210
7	5973240	5973540	Deletion	300b	1	6	-1.10507	High	CDS (target)	G	8210
16	2102058	2102258	Deletion	200b	1	4	-0.7722	High	CDS (target)	G	8210
7	94597891	94598041	Deletion	150b	1	3	-1.28924	Low	CDS (target)	G	8210
4	44967	3843401	Deletion	3.8Mb	1	399	-1.01016	High	CDS (target)	GRCh38	8210

Figure 162: Accessing the shortlist of selected variants

The shortlist opens in a separate view. A variant can be removed from the shortlist by clicking on the red bin icon in the shortlist view.

Chr	Start	End	Type	Length	Copy Number	#Markers	Mean	Confidence	Overlaps	
1	152305335	152305635	DEL	300	0	6	-4.20985	High	CDSTarget	
21	14110138	17533897	DEL	3423759	1	49	-1.01204	High	CDSOther	

Figure 163: Viewing the shortlist of CNV or LOH variants

Alternatively, a variant can be removed from the shortlist using the CNV options menu.

Chromosome	Start	End	Type	Length	Copy Number	# Markers	Mean	Confidence	Overlap	Genome Build	Classification	Depth	Frequency	Sample
6	257160	335276	Deletion	78.12kb	1	8	-0.80957	High	CDS (other)	GRCh38	Unclassified	229		8210
1	152305335	152305635	Deletion	300b	0	6	-4.20985	High	CDS (target)	GRCh38	Unclassified	31		8210
21	14110138	17533897	Deletion	3.42Mb	1	49	-1.01204	High	CDS (other)	GRCh38	Unclassified	123		8210
7	5973240	5973540	Deletion	300b	1	6	-1.10507	High	CDS (target)	GRCh38	Unclassified	90		8210
16	2102058	2102258	Deletion	200b	1	4	-0.7722	High	CDS (target)	GRCh38	Unclassified	206		8210
7	94597891	94598041	Deletion	150b	1	3	-1.28924	Low	CDS (target)	GRCh38	Unclassified	13		8210
4	44967	3843401	Deletion	3.8Mb	1	399	-1.01016	High	CDS (target)	GRCh38	Unclassified	182	0.62	8210
17	46170855	46225432	Deletion	54.58kb	1	43	-0.714658	High	CDS (target)	GRCh38	Unclassified	579	0.47	8210

Figure 164: Deleting a variant from the shortlist

The shortlist will be updated to reflect the removal of a variant.

Chr	Start	End	Type	Length	Copy Number	#Markers	Mean	Confidence	Overlaps	
1	152305335	152305635	DEL	300	0	6	-4.20985	High	CDSTarget	

Figure 165: The shortlist showing that the variant has been removed

Variant Classification

A variant can be classified from the list that is included by default. These are:

- Benign
- Uncertain significance, likely benign
- Uncertain significance
- Uncertain significance, likely pathogenic
- Pathogenic

Additional classifications can be added in the Admin Controls section of the software (Admin Controls > Analysis > Classifications)

Chromosome	Start	End	Type	Length	Copy Number	# Markers	Mean	Confidence	Overlap	Genome Build	Classification	Depth	Frequency	Sample
21	14110138	17533897	Deletion	3.42Mb	1	49	-1.01204	High	CDS (other)	GRCh38	Unclassified	123		8210
7	5973240	5973540	Deletion	300b	1	6	-1.10507	High	CDS (target)	GRCh38	Unclassified	90		8210
16	2102058	2102258	Deletion	200b						GRCh38	Unclassified	206		8210
7	94597891	94598041	Deletion	150b						GRCh38	Unclassified	13		8210
4	44967	3843401	Deletion	3.8Mb						GRCh38	Unclassified	182	0.62	8210
17	46170855	46225432	Deletion	54.58kb						GRCh38	Unclassified	579	0.47	8210
6	257160	335276	Deletion	78.12kb	1	8				GRCh38	Unclassified	229		8210
1	152305335	152305635	Deletion	300b	0	6				GRCh38	Unclassified	31		8210
1	152308235	152309085	Duplication	850b	3	17				GRCh38	Unclassified	1386	0.46	8210
22	50720285	50720935	Duplication	650b	3	12				GRCh38	Unclassified	302		8210

Figure 166: Default classifications

Chromosome	Start	End	Type	Length	Copy Number	# Markers	Mean	Confidence	Overlap	Genome Build	Classification	Depth	Frequency	Sample
21	14110138	17533897	Deletion	3.42Mb	1	49	-1.01204	High	CDS (other)	GRCh38	Pathogenic	123		8210
7	5973240	5973540	Deletion	300b	1	6	-1.10507	High	CDS (target)	GRCh38	Unclassified	90		8210
16	2102058	2102258	Deletion	200b						GRCh38	Unclassified	206		8210
7	94597891	94598041	Deletion	150b						GRCh38	Unclassified	13		8210
4	44967	3843401	Deletion	3.8Mb						GRCh38	Unclassified	182	0.62	8210
17	46170855	46225432	Deletion	54.58kb						GRCh38	Unclassified	579	0.47	8210
6	257160	335276	Deletion	78.12kb	1	8				GRCh38	Unclassified	229		8210
1	152305335	152305635	Deletion	300b	0	6				GRCh38	Unclassified	31		8210
1	152308235	152309085	Duplication	850b	3	17				GRCh38	Unclassified	1386	0.46	8210
22	50720285	50720935	Duplication	650b	3	12				GRCh38	Unclassified	302		8210

Figure 167: Classifying a CNV deletion as pathogenic

Chromosome	Start	End	Type	Length	Copy Number	# Markers	Mean	Confidence	Overlap	Genome Build	Classification	Depth	Frequency	Sample
21	14110138	17533897	Deletion	3.42Mb	1	49	-1.01204	High	CDS (other)	GRCh38	Pathogenic	123		8210
7	5973240	5973540	Deletion	300b	1	6	-1.10507	High	CDS (target)	GRCh38	Unclassified	90		8210

Figure 168: Updating of the variant to show the new classification

Chromosome	Start	End	Type	Length	Copy Number	# Markers	Mean	Confidence	Overlap	Genome Build	Classification	Depth	Frequency	Sample
21	14110138	17533897	Deletion	3.42Mb	1	49	-1.01204	High	CDS (other)	GRCh38	Pathogenic	123		8210
7	5973240	5973540	Deletion	300b	1	6	-1.10507	High	CDS (target)	GRCh38	Unclassified	90		8210
16	2102058	2102258	Deletion	200b						GRCh38	Unclassified	206		8210
7	94597891	94598041	Deletion	150b						GRCh38	Unclassified	13		8210
4	44967	3843401	Deletion	3.8Mb						GRCh38	Unclassified	182	0.62	8210
17	46170855	46225432	Deletion	54.58kb						GRCh38	Unclassified	579	0.47	8210
6	257160	335276	Deletion	78.12kb	1	8				GRCh38	Unclassified	229		8210
1	152305335	152305635	Deletion	300b	0	6				GRCh38	Unclassified	31		8210
1	152308235	152309085	Duplication	850b	3	17				GRCh38	Unclassified	1386	0.46	8210
22	50720285	50720935	Duplication	650b	3	12				GRCh38	Unclassified	302		8210

Figure 169: Removing a variant classification


View Classification History

User can review the classification of a variant by selecting that option in the menu.

Sample	Chromosome	Start	End	Type	Length	Copy Number	# Markers	Mean	Confidence	Overlap	Genome Build	Bands	Classification
Sample-1	13	50566440	51782998	Duplication	1.22Mb	3	70	0.203889	High	CDS (target)	GRCh37	q14.2-q14.3	Unclassified
Sample-1	X	100611029	100611266	Duplication	237b	3	3	0.632296	High	CDS (target)	GRCh37	q22.1	Unclassified

Figure 170: Selecting view classification history option

When chosen a table appears displaying how a variant has been classified, who made the classification and when any changes were made



Classification	User	Date
Unclassified	admin	Oct 21, 2021 4:13:44 PM
Pathogenic	admin	Oct 21, 2021 3:46:36 PM
Unclassified	admin	Oct 15, 2021 9:18:38 AM
Pathogenic	admin	Sep 21, 2021 2:27:30 PM

Figure 171: An example of a variant's classification history

View in CytoSure Interpret

CytoSure Interpret is OGT's class-leading microarray software analysis platform. For existing microarray customers, CNV and LOH events can be loaded into CytoSure Interpret.

<input type="checkbox"/>	Chromosome	Start	End	Type	Length	Copy Number	# Markers	Mean	Confidence	Overlap	Genome Build	Classification	Depth	Frequency	Sample
<input checked="" type="checkbox"/>	21	14110138	17533897	Deletion	3.42Mb	1	45	-0.01204	High	CDS (other)	GRCh38	Unclassified	123		8210
<input type="checkbox"/>	7	5973240	5973540	Deletion	300b	1	1	-0.10507	High	CDS (target)	GRCh38	Unclassified	90		8210
<input type="checkbox"/>	16	2102058	2102258	Deletion	200b	1	1	-0.07722	High	CDS (target)	GRCh38	Unclassified	206		8210
<input type="checkbox"/>	7	94597891	94598041	Deletion	150b	1	1	-0.026924	Low	CDS (target)	GRCh38	Unclassified	13		8210
<input type="checkbox"/>	4	44967	3843401	Deletion	3.8Mb	1	1	-0.01016	High	CDS (target)	GRCh38	Unclassified	182	0.62	8210
<input type="checkbox"/>	17	46170855	46225432	Deletion	54.58Kb	1	1	-0.01016	High	CDS (target)	GRCh38	Unclassified	579	0.47	8210
<input type="checkbox"/>	6	257160	335276	Deletion	78.12Kb	1	8	-0.80957	High	CDS (other)	GRCh38	Unclassified	229		8210
<input type="checkbox"/>	1	152305335	152305635	Deletion	300b	0	6	-4.20985	High	CDS (target)	GRCh38	Unclassified	31		8210
<input type="checkbox"/>	1	152308235	152309085	Duplication	850b	3	17	0.651716	High	CDS (target)	GRCh38	Unclassified	1386	0.46	8210
<input type="checkbox"/>	22	50720285	50720935	Duplication	650b	3	12	0.401733	High	CDS (target)	GRCh38	Unclassified	302		8210

Figure 172: Selecting to view a CNV deletion in CytoSure Interpret microarray software

It is necessary to have CytoSure Interpret open prior to selecting this option. If it is not yet running Interpret will issue a prompt to the user.





Figure 173: Prompt from Interpret if trying to load data in CytoSure Interpret when it is not running

Variant Links

The software allows users to link out to external sources of documentation. Currently included are:


- Ensembl

Additional resources can be added in the Admin Controls (Admin Controls > Analysis > Manage Links).



	Chromosome	Start	End	Type	Length	Copy Number	# Markers	Mean	Confidence	Overlap	Genome Build	Classification	Depth	Frequency	Sample
<input checked="" type="checkbox"/>	21	14110138	17533897	Deletion	3.42Mb	1	4	-1.01204	High	CDS (other)	GRCh38	Unclassified	123		8210
<input type="checkbox"/>	7	5973240	5973540	Deletion	300b			-1.0507	High	CDS (target)	GRCh38	Unclassified	90		8210
<input type="checkbox"/>	16	2102058	2102258	Deletion	200b			-0.7722	High	CDS (target)	GRCh38	Unclassified	206		8210
<input type="checkbox"/>	7	94597891	94598041	Deletion	150b			-1.26924	Low	CDS (target)	GRCh38	Unclassified	13		8210
<input type="checkbox"/>	4	44967	3843401	Deletion	3.8Mb			-1.01016	High	CDS (target)	GRCh38	Unclassified	182	0.62	8210
<input type="checkbox"/>	17	46170855	46225432	Deletion	54.58Kb			-1.01016	High	CDS (target)	GRCh38	Unclassified	579	0.47	8210
<input type="checkbox"/>	6	257160	335276	Deletion	78.12Kb	1	8	-4.20985	High	CDS (other)	GRCh38	Unclassified	229		8210
<input type="checkbox"/>	1	152305335	152305635	Deletion	300b	0	6	-1.01016	High	CDS (target)	GRCh38	Unclassified	31		8210
<input type="checkbox"/>	1	152308235	152309085	Duplication	850b	3	17	0.651716	High	CDS (target)	GRCh38	Unclassified	1386	0.46	8210
<input type="checkbox"/>	22	50720285	50720935	Duplication	650b	3	12	0.401733	High	CDS (target)	GRCh38	Unclassified	302		8210

Figure 174: Accessing variant links




	Chromosome	Start	End	Type	Length	Copy Number	# Markers	Mean	Confidence	Overlap	Genome Build	Classification	Depth	Frequency	Sample
<input checked="" type="checkbox"/>	21	14110138	17533897	Deletion	3.42Mb	1	4	-1.01204	High	CDS (other)	GRCh38	Unclassified	123		8210
<input type="checkbox"/>	7	5973240	5973540	Deletion	300b			-1.0507	High	CDS (target)	GRCh38	Unclassified	90		8210
<input type="checkbox"/>	16	2102058	2102258	Deletion	200b			-0.7722	High	CDS (target)	GRCh38	Unclassified	206		8210
<input type="checkbox"/>	7	94597891	94598041	Deletion	150b			-1.26924	Low	CDS (target)	GRCh38	Unclassified	13		8210
<input type="checkbox"/>	4	44967	3843401	Deletion	3.8Mb			-1.01016	High	CDS (target)	GRCh38	Unclassified	182	0.62	8210
<input type="checkbox"/>	17	46170855	46225432	Deletion	54.58Kb			-1.01016	High	CDS (target)	GRCh38	Unclassified	579	0.47	8210
<input type="checkbox"/>	6	257160	335276	Deletion	78.12Kb	1	8	-4.20985	High	CDS (other)	GRCh38	Unclassified	229		8210
<input type="checkbox"/>	1	152305335	152305635	Deletion	300b	0	6	-1.01016	High	CDS (target)	GRCh38	Unclassified	31		8210
<input type="checkbox"/>	1	152308235	152309085	Duplication	850b	3	17	0.651716	High	CDS (target)	GRCh38	Unclassified	1386	0.46	8210
<input type="checkbox"/>	22	50720285	50720935	Duplication	650b	3	12	0.401733	High	CDS (target)	GRCh38	Unclassified	302		8210

Figure 175: Accessing EnSEMBL as an external source for data annotation

Adding Notes to CNVs

Users can add notes to CNVs



	Sample	Chromosome	Start	End	Type	Length	Copy Number	# Markers	Mean	Confidence	Overlap	Genome Build	Bands	Classification	# Samples	% Sat
<input checked="" type="checkbox"/>	Sample-1	13	50566440	51782998	Deletion	1.21Mb	3	70	0.203889	High	CDS (target)	GRCh37	q14.2;q14.3	Unclassified	5	16.22
<input type="checkbox"/>	Sample-1	X	100611029	100611266			3	3	0.632296	High	CDS (target)	GRCh37	q22.1	Unclassified	5	35.14

Figure 176: Selecting the option to add notes to a CNV

When chosen a text editor is displayed as well any pre-existing notes. At the top is an option to choose a file and to then upload it. Below is the text box where details can be entered.

Figure 177: A blank template for creating a note

Interpret User Guide v96-20250820121007

Figure 178: Addition of text and a file to a note

Selecting + Add Note completes creation of the note and it is added to the existing notes section. Any file that has been uploaded with the note is shown and can be downloaded if required.




Figure 179: A new note is shown

Now when the user accesses the menu for the variant there is an additional option providing the display of any notes.




Figure 180: Existing notes are available to view

If selected, the note(s) are shown in a separate box.




Figure 181: Viewing an existing note

Additional notes can be added as shown below.




Figure 182: A blank note template with two existing notes

Manual Creation of CNVs

It may be that the user believes, based on the visual representation of the CNV data, that the software has missed a CNV call and would like to manually generate it. For example, a user may believe that the region highlighted in the screenshot below represents a CNV, but it has not been automatically detected by the software.




Figure 183: A region, not called by the software as a CNV, that the user wants to manually define as a CNV

In order to manually create the CNV call, the user defines the CNV region by using the mouse to select the region in the chromosome ruler track. Alternatively, the coordinates can be provided in the text box in the menu bar of IGV.




Figure 184: Using the ruler region in IGV to select the region to display

The IGV window resets to the size of the required region



Figure 185: IGV set to the boundaries of the region to be defined as CNV

Select the Add CNV option in the Actions menu.




Figure 186: Selection of the Add CNV option from the Actions menu in the variant table header

Users have the option to define the entire region in IGV as the CNV or the software will snap to the nearest probe at each end.




Figure 187: The option to snap to the nearest probe or to keep the region as that displayed in IGV

! Once a CNV has been created manually it is **NOT** automatically displayed in the software. The creation of the CNV can be confirmed by the number of CNVs detected that feed into the protocol filter; in the figure below the number has incremented by one to five.




Figure 188: The number of CNVs has been incremented by one

However, in order to display it in the variant table additional steps need to be taken. In the original analysis the protocol and filters did not select the manually defined region as being a CNV and so in order to include it the default protocol CNV filter need to be updated. This can be done in the Admin Controls > Analysis > Protocols section.




Figure 189: The default CNV filter in the default analysis protocol

Users need to create a filter that allows manually created CNVs to be included and this can be added to the default CNV filter. This is shown in the figure below.




Figure 190: The CNV filter in the Default Protocol has been edited to include OR Source is Manual

Repeating the analysis with the updated filter will result in any manual CNVs being added to the sample variant list.

Merging CNV calls

There are occasions when CNVs are called with small regions in between that the user would like to combine into a single larger CNV.

In order to do this, adjust the scaling in the IGV window such that both CNVs are visible, then right clicking between them in the track will generate a popup menu with the option to Merge Displayed provided.




Figure 191: Selecting the option to merge displayed CNVs

If selected the software will request confirmation of the merge option.




Figure 192: Confirmation of merge option

Following confirmation of the merge option the variant table will be updated. There will be a single row containing the new merged CNV that spans the two previously separate calls. Additionally, the variant counts above the table will be updated.

<div><div><div>All (5)</div><div>CNVs (3)</div><div>LOH Calls (2)</div></div><div>Page 1 of 1 (1 - 3 of 3)</div><div>Page Size: 20</div><div>Actions</div></div>													
<input type="checkbox"/>	Sample	Chromosome	Start	End	Type	Length	Copy Number	# Markers	Mean	Confidence	Overlap	Genome Build	Base
<input type="checkbox"/>	Sample-1	X	100611029	100611266	Duplication	237b	3	3	0.632296	High	CDS (target)	GRCh37	q22.1
<input type="checkbox"/>	Sample-1	13	50566440	51782998	Duplication	1.22Mb	3	70	0.203889	High	CDS (target)	GRCh37	q14.2-q1
<input type="checkbox"/>	Sample-1	13	49384108	49711010	Duplication	326.9Kb	2	21	0.0431704			GRCh37	q14.2

All (4)

CNVs (2)

LOH Calls (2)

Page 1 of 1 (2 - 2 of 2)

Page Size: 20

Actions

<input type="checkbox"/>	Sample	Chromosome	Start	End	Type	Length	Copy Number	# Markers	Mean	Confidence	Overlap	Genome Build	Base
<input type="checkbox"/>	Sample-1	X	100611029	100611266	Duplication	237b	3	3	0.632296	High	CDS (target)	GRCh37	q22.1
<input type="checkbox"/>	Sample-1	13	49384108	51782998	Duplication	2.4Mb	2	121	0.0870438	Low	CDS (target)	GRCh37	q14.2-q1

Figure 193: An updated variant table showing the merging of 2 CNVs to a single row in the table as well as the decrease in the number of All variants and CNV variants

Likewise, in the IGV window the two calls are now combined.



Figure 194: Following the merger of 2 CNVs a single CNV is now displayed

Separating Merged CNV calls

Having been created, users are able to dissolve a merged CNV. Right clicking on the CNV row in the variant table will display the standard popup menu but now with an additional option of Dissolve which will split the merged CNV back into the original separate calls.




Figure 195: Selecting the dissolve option in CNV variant table

Aneuploidy Plots

Interpret is able to provide aneuploidy plots in order for the user to assess whether there is a difference in chromosome number in a set of patients.

This functionality is accessed through the Tools sub-menu in the software dashboard as shown below.




Figure 196: The create aneuploidy plot option

The aneuploidy plot option can only be used when the user is viewing CNV data. If this is not the case then the following error message will be displayed. As with all error messages in the software it can be removed by clicking on it.

Please view the CNV data for which you like to generate an Aneuploidy plot first

Figure 197: Error message from aneuploidy plot option

The figure below illustrates the correct view from which to launch an aneuploidy plot.

 The screenshot shows the 'CNVs + LOH Calls' view. At the top, there's a navigation bar with 'Tools' selected, showing 'Create Aneuploidy Plot'. Below the navigation bar, there's a table of genomic data. The table has columns: Sample, Chromosome, Start, End, Type, Length, Copy Number, # Markers, Mean, Confidence, Overlap, Depth, Frequency, and Inheritance. The table is filtered to show 'CNVs (47)' out of 361 total items. The 'Create Aneuploidy Plot' button is visible above the table.

Sample	Chromosome	Start	End	Type	Length	Copy Number	# Markers	Mean	Confidence	Overlap	Depth	Frequency	Inheritance
10847	X	154722218	154722368	Deletion	150b	0	3	-6.96929	High	CDS (target)	1		Not Tested
12878	X	154722218	154722368	Deletion	150b	0	3	-7.24167	High	CDS (target)	1		Not Tested
11516	1	152278342	152279231	Deletion	889b	1	7	-0.793907	Moderate	CDS (target)	403		Not Tested
11516	15	22772721	28708642	Deletion	5.94Mb	1	379	-1.00102	High	CDS (target)	165	0.34	Not Tested
11516	16	89853091	89866067	Deletion	12.98Kb	1	5	-0.830938	High	CDS (other)	175		Not Tested
12878	19	56531693	56531843	Deletion	150b	1	3	-0.934804	High	CDS (target)	220		Not Tested
12878	22	22395009	23011626	Deletion	616.62Kb	1	12	-0.854273	High	CDS (other)	152		Not Tested
12878	22	23055395	23197546	Deletion	142.15Kb	0	6	-4.45401	High	CDS (other)	12		Not Tested
14130	5	78928	33885969	Deletion	33.81Mb	1	688	-0.993025	High	CDS (target)	175	0.64	Not Tested
14130	14	106321887	106484831	Deletion	162.94Kb	0	4	-6.87598	High	CDS (other)	2		Not Tested
5881	16	2152051	2152267	Deletion	216b	1	3	-0.957263	High	CDS (target)	167		Not Tested
5881	19	43353946	43529051	Deletion	175.1Kb	1	6	-0.817712	High	CDS (other)	123		Not Tested

Figure 198: Selection of aneuploidy plot from menu bar Tools

The user will then be asked to provide the region list that needs to be evaluated for plotting. Information on creating a region list is documented in this manual in the Administration Controls > Analysis > Region lists.




Figure 199: Select a region list pop-up

The available regions lists will be displayed in a drop-down menu; in this example there is a single region list created in the software.




Figure 200: Selection of a region from the drop-down list

Following selection of the required region list the user clicks on the View button.




Figure 201: Selection of the region Chromosome_1

An example aneuploidy plot is displayed below.




Figure 202: An example of an aneuploidy plot for the chosen region across the samples in the batch

Viewing Translocation Events

The variant table has a column selector icon  allowing user to configure which columns are displayed. The figure below shows the columns available for display.




Figure 203: Columns available to select for display in the translocations variant page




Figure 204: Example of a translocation

Translocation Options

As with the page displaying SNV and Indel calls there are options available for each translocation variant called by the software,

Right clicking on a variant will provide a menu of the possible options.




Figure 205: Translocation options

Adding to Shortlist




Figure 206: Adding a translocation to the shortlist

Once a variant has been added to a shortlist the available option is updated to now allow that variant to be deleted from the shortlist.




Figure 207: Selecting to delete a variant from the shortlist

Variant Classification

A variant can be classified from the list that is included by default. These are:

- Benign
- Uncertain significance, likely benign
- Uncertain significance
- Uncertain significance, likely pathogenic
- Pathogenic

Additional classifications can be added in the Admin Controls section of the software (Admin Controls > Analysis > Classifications)




Figure 208: Classify the translocation

A variant classification may change over time and it is possible to track the changes and view the classification history.




Figure 209: Viewing a variant's classification history

Initially, the classification will be blank.




Figure 210: A variant with no classification history

When a classification is made the history table will show the classification type, who made it and when it was made.





Figure 211: Example of a benign classification

Any updates to the classification will be recorded with previous designations retained.



Classification History

Classification history for 22:23290555 > chr9:130731760

Classification	User	Date
Unclassified	admin	21-Apr-2020 16:10:41
Uncertain significance	admin	21-Apr-2020 16:10:08
Benign	admin	21-Apr-2020 16:08:30

Figure 212: Example of a tracking a translocation classification change

Display with Flanking

Users can select to view translocations with flanking sequence




Figure 213: Selecting to show a translocation with flanking sequence

Variant Links

Links to external data sources are available; these are managed in Admin Controls > Analysis > Manage Links




Figure 214: Linking out to external data sources

Variant Table Options

Column Sorting

Rows in the variant table can be sorted using the column header. In the example below the results have been sorted by decreasing and increasing allele frequency. Currently, data can only be sorted by one column.

Alt	Allele Frequency	Type
G	100%	SNV
C	99.57%	SNV
A	48.9%	SNV
T	43.24%	SNV
C	51.07%	SNV
C	100%	SNV
CT	86.16%	Insertion
T	40.99%	SNV
C	47.55%	SNV
A	51.5%	SNV

Alt	Allele Frequency	Type
A	100%	SNV
T	100%	SNV
G	100%	SNV
A	100%	SNV
A	100%	SNV
T	100%	Deletion
G	100%	SNV
G	100%	SNV
T	100%	SNV
G	100%	SNV

Alt	Allele Frequency	Type
A	21.98%	SNV
G	22.95%	SNV
G	24%	Deletion
C	24.27%	Deletion
G	24.69%	Deletion
A	26.3%	SNV
C	26.8%	Deletion
AAAACA	27.27%	Complex
G	27.86%	Deletion
G	28.4%	Deletion

Figure 215: Sorting by Allele Frequency

Dynamic Filtering

As shown previously the variants page displays the Protocol Filter, the number of variants detected by the pipeline and presented to the filter is depicted in a red box and the number remaining in a green box.

In the image below you can see that there are 2946 variants (in the red box) detected by the pipeline that are to be filtered based on the settings in the protocol. Subsequently there are 2754 remaining (as shown in the green box).



Figure 216: The filter used by the protocol in the analysis of the sample displayed in the Variants page

However, the user is able to implement additional filtering dynamically. Any column header with the funnel icon  can be used as a filter. For example, a user may want to filter on Gene Symbol

	HGVSc (Gene Symbol)	Chromosome	Start	End	Ref	Alt	Allele Frequency	Type	Ref Depth	Alt Depth	Total Depth	Ref Quality	Alt Quality	Sample	
---	---------------------	------------	-------	-----	-----	-----	------------------	------	-----------	-----------	-------------	-------------	-------------	--------	---

Figure 217: Selection of the funnel icon for the Gene Symbol column

In this case they want to see only variants found when the total depth is greater than 200.




Figure 218: Dynamic filtering of the variants using the Total Depth column

After updating the Variants view now shows a Dynamic Filter window and within it is the "Total Depth > 200" filter.

From the 2754 variants generated by the protocol using the default filter, it can be seen that a further 265 have been removed filtered with 2489 remaining.




Figure 219: Example of Variants filtered by the protocol filter and a dynamic filter

Dynamic filters can be chained together so additional filters can be added for instance an Allele Frequency greater than 80%




Figure 220: Selection of another dynamic filter to start creating combinations

Now the Dynamic Filter shows "Total Depth > 200 and Allele Frequency > 0.8" and there are now 1457 variants remaining from the input of 2754.





Figure 221: Variants being filtered by a compound dynamic filter

There is no requirement for the user to have to repeat the setting of dynamic filters every time they use the software, there is the option to name the filter and pressing



to retain for re-use.

Alternatively, all dynamic filters can be removed from the display by selecting to clear the filter



Viewing a Sample in IGV

Selection of a variant in the Variants Table causes it to be displayed in the embedded IGV.




Figure 222: A variant selected and the aligned displayed in IGV

Within the IGV window there are several options for modifying the data being displayed.




Figure 223: Display options for IGV

By default, the sequence viewer is centered upon the selected variant but users can drag the display upstream and downstream of the variant position. Also, it is possible to zoom in and out via the magnification slider at the top of the window.

Additionally, the tracks displayed can also be modified via the setting options available on right hand side of the viewer .

For example:




Figure 224: Display options for the integrated IGV browser, firstly using a left click on the mouse and secondly using a right click




Figure 225: Display options available for sample reads, firstly with a left click and secondly with a right click

Using Tracks

Users can add or remove data tracks to the IGV view. This can be from publicly available sources or from proprietary internal or subscription-based sources.

Tracks can be added in the Software section of the Admin Controls (Admin Controls > Software > Annotation) and documentation of how to do this is in this section of the user guide.

To use this functionality, users need to access the Tools tab of IGV




Figure 226: The Tools tab for adding data tracks to annotate an alignment displayed in IGV

Once accessed selecting the drop-down arrow will list the available tracks.




Figure 227: The drop-down list of data tracks available

Select the data track to be added.




Figure 228: Selection of a track

And then click on the  icon to add it to the set of tracks for the software to display




Figure 229: Click on the + icon to add the data track to the display

The selected track will be displayed. It can be removed by clicking on the minus icon




Figure 230: The new track is now loaded

Select any further tracks to add to the view




Figure 231: Addition of the required tracks

Finally, close the Tools tab and the data tracks will be displayed.




Figure 232: Display of the selected data tracks following closure of the tools tab

'Popping out' of the IGV display

There is, potentially, a substantial amount of information that can be displayed in the IGV view. To accommodate the information and make it easier for the user it is possible to 'pop out' the IGV view into a new browser tab.

This is accomplished using the button in the display

Page 1 of 1 (1 - 1 of 1) Page Size: 10 Actions									
<input checked="" type="checkbox"/>	Sample	Donor Gene	Recipient Gene	Donor Chromosome	Donor Breakpoint	Recipient Chromosome	Open IGV in a separate window		
<input checked="" type="checkbox"/>	t010	BCR	ABL1	22	23290555	9	130731760	3340	2758

Figure 233: Button to allow display of IGV into a new tab in the browser

Selecting Multiple Samples

As discussed above users can opt to view multiple samples simultaneously by selecting them in the batch view.

For selected, view:
SNVs/Indels
CNVs/LOH Calls

Completed Samples

Sample	View	# SNVs	# CNVs	# LOH	Report	QC
5881	VCF Logs	2,754	8	16		View
6937	VCF Logs	2,695	13	15		View
7408	SNVs CNVs/LOH VCF Logs	2,740	7	12		View
8210	SNVs CNVs/LOH VCF Logs	2,666	10	16		View
10384	SNVs CNVs/LOH VCF Logs	2,650	4	17		View
10847	SNVs CNVs/LOH VCF Logs	2,669	5	13		View
11516	SNVs CNVs/LOH VCF Logs	2,571	7	16		View
12878	SNVs CNVs/LOH VCF Logs	2,627	14	18		View
14130	SNVs CNVs/LOH VCF Logs	2,614	18	14		View

Figure 234: Selecting multiple samples to view in the Variants page

When multiple samples are selected there will be separate tracks for each sample in IGV. This makes it possible to compare the same variant in different samples.




Figure 235: An example of two samples sharing a variant as displayed in the integrated IGV browser

Variant Table Options

There are options within the Variant table accessed using the Actions drop down menu.




Figure 236: Accessing the options in the variant table

Reporting

Results can be exported by clicking on the Generate Report button below the variant table.




Figure 237 shows the IGV interface with a variant table. The table has columns: Sample, HGVS (Gene Symbol), Chromosome, Start, End, Ref, Alt, Allele Frequency, Type, Total Depth, Ref Quality, Alt Quality, and Quality Score. A dropdown menu is open over the 'Actions' header, showing options: Generate Report, View Shortlist, Filter, and With selected. The 'Generate Report' option is highlighted.

Sample	HGVS (Gene Symbol)	Chromosome	Start	End	Ref	Alt	Allele Frequency	Type	Total Depth	Ref Quality	Alt Quality	Quality Score
5881	FGFR3:c.882T>C	4	1801977	1801977	T	C	48.28%	SNV	406	7336.00	6799.00	4792.21
5881	FGFR3:c.1959G>A	4	1806167	1806167	G	A	100%	SNV	317	0.00	10840.00	9706.90
5881	FGFR3:c.*156_*157del	4	1807400	1807402	CGT	C	99.26%	Deletion	314	62.00	15921.00	13518.20
5881	FGFR3:c.*798C>T	4	1808060	1808060	C	T	100%	SNV	358	0.00	12200.00	10870.20
5881	SH3BP2:c.123G>T	4	2820740	2820740	G	T	99.16%	SNV	238	72.00	8128.00	7168.91
5881	SH3BP2:c.300T>C	4	2824673	2824673	T	C	99.42%	SNV	2	346	348	10668.70
5881	SH3BP2:c.*634C>A	4	2834468	2834468	C	A	99.1%	SNV	4	441	445	13622.30
5881	SH3BP2:c.*1067C>T	4	2834901	2834901	C	T	100%	SNV	0	377	377	11798.10
5881	SH3BP2:c.*2202G>C	4	2836036	2836036	G	C	99.71%	SNV	1	349	350	10774.50
5881	SH3BP2:c.*3877G>C	4	2837711	2837711	G	C	99.02%	SNV	3	304	307	9255.40

Figure 237: Selecting the Generate Report option from the variant table header menu

Interpret provides multiple types of report and for each of these types there are templates. These are highly customisable and updates can be easily applied in Admin Controls-Analysis-Reports.

When selecting to generate a report, the initial window allows users to select the type of report to generate.




Figure 238 shows the initial report option window. It has a 'Report Type' dropdown menu, a checkbox for 'Just report shortlisted variants', and a 'Columns' section with radio buttons for 'Include all columns' (selected) and 'Only include displayed columns'. A green 'Generate Report' button is at the bottom.

Figure 238: Initial report option

Default reports supplied with the software are listed.




Figure 239 shows the selection of PDF report type window. The 'Report Type' dropdown menu is open, showing options: PDF Report with Template, VCF Report, and Report with Template. A hand cursor is pointing at 'PDF Report with Template'. The 'Generate Report' button is at the bottom.

Figure 239: Selection of PDF report type

Once the report type has been selected the user needs to specify the template to use.




Figure 240: Selection of the template to use with the PDF report type

Once all options are chosen, pressing Generate Report will create the PDF file and the web browser will download it.

Sample	Chromosome	Start	End	Length	Genome Build	Ref	Alt	Type
5881	4	1801977	1801977	0b	GRCh38	T	C	SNV
5881	4	1806167	1806167	0b	GRCh38	G	A	SNV
5881	4	1807400	1807402	2b	GRCh38	CGT	C	Deletion
5881	4	1808060	1808060	0b	GRCh38	C	T	SNV
5881	4	2820740	2820740	0b	GRCh38	G	T	SNV
5881	4	2824673	2824673	0b	GRCh38	T	C	SNV
5881	4	2834468	2834468	0b	GRCh38	C	A	SNV

Figure 241: An example of the PDF report generated

Other options are included, for instance an HTML based report.




Figure 242: Selecting a template type report

Again, a template needs to be chosen




Figure 243: Selection of a HTML format report

The HTML formatted report is then generated and available,

Sample	Chromosome	Start	End	Length	Genome Build	Ref	Alt	Type	Genomic Context	Contig Length	HGVSc	HGVSp	HGVSc (Gene Symbol)	Classification	Genotype	Zygosity	Inheritance	Total Depth	Ref Depth	Alt Depth	Allele Frequency
5881	4	1801977	1801977	0b	GRCh38	T	C	SNV	NaN		ENST00000340107.8:c.882T>C	ENSP00000339824.4:p.Asn294=	FGFR3:c.882T>C	Unclassified	0/1	Heterozygous	Not Tested	406	210	196	48.28%
5881	4	1806167	1806167	0b	GRCh38	G	A	SNV	NaN		ENST00000340107.8:c.1959G>A	ENSP00000339824.4:p.Trp553=	FGFR3:c.1959G>A	Unclassified	1/1	Homozygous	Not Tested	317	0	317	100%
5881	4	1807400	1807402	2b	GRCh38	CGT	C	Deletion	Simple_repeat	36	ENST00000340107.8:c.*156_*157del		FGFR3:c.*156_*157del	Unclassified	1/1	Homozygous	Not Tested	314	2	268	89.26%
5881	4	1808060	1808060	0b	GRCh38	C	T	SNV	NaN		ENST00000340107.8:c.*798C>T		FGFR3:c.*798C>T	Unclassified	1/1	Homozygous	Not Tested	358	0	358	100%
5881	4	2820740	2820740	0b	GRCh38	G	T	SNV	NaN		NM_001126891.2:c.123G>T	NP_001116153.1:p.Leu41=	SH3BP2:c.123G>T	Unclassified	1/1	Homozygous	Not Tested	238	2	236	99.16%
5881	4	2824873	2824873	0b	GRCh38	T	C	SNV	NaN		NM_001126891.2:c.300T>C	NP_001116153.1:p.His100=	SH3BP2:c.300T>C	Unclassified	1/1	Homozygous	Not Tested	348	2	346	99.42%
5881	4	2834468	2834468	0b	GRCh38	C	A	SNV	NaN		NM_001126891.2:c.*634C>A		SH3BP2:c.*634C>A	Unclassified	1/1	Homozygous	Not Tested	445	4	441	99.1%

Figure 244: An example of the HTML report

Actions




Figure 245: Options available for configuring the view in IGV

Display Flanking

Users can choose whether or not to display flanking sequence in the IGV display.




Figure 246: Selecting the Display flanking option in the Actions menu

Manage Tracks

Users can add or remove data tracks to the IGV view. This can be from publicly available sources or from proprietary internal or subscription-based sources.

Tracks can be added in the Software section of the Admin Controls (Admin Controls > Software > Annotation) and documentation of how to do this is in this section of the user guide.




Figure 247: Selecting the manage tracks options

The available tracks will be displayed in a pop-up window and users can select the tracks that they want to add to the display.




Figure 248: Tracks available to display

Once the required tracks are selected, users can press Update to update the IGV display.




Figure 249: Selecting tracks to add to the IGV display




Figure 250: Displaying of tracks in the IGV display

Viewing Analysis Results By Variant

As results of samples are generated, they are stored in the Interpret database and can be analysed from a variant-centric point of view

Accessing of this viewpoint is via the Variants button on the dashboard menu bar shown in the figure below.




Figure 251: Selection of Variants from the Dashboard menu bar

Selecting the Variants tab in the menu bar opens up a new page to display all the variants recorded in the database.




Figure 252: The start page for viewing variants

There is a substantial amount of information available in the variants page and the different sections are highlighted in the figure below.




Figure 253: The different sections of the variant page

There a number of active regions

- Filter Pane



The filter pane allows for dynamic filtering of the variants. By default the filter is set to All Variants, so all variants are displayed in the variants table, however these can be refined according to your specific requirements.

- **Variant Table**

This displays all variants in the database that meet the filtering requirements of the dynamic filter. By default this is for displaying all variants.

- **Sample Table**

When a row in the variant table is selected all samples that contain the selected variant will be displayed in this table.

- **IGV Pane**

Selection of samples in the


- **Notes**

Users can add notes to variants.

- **Quick Filters**

These are selection of options that allow users to rapidly filter variants on the basis of some general conditions.

Clicking on the  icon on far right of the column headers in the variant table will display all the columns that can be selected for display in the variant table.



HGVSc (Gene Symbol)	HGVSc
Chromosome	Canonical?
Start	rsID
End	Minor Allele Frequency
Ref	Minor Allele
Type	gnomAD - Total
Genome Build	gnomAD - African
Alt	gnomAD - Latino
Context Length	gnomAD - Ashkenazi Jewish
Genomic Context	gnomAD - East Asian
HGVSp	gnomAD - European (Finnish)
HGVSc	gnomAD - European (non-Finnish)
Classification	gnomAD - South Asian
# Samples	gnomAD - Other
% Samples	ClinVar Significance
#	Gene ID
Most Severe Consequence	Gene Symbol
Impact	Transcript ID
Consequence Terms	Exon Number
PolyPhen Prediction	Protein ID
PolyPhen Score	Length
SIFT Prediction	Transcript Resolution Method
SIFT Score	Exon ID

Figure 254: Columns available to select for display in the variant table


Similarly clicking on the same icon on the sample table provides a series of column options.

Sample ID	Ref Strand Bias
Allele Frequency	Alt Strand Bias
Ref Depth	Reads Placed Left
Alt Depth	Reads Placed Right
Total Depth	Sex
Quality Score	Homozygosity
Ref Quality	Read 1
Alt Quality	Read 2
Genotype	Read 1 Size
Zygosity	Read 2 Size
Inheritance	Batch Name
Log Ratio	Batch Date
# Ref Reads (+)	User
# Ref Reads (-)	Protocol
# Alt Reads (+)	Panel
# Alt Reads (-)	

Figure 255: Columns available to select for display in the sample table

Dynamic Filtering


The dynamic filters are the same set of options discussed in the previous section. They can be accessed from the Actions drop down menu shown below.



	HGVSc (Gene Symbol)	Chromosome	Start	End	Ref	Type	Genome Build	Alt	Context Length	Genomic Context	HGVSp	HGVSc	Generate Report	# Samples	% Samples	#	Most Severe	Impact	Consequence	Polyphen Prediction
<input type="checkbox"/>	DVL1c.*627C>T	1	1270895	1270895	G	SNV	GRCh37	A		NaN		XM_005244731.1:c.*62	Filter	1	10%	5132	3_prime_UTR	MODIFIER	3 Prime UTR Va	
<input type="checkbox"/>	DVL1c.*630C>T	1	1271059	1271059	G	SNV	GRCh37	A		NaN		XM_005244731.1:c.*62		1	10%	2922	3_prime_UTR	MODIFIER	3 Prime UTR Va	
<input type="checkbox"/>	DVL1c.*583C>T	1	1271139	1271139	G	SNV	GRCh37	A		NaN		XM_005244731.1:c.*583C>T	Unclassified	1	10%	5133	3_prime_UTR	MODIFIER	3 Prime UTR Va	
<input type="checkbox"/>	DVL1c.*547T>C	1	1271175	1271175	A	SNV	GRCh37	G		NaN		XM_005244731.1:c.*547T>C	Unclassified	10	100%	1	3_prime_UTR	MODIFIER	3 Prime UTR Va	

Figure 256: Accessing the dynamic filtering options

The dynamic filtering window provides a detailed set of options for investigating the variants stored within the Interpret database.



Filter Variants

Create New Filters

Category

- Basic Variant Attributes
- Consequence/Severity Predictions
- Population Frequencies
- Classification/Pathogenicity
- Group Frequencies and Odds

Type

Settings

Current Filter

Filter Name: All Variants

Existing Filters

Include combinations?

Recent

Favourites

All

AND OR

Cancel Update Reset All

Figure 257: Full filtering options available to filter variants

Filtering by Quick Filters

Selecting the Quick Filters tab on the side of the variants page opens a tab that provides some options for quickly drilling down into variants of interest. The options currently available are to select based on a classification type, the NGS panel or gene.




Figure 258: Quick filter options

Classification or Panel can be selected by pressing the corresponding buttons with multiple selections allowed. To filter by genes start typing the gene name in the text box matching values will be displayed.



Gene

SKI

SKI - ENSG00000157933

SKI

SKINTL - ENSG00000242267

SKINTL

SKIL - ENSG00000136603

SKIL

SKIV2L2

SKIV2L2 - ENSG00000039123

SKIV2L

SKIV2L - ENSG00000204351

Figure 259: Using quick filters to select for the gene SKI

Once a gene is selected it will be displayed as below and can be removed by clicking on the x next to the gene name.

Gene

SKI x

Figure 260: Using quick filters to filter by the SKI gene

The dynamic filter is now updated and shows that, from the input of 8432 variants, there are only 10 found within the SKI gene.

Dynamic Filters: In SKI

8432

In SKI

10

Page 1 of 1 (1 - 10 of 10)

Page Size: 30

Actions

Figure 261: Displaying only variants in the SKI gene

Additional genes can be selected and these will be displayed in the same way.

Gene

SKI x DVL1 x

Figure 262: Using quick filters to select variants in the SKI or DVL1 genes

When the 2-gene filter is applied the output now increases to 19 variants being displayed.



Each row of the variant table represents a variant that has been detected in at least one sample. Selecting a variant displays all the samples in which the variant is present. In the figure below, the variant DVL1:c.*347T>C is present in 10 samples.



Subsequently, selecting any of the sample or samples rows will display the alignment for the variant in the corresponding sample.





Adding Notes to Variants

It is possible for users to add annotations to variants through the notes function. When a variant has been selected and there are rows populated in the sample table, the user can make a right click on one of these. From the popup menu select the Notes > Add Note options




Figure 266: Selecting the Add Note option

A window is displayed with a text box where up to 250 characters can be used. Any other pre-existing notes will also be shown.




Figure 267: Note creation template

The user can enter the required notation.

+Add Note To Variant:

Notes for DVL1:c.*347T>C

Characters remaining:194

This is an important variant that needs to be monitored.

Existing Notes (0)

Cancel + Add Note

Figure 268: Note creation

Then selecting + Add Note completes the process and the existing notes section is updated to include the newly created note.

+Add Note To Variant:

Notes for DVL1:c.*347T>C

Characters remaining:250

Existing Notes (1)

Admin User 2021-01-19 21:25:00.0

This is an important variant that needs to be monitored.

Cancel + Add Note

Figure 269: Note generation

Once all changes have been made, the notes window can be closed and the view will return the normal variant display with the note now being displayed in the Notes panel.

DVL1:c.*347T>C															Notes (1)	
Sample ID	Allele Frequency	Ref Depth	Alt Depth	Total Depth	Quality Score	Ref Quality	Alt Quality	Genotype	Zygosity	Inheritance	Log Ratio	# Ref Reads (+)	# Ref Reads (-)	# Alt Reads (+)	# Alt Reads (-)	Ref Strand Bias
10847	100%	0	118	118	3,508.4199	0	4.038	1/1	Homozygous	Not Tested	86.3965	0	0	87	31	?
11516	99.40%	1	183	184	5,470.4902	36	6.324	1/1	Homozygous	Not Tested	125.598	1	0	129	34	100%
12878	100%	0	176	176	5.265	0	6.046	1/1	Homozygous	Not Tested	126.599	0	0	123	53	?
14130	100%	0	174	174	5,252.23	0	6.024	1/1	Homozygous	Not Tested	125.213	0	0	121	53	?
4315	100%	0	180	180	5,966.6899	0	6.217	1/1	Homozygous	Not Tested	125.372	0	0	113	67	?
5681	100%	0	138	138	4,122.5898	0	4.771	1/1	Homozygous	Not Tested	100.259	0	0	100	38	?
6937	100%	0	212	212	6,410.6201	0	7.330	1/1	Homozygous	Not Tested	151.552	0	0	145	67	?
7408	100%	0	176	176	5,277.2598	0	6.051	1/1	Homozygous	Not Tested	126.599	0	0	123	53	?
10384	100%	0	130	130	3,867.28	0	4.443	1/1	Homozygous	Not Tested	94.7143	0	0	90	40	?
8210	100%	0	141	141	4,273.6299	0	4.872	1/1	Homozygous	Not Tested	102.339	0	0	84	57	?

Figure 270: Displaying a note in the note panel

A note can be deleted by pressing the red rubbish bin icon. If the Confirm Delete option is then selected the note will be removed.




Figure 271: Deleting a note

Users can also edit a note by clicking on the pen icon; which will show the note in a text box where changes can be made. The update is confirmed by the pressing the Apply button.




Figure 272: Update an existing note

The notes panel is subsequently updated to show the revised note.

DVL1:c.*347T>C																	Notes (1)	
Page 1 of 1 (1 of 10) of 10 Page Size: 10 Actions Filter																		
Sample ID	Allele Frequency	Ref Depth	Alt Depth	Total Depth	Quality Score	Ref Quality	Alt Quality	Genotype	Zigosity	Inheritance	Log Ratio	# Ref Reads (+)	# Ref Reads (-)	# Alt Reads (+)	# Alt Reads (-)	Ref Strand Bias		
10847	100%	0	118	118	3,508,419.9	0	4,038	1/1	Homozygous	Not Tested	86,396.5	0	0	87	31	?		
11516	99.46%	1	183	184	5,470,490.2	36	6,324	1/1	Homozygous	Not Tested	125,938	1	0	129	54	100%		
12878	100%	0	176	176	5,265	0	6,046	1/1	Homozygous	Not Tested	126,599	0	0	123	53	?		
14130	100%	0	174	174	5,262.23	0	6,024	1/1	Homozygous	Not Tested	126,215	0	0	121	53	?		
4315	100%	0	180	180	5,496,689.9	0	6,217	1/1	Homozygous	Not Tested	129,372	0	0	113	67	?		
5881	100%	0	138	138	4,122,589.8	0	4,771	1/1	Homozygous	Not Tested	100,259	0	0	100	38	?		
6937	100%	0	212	212	6,416,620.1	0	7,350	1/1	Homozygous	Not Tested	151,552	0	0	145	67	?		
7408	100%	0	176	176	5,277,295.8	0	6,051	1/1	Homozygous	Not Tested	126,599	0	0	123	53	?		
10384	100%	0	130	130	3,867.28	0	4,443	1/1	Homozygous	Not Tested	94,714.3	0	0	90	40	?		
8210	100%	0	141	141	4,275,629.9	0	4,872	1/1	Homozygous	Not Tested	102,339	0	0	84	57	?		

This is an important variant that needs to be monitored.
Update: this variant is now classified as pathogenic.
by admin @ 19 Jan 2021 21:25:00

Figure 273: Display of an updated note

Administration Controls

These are accessed by selecting the Admin Controls options in the user drop down menu.




Figure 274: Accessing the Admin Controls section

There are 4 main parts to the admin controls:

1. Overview
2. User Controls
 - Current Users
 - Add Users
3. Analysis
 - Manage Samples
 - Current Analyses
 - Protocols
 - Panels
 - Region Lists
 - Variant Lists
 - Classifications
 - Metric Sets
 - Manage Links
 - Filters
 - Preferred Transcripts
 - Reports

- Guidelines
- #### 4. Software
- Software Overview
 - Annotation
 - Advanced Settings
 - Plug-ins

Administration Overview

The administration overview provides a view of the latest activity in Interpret

Admin Controls			
Overview			
Overview			
User Controls	Date	Action	User
Analysis	08-Jan-2020 12:17:35	Started analysis of CytoSure NGS Panel - 298678 - 8 Jan 2020 12:17:35	admin
Software	08-Jan-2020 11:56:21	Started analysis of deleted batch	admin

Figure 275: The overview window showing the latest activity in Interpret

User Controls

Current User

Current Users shows a list of all the current users.


 Home Batches Samples Order Help & Support Tools v3.2.5				
Admin Controls				
Overview				
User Controls				
User Controls	Username	First Name	Last Name	Roles
Current Users	admin	Admin	User	Viewer, Administrator
Add Users				
Analysis				
Software				

Figure 276: Front page for current users in user controls



Selecting a user provides an overview of the user account

Home

Batches

Order

Help & Support

Filters

User Account - admin

Username

admin

Firstname

Surname

Last Logged In

11 Jan 2018 15:36:16

Roles

☒ Administrator

☐ Analyst

☐ Director

☐ Technician

☒ Viewer

View Activity

Modify Permissions

Change Password

Figure 277: An overview of a user account

Selecting the "View Activity" button will give an overview of what the user has done and when.

User Account - admin

Activity		
Date	Action	User
08-Jan-2020 12:17:35	Started analysis of CytoSure NGS Panel - 298678 - 8 Jan 2020 12:17:35	admin
08-Jan-2020 11:56:21	Started analysis of deleted batch	admin

Figure 278: Viewing a user's activity

Selecting "Modify Permissions" will display the current permissions available and those with a tick in the adjacent checkbox have been enabled for the user.

As the current user is admin, as expected, all permissions are enabled.




Figure 279: Permission options that can be modified

Finally, selecting "Change Password" gives a pop-up window that allows the user to change their password.

Figure 280: Pop-up menu to change a user's password

Add User

Add User provides a means to add new users. The following information needs to be supplied or selected:

- User details
- Login details
- Roles

Following correct completion of the form the 'Create User' button becomes active allowing the process of adding a user to be completed.

Admin Controls

Overview | **Add new user**

User Controls

Current Users | Add Users | Analysis | Software

User details

First Name: *

Last Name: *

Login details

Username: *

Password: *

Confirm Password: *

Roles

☐ Administrator

☐ Analyst

☐ Director

☐ Technician

☒ Viewer

[Create User](#)

Figure 281: The add user start page

When the details are syntactically correct there will be a tick in the corner of the window.

User details

First Name: *

Last Name: *

Field data validated

Figure 282: In-form validation of the new user details

Similarly, there are controls in place for the login details. A username that already exists cannot be used.

User details

First Name: *

Last Name: *

Login details

Username: *

Password: *

Confirm Password: *

Figure 283: In-form validation of the new user login details

An existing user name cannot be used.

Login details

Username: *

Password: *

Confirm Password: *

Figure 284: Example of using a pre-existing username

The form also checks to ensure the password is entered correctly.

Login details

Username: *

Password: *

Confirm Password: *

Figure 285: Example of entering non-matching passwords

Lastly, roles are assigned from the list using the checkboxes.

Roles

☒ Administrator

☐ Analyst

☒ Director

☐ Technician

☒ Viewer

[Create User](#)

Figure 286: Selection of the roles for the new user

There are currently 5 different roles defined within the software and these are described in the following table.



Role	Abilities
Administrator	All
Director	All
Analyst	Starting an analysis and viewing results
Technician	Loading FASTQs and Starting an analysis
Viewer	Basic permissions required for all users

Table 4: Roles available to assign

In order to login to Interpret, all users are required to be a Viewer. Other roles can be assigned as required.

Selecting "Create User" processes the form and there is a popup to confirm that the new user has been created.




Figure 287: Popup menu confirming creating of the new user

Now looking at the Current Users in the User Controls you can see that the new user is listed.

User Controls

Username	First Name	Last Name	Roles
admin	Admin	User	Viewer, Administrator
peters	Peter	Smith	Viewer, Administrator, Director

New user

Figure 288: Display of the newly created user in the user table

Analysis

Manage Samples

Interpret allows users to manage samples and the data associated with them.



Figure 289: The Manage Samples start page

There are two tabbed panes displayed, an overview and a variables table.




Figure 290: The two sections of the manage samples page

Selecting a sample in the overview tab brings up a series of sub-pages.
There is a batch history showing when the sample had been used in an analysis.




Figure 291: The batch history for a sample

There is a table showing all data associated with a sample name, including duplicate data if the sample data has been uploaded more than once.

These data can be deleted by selecting the Delete Data Files option.




Figure 292: The associated raw data for a sample

If a user tries to delete data files, there will be a popup menu asking for confirmation prior to any data being deleted.




Figure 293: Popup menu asking for deletion confirmation

The user can see the variant and QC results associated with a sample.




Figure 294: The results generated for a sample

Lastly, there are actions available for a sample. Currently, this is limited to updating a sample name.




Figure 295: Actions available for a sample

Entering a new name and selecting Update Sample Name will update the sample name.




Figure 296: Updating a sample name

In the variables tab, users can modify existing variables or create new ones.

Manage Samples

OverviewVariables

Samples Variables

Sample	Edit
5881	Edit
6937	Edit
7408	Edit
8210	Edit
10384	Edit
10847	Edit
11516	Edit
12878	Edit
14130	Edit
4315	Edit
Test Sample	Edit

Manage Variables

Figure 297: Sample variables

Selecting Edit allows the user to add variables for a sample or a new variable category.

Manage Samples

OverviewVariables

Manage Sample Variables

+ Add New Variable

+ Add New Variable Category

Save Changes

Discard Changes

return

Figure 298: Managing variables for a sample

Current Analyses

This provides a means to view all current analyses.

Admin Controls

Overview

User Controls

Analysis

Manage Samples

Current Analyses

Protocols

Panels

Region Lists

Variant Lists

Classifications

Metric Sets

Manage Links

Filters

Preferred Transcripts

Reports

Guidelines

Software

Analyses

	Name	Panel	Protocol	Number of Samples	Status	Date	User
	CytoSure NGS Batch 0002	CytoSure NGS Panel	Default Protocol	3	Completed	10 Jan 2020 11:38:27	admin
	CytoSure NGS Batch 0001(1)	CytoSure NGS Panel	Default Protocol	10	Completed	8 Jan 2020 14:51:53	admin

Figure 299: The Current Analyses start page



As with other pages in the software where there is a column selector icon

Column Selector

Analyses							
<input type="checkbox"/>	Name	Panel	Protocol	Number of Samples	Status	Date	User
<input type="checkbox"/>	CytoSure NGS Batch 0002	CytoSure NGS Panel	Default Protocol	3	Completed	10 Jan 2020 11:38:27	admin
<input type="checkbox"/>	CytoSure NGS Batch 0001	CytoSure NGS Panel	Default Protocol	10	Completed	8 Jan 2020 14:51:53	admin

Figure 300: Current analyses with the column selector highlighted

Columns can be added or removed as required from the popup menu.

Analyses							
<input type="checkbox"/>	Name	Panel	Protocol	Number of Samples	Status	Date	User
<input type="checkbox"/>	CytoSure NGS Batch 0002	CytoSure NGS Panel	Default Protocol	3	Completed	10 Jan 2020 11:38:27	admin
<input type="checkbox"/>	CytoSure NGS Batch 0001(1)	CytoSure NGS Panel	Default Protocol	10	Completed	8 Jan 2020 14:51:53	admin

Select Displayed Columns

- ☒ Name
- ☒ Panel
- ☒ Protocol
- ☒ Number of Samples
- ☒ Status
- ☒ Date
- ☒ User
- ☐ Samples

Figure 301: Columns available to select for display


Additionally columns can be sorted using column filters.

Column Filters

<input type="checkbox"/>	Name	Panel	Protocol	Number of Samples	Status	Date	User
<input type="checkbox"/>	CytoSure NGS Batch 0002	CytoSure NGS Panel	Default Protocol	3	Completed	10 Jan 2020 11:38:27	admin
<input type="checkbox"/>	CytoSure NGS Batch 0001	CytoSure NGS Panel	Default Protocol	10	Completed	8 Jan 2020 14:51:53	admin

Figure 302: Column filtering options available

For example, in order to view all analysis with the Default Protocol this can be selected from the drop down in the protocol column




Name	Panel	Protocol	Number of Samples	Status	Date	User
CytoSure NGS Batch 0002	CytoSure NGS Panel	Default Protocol	3	Completed	10 Jan 2020 11:38:27	admin
CytoSure NGS Batch 0001	CytoSure NGS Panel	Mosaic Protocol	10	Completed	8 Jan 2020 14:51:53	admin

Figure 303: Selection of analyses only performed using the Default Protocol

Protocols

A Protocol defines how a sample is analysed.



Admin Controls

Overview

User Controls

Analysis

Manage Samples

Current Analyses

Protocols

Panels

Region Lists

Variant Lists

Classifications

Metric Sets

Manage Links

Filters

Preferred Transcripts

Reports

Guidelines

Software

Protocols

Default Protocol

Mosaic Protocol

Name Default Protocol

Version 1

Metric Set Default

Created By admin

Created On 8 Jan 2020

Protocol SNV Filter Default SNV Filter [Show](#)

Protocol CNV Filter Default CNV and LOH Filter [Show](#)

Basic Pipeline Configuration **Advanced Pipeline Configuration** **Germline Configuration**

Pipeline Version CytoSure NGS v0.1.9_patch1

Enable CNV and LOH Calling Yes

Import SNV Results Yes

Reference Samples for CNV detection 5881 (5881_7_L001001), 6937 (6937_7_L001001), 7408 (7408_7_L001001), 8210 (8210_7_L001001), 10384 (10384_7_L001001), 10847 (10847_7_L001001), 11516 (11516_7_L001001), 12878 (12878_7_L001001), 14130 (14130_7_L001001)

[+ New Protocol](#) [Edit](#) [Copy](#) [Delete](#)

Figure 304: The Protocols start page

There are 4 components in a Protocol:

- 1. Pipeline Type (and Pipeline Capabilities)**– the analysis pipeline type with which the protocol is compatible, and the functions of that pipeline that should be run as part of the analysis of samples in batches which use the protocol (its "capabilities"). When a batch is created, Interpret matches pipeline type and



capabilities with those supported by the selected panel.

Pipeline Type

Somatic

Pipeline Capabilities

- ☒ CNV Calling
- ☒ Translocation Calling
- ☒ PTD Calling
- ☒ ITD Calling
- ☒ UMI Analysis

2. **Metric Set** – settings with which to qualitatively assess the run data.
3. **Protocol Filter(s)** – a filter with which to process all variants produced by the analysis pipeline. Currently there are filters for the following:
 - a. SNV filter
 - b. CNV & LOH filter
 - c. Translocation filter
4. **Pipeline Configurations** – specific configurable settings used in the pipeline.
 - a. Basic pipeline configuration
 - b. Advanced pipeline configuration
 - c. Germline configuration – if germline mode is selected in the basic configuration settings
 - d. Somatic configuration – if somatic mode is selected in the basic configuration settings

The screenshot displays the configuration interface for an analysis protocol. At the top, a header bar contains the text "Name Default Protocol", "Version 1", "Metric Set Default", "Created By admin", and "Created On 7 Jan 2020". Below this, the "Protocol SNV Filter: Default SNV Filter" section is expanded, showing a list of filters: "Alt Strand Bias < 80%", "Reads Placed Left > 4", "Reads Placed Right > 4", and "Quality Score > 1.0". The "Protocol CNV Filter: Default CNV and LOH Filter" section is also expanded, showing a complex logical expression: "NOT Type is LOH" OR "Mean >= 0.4" OR "Mean <= -0.6" OR "Copy Number = 0" OR "Type is LOH" OR "Mean Standard Error <= 0.12" OR "Mean <= -2" OR "# Markers >= 3". Below the filters, there are three tabs: "Basic Pipeline Configuration", "Advanced Pipeline Configuration", and "Germline Configuration". The "Basic Pipeline Configuration" tab is selected, showing settings for "Pipeline Version" (CytoSure NGS v0.1.9_patch1), "Enable CNV and LOH Calling" (Yes), "Import SNV Results" (Yes), and "CNV Calling Mode" (Germline). At the bottom, there are three buttons: "Edit", "Copy", and "Delete".

Metric Set

SNV Filter

CNV Filter

Pipeline Configuration

Figure 305: The sections of an analysis protocol

At the bottom of the protocol are a set of tabbed pages with different configuration settings.

There are basic pipeline settings:



<div>Basic Pipeline Configuration Advanced Pipeline Configuration Germline Configuration</div>	
Pipeline Version	CytoSure NGS v0.2.1
Enable CNV and LOH Calling	Yes
Import SNV Results	Yes
CNV Calling Mode	Germline
Enable Translocation Detection	Yes
Reference Samples for CNV detection	All batch samples

Figure 306: Default basic pipeline configurations

There are advanced pipeline settings:



Basic Pipeline Configuration		Advanced Pipeline Configuration	Germline Configuration
Flanking Region (bp)	0		
Allele Balance Priors Off	Yes		
Minimum Alt Count	5		
Minimum Base Quality	20		
Minimum Total Read Count	20		
Minimum Mapping Quality	30		
Max Read Mismatch Fraction	0.04		
Threshold for Targets Not Covered	1		
Min Depth for CNV segment BAF	20		
Drop Outliers Threshold	10		
Use CBS with Smoothing	Yes		
# CPUs for Processing	12		
RAM (GB)	2		
Call PTDs	Yes		
Enable Enhanced Structural Variant Calling	No		
Minimum Reads for Translocation	3		
Mosaicism Correction Factor (Duplications)	0.1		
Mosaicism Correction Factor (Deletions)	0		

Figure 307: Default advanced pipeline configurations

Depending on the CNV calling mode selected in the basic pipeline configuration there will be either a Germline Configuration tab



Basic Pipeline Configuration	Advanced Pipeline Configuration	Germline Configuration
Germline Minimum Alt Fraction	0.2	
Germline Pooled Continuously	No	
Germline Pooled Discretely	No	
Germline Segmentation p-value Threshold	0.01	
Germline Copy Number Thresholds	-2.0,-0.25,0.2,0.7	

Figure 308: Germline configuration settings

Or a Somatic Configuration tab.

Basic Pipeline Configuration	Advanced Pipeline Configuration	Somatic Configuration
Somatic Minimum Alt Fraction	0.01	
Somatic Pooled Continuously	Yes	
Somatic Pooled Discretely	Yes	
Somatic Segmentation p-value Threshold	0.05	
Somatic Copy Number Thresholds	-2.0,-0.07,0.07,0.7	
Tumour Content Estimation	Yes	

Figure 309: Somatic configuration settings

Modifying a Protocol

To modify a protocol users select the Edit button at the bottom of the protocol. A new version of the protocol is displayed which is the same, with the exception of the Version number that is incremented.

The user can retain the same metric set and filters or select alternatives from the drop-down lists.



Name

Default Protocol

Version

2

Metric Set

Default

▼

Protocol SNV Filter

Default SNV Filter [16]

▼

▼ Protocol SNV Filter: Default SNV Filter

Hide

Alt Strand Bias < 80%

+

Reads Placed Left > 4

+

Reads Placed Right > 4

+

Quality Score > 1.0

Protocol CNV Filter

Default CNV and LOH Filter [30]

▼

▼ Protocol CNV Filter: Default CNV and LOH Filter

Hide

NOT

Type is LOH

+

Mean >= 0.4

OR

Mean <= -0.6

+

Mean Standard Error <= 0.12

OR

Copy Number = 0

+

Markers >= 3

OR

Mean <= -2

Type is LOH

OR

Mean >= 0.4

+

Mean Standard Error <= 0.12

+

Markers >= 3

Figure 310: Editing an existing protocol

Protocol configurations can also be set; with the default values being updated. Software update does not automatically change the pipeline version for the protocol. User must modify protocol to modify desired pipeline version.

Basic Pipeline Configuration

Advanced Pipeline Configuration

Germline Configuration

Pipeline Version

CytoSure NGS v0.1.9_patch1

▼

Enable CNV and LOH Calling

☒ Yes ☐ No

Import SNV Results

☒ Yes ☐ No

Reference Samples for CNV detection

▼

CNV Calling Mode

☒ Germline ☐ Somatic

Figure 311: Setting basic pipeline configurations

Similarly with advanced pipeline configurations.



Basic Pipeline Configuration **Advanced Pipeline Configuration** Germline Configuration

Flanking Region (bp)

0

Allele Balance Priors Off

☒ Yes ☐ No

Minimum Alt Count

5

Minimum Alt Fraction

0.2

Minimum Base Quality

20

Minimum Total Read Count

20

Minimum Mapping Quality

30

Pooled Continuously

☐ Yes ☒ No

Pooled Discretely

☐ Yes ☒ No

Max Read Mismatch Fraction

0.04

Threshold for Targets Not Covered

1

Min Depth for CNV segment BAF

20

Segmentation p-value Threshold

0.01

Drop Outliers Threshold

10

Use CBS with Smoothing

☒ Yes ☐ No

Copy Number Thresholds

-2.0,-0.25,0.2,0.7

CPUs for Processing

12

RAM (GB)

2

Call PTDs

☒ Yes ☐ No

Mosaicism Correction Factor (Duplications)


0.1

Mosaicism Correction Factor (Deletions)

0

Figure 312: Setting advanced pipeline configurations

And lastly with the germline configurations.



Basic Pipeline Configuration Advanced Pipeline Configuration **Germline Configuration**

Germline Minimum Alt Fraction

Germline Pooled Continuously ☐ Yes ☒ No

Germline Pooled Discretely ☐ Yes ☒ No

Germline Segmentation p-value Threshold

Germline Copy Number Thresholds

Copy Number	Log Ratio
0	<input type="text" value="-2"/>
1	<input type="text" value="-0.25"/>
3	<input type="text" value="0.2"/>
4	<input type="text" value="0.7"/>

Figure 313: Setting germline configurations

Creating a Reference Pool for calling CNVs


In order to call CNVs, the software will use a set of samples to curate a pool of references against which variants can be detected.

Samples that will be part of the reference pool to be used for calling CNVs can be selected from all the samples that have previously been loaded in the software.



It is important that samples used for the reference have been processed using the NGS panel with baits from the same lot that has been used for processing the test samples.

OGT can provide a set of suitable reference sample data for each of the bait lots that can be purchased. These can be used by the user as a starting point for building a reference of control samples suitable for CNV calling.



Basic Pipeline Configuration Advanced Pipeline Configuration

Pipeline Version

Enable CNV and LOH Calling

Import SNV Results

Reference Samples for CNV detection

- 5881 (5881_7_L001001)
- 6937 (6937_7_L001001)
- 7408 (7408_7_L001001)
- 8210 (8210_7_L001001)
- 10384 (10384_7_L001001)
- 10847 (10847_7_L001001)
- 11516 (11516_7_L001001)
- 12878 (12878_7_L001001)
- 14130 (14130_7_L001001)

1-9/10

Figure 314: Selecting samples for the reference pool

Once selected the sample will be displayed, highlighted in blue to denote pending status.

Reference Samples for CNV detection

5881 (5881_7_L001001) ×

Figure 315: A sample selected for the reference pool

Samples can be added as required.

Reference Samples for CNV detection

5881 (5881_7_L001001) ×

6937 (6937_7_L001001) ×

7408 (7408_7_L001001) ×

8210 (8210_7_L001001) ×

10384 (10384_7_L001001) ×

10847 (10847_7_L001001) ×

11516 (11516_7_L001001) ×

12878 (12878_7_L001001) ×

14130 (14130_7_L001001) ×

Figure 316: A reference pool of samples for a protocol

Once all required samples have been selected for the pool, the configuration can be saved.

Basic Pipeline Configuration

Advanced Pipeline Configuration

Germline Configuration

Pipeline Version

CytoSure NGS v0.1.9_patch1

Enable CNV and LOH Calling

Yes

Import SNV Results

Yes

Reference Samples for CNV detection

5881 (5881_7_L001001), 6937 (6937_7_L001001), 7408 (7408_7_L001001), 8210 (8210_7_L001001), 10384 (10384_7_L001001), 10847 (10847_7_L001001), 11516 (11516_7_L001001), 12878 (12878_7_L001001), 14130 (14130_7_L001001)


CNV Calling Mode

Germline

Figure 317: The basic pipeline configuration showing a a reference pool

Basic Pipeline Configuration

Setting	Description	Default Setting	Other Options
Pipeline Version	Select the version of the pipeline to use in the protocol.		
Enable CNV and LOH Calling	Enable or disable calling of CNVs and LOH by the protocol.	Yes	No
Import SNV Results	Enable importing of the SNV calls into the database.	Yes	No
CNV Calling Mode (CytoSure NGS (< 0.2.13) only)	Select the calling mode to either germline or somatic variants.	Germline	Somatic
Enable Translocation Detection	Enable or disable detection of translocations.	Yes	No




Setting	Description	Default Setting	Other Options
Reference Samples for CNV Detection	Specify to use all samples in the batch, select specific samples from all samples in the system, or indicate that the user should select one or more of the samples in the batch when the batch is started.	All batch samples	Specific samples, Samples from batch
Call PTDs	Enable calling of partial tandem duplications	Yes	No
Call ITDs	Enable calling of internal tandem duplications	Yes	No
Run UMI Processing (Somatic only)	Enable analysis with unique molecular identifiers (UMIs)	Yes	No
Hotspots (Somatic only)	Indicate which variants should be specifically analysed by the pipeline for monitoring purposes	(empty)	


Table 5: Basic Pipeline Configuration Settings

Advanced Pipeline Configuration

Category	Setting	Description	Default Setting	Other Options
SNV Detection	Flanking Region (bp)	The amount of flanking sequence to include in the analysis.	0	Any value from 0 to 60
SNV Detection	Allele Balance Priors Off	Disable use of aggregate probability of observation balance between alleles as a component of the priors.	Yes	No
SNV Detection	Minimum Alt Count	The minimum alternative allele read count	5	
SNV Detection	Minimum Alt Fraction	The minimum alternative allele read fraction	0.2	
SNV Detection	Minimum Base Quality	The minimum allowed base quality	20	
SNV Detection	Minimum Total Read Count	The minimum total read count	20	
SNV Detection	Minimum Mapping Quality	The minimum mapping quality of the reads	30	



Category	Setting	Description	Default Setting	Other Options
SNV Detection	Pooled Continuously	Output all alleles which pass input filters, regardless of genotyping outcome or model.	No	Yes
SNV Detection	Pooled Discretely	Assume that samples result from pooled sequencing and model pooled samples using discrete genotypes across pools.	No	Yes
SNV Detection	Max Read Mismatch Fraction	The maximum fraction of mismatches in the read	0.04	
Quality Control	Threshold for Targets Not Covered	The coverage threshold to consider a target as not covered	1	
CNV Detection	Min Depth for CNV Segment BAF	The minimum depth required calculating the CNV segment B-allele frequency	20	
CNV Detection	Segmentation p-value Threshold	The significance threshold to accept segment	0.01	
CNV Detection	Drop Outliers Threshold	Drop bins that lie more than this many multiples of the 95th quantile away from the average within a rolling window as they are considered as outliers.	10	
CNV Detection	Use CBS with Smoothing	Use smoothing with the circular binary segmentation algorithm.	Yes	No
CNV Detection	Copy Number Thresholds	Thresholds for calling each integer copy number, i.e. 0,1,3,4	-2.0, -0.25, 0.2, 0.7	
Hardware Resources	# CPUs for Processing	The number of CPUs to be used by the protocol	12	
Hardware Resources	RAM(GB)	The amount of memory to be used by the protocol	2	
Structural Variant Detection	Threshold for PTDs	Exon-to-control ratio threshold to indicate a PTD	1.9	
Structural Variant Detection	Minimum Reads for Translocation	The minimum number of reads required to call translocations	3	
CNV Detection	Mosaicism Correction Factor (Duplications)	The mosaicism correction factor for duplications	0.1	



Category	Setting	Description	Default Setting	Other Options
CNV Detection	Mosaicism Correction Factor (Deletions)	The mosaicism correction factor for deletions	0	
UMI Processing	Minimum Input Base Quality	Minimum base quality to be used to call molecular consensus reads	10	
UMI Processing	Minimum Read Mapping Quality	Minimum read mapping quality to be used in reads grouping by UMI	30	
UMI Processing	Minimum Reads Supporting Consensus	Minimum number of reads supporting a consensus base/read. Used in consensus reads filtering. 0 indicates no filtering.	0	
Fusion Detection	Supporting Threshold	Minimum number of reads in support of the fusion call	5	

Table 6: Advanced Pipeline Configuration Settings

Germline Configuration

Setting	Description	Default Setting	Other Options
Germline Minimum Alt Fraction	The minimum alternate allele fraction	0.2	
Germline Pooled Continuously	Output all alleles which pass input filters, regardless of genotyping outcome or model.	No	Yes
Germline Pooled Discretely	Assume that samples result from pooled sequencing and model pooled samples using discrete genotypes across pools.	No	Yes
Germline Segmentation p-value Threshold	The significance threshold to accept segment	0.01	

Setting	Description	Default Setting		Other Options
Germline Copy Number Thresholds	Set the log ratio thresholds for each of the copy number variations.	Cop y Nu mb er	Lo g Rat io	
		0	-2	
		1	-0.25	
		3	0.2	
		4	0.7	

Table 7: Germline Configuration Settings

Somatic Configuration

Setting	Description	Default Setting	Other Options
Somatic Minimum Alt Fraction	The minimum alternate allele fraction	0.01	
Somatic Pooled Continuously	Output all alleles which pass input filters, regardless of genotyping outcome or model.	Yes	No
Somatic Pooled Discretely	Assume that samples result from pooled sequencing and model pooled samples using discrete genotypes across pools.	Yes	No
Somatic Segmentation p-value Threshold	The significance threshold to accept segment	0.05	



Setting	Description	Default Setting		Other Options
Somatic Copy Number Thresholds	Set the log ratio thresholds for each of the copy number variations.	Cop y	Lo g	
		Nu mb er	Rat io	
		0	-2	
		1	-0.07	
		3	0.07	
		4	0.7	
Tumour Content Estimation	Enable the software to estimate the percentage tumour content in the sample	Yes		No

Table 8: Somatic Configuration Settings

Panels

The **Panels** page displays information about the panels currently loaded into the software.

Home

Batches

Samples

Variants

Help & Support

Tools

GRCh37v3.5.21

Admin Controls

Overview

User Controls

Analysis

Manage Samples

Current Analyses

Protocols

Panels

Region Lists

Variant Lists

Classifications

Metric Sets

Manage Links

Filters

Preferred Transcripts

Reports

Guidelines

Software

Name	Product Code	Deprecated	Types	Capabilities	
CytoSure NGS Panel	502003	false	[CytoSure NGS (< 0.2.13)]	[ITD Calling, Translocation Calling, PTD Calling, CNV Calling]	
SureSeq CLL + CNV Panel	602022	false	[CytoSure NGS (< 0.2.13)]	[ITD Calling, Translocation Calling, PTD Calling, CNV Calling]	
SureSeq Ovarian Cancer Panel	600073	false	[CytoSure NGS (< 0.2.13)]	[ITD Calling, Translocation Calling, PTD Calling, CNV Calling]	
SureSeq Myeloid Panel	600075	false	[CytoSure NGS (< 0.2.13)]	[ITD Calling, Translocation Calling, PTD Calling, CNV Calling]	
SureSeq Comprehensive FH Panel	601004	false	[CytoSure NGS (< 0.2.13)]	[ITD Calling, Translocation Calling, PTD Calling, CNV Calling]	
SureSeq Core MPN Panel	602001	false	[CytoSure NGS (< 0.2.13)]	[ITD Calling, Translocation Calling, PTD Calling, CNV Calling]	

+ Add Panel

Figure 318: The panel start page

Adding New Panels

New panels can be added via the **Add Panel** button.




Figure 319: The Add Panel window

To add a new panel to the system, users must indicate with which of the available **Pipeline Types**, and which of the **Optional Capabilities** of those pipeline types, the panel is compatible (contact OGT support if unsure). Assign a unique **Name** and **Product Code**, select the correct **Genome Build** for the **Targets File**, then select the file from the file system by clicking the **Select Targets File** button, and finally click **Add Panel**.



Targets File Format

Targets Files should be provided by OGT. Standard BED files are not compatible with Interpret and their use will result in pipeline failure. Contact OGT for the correct Targets File for your panel and build.

The **Add Panel** window may also be used to add a targets file to an existing panel for different genome build – select the **Existing Panel** radio button and the select the appropriate panel from the **Existing Panels** drop-down list in order to do this.

Modifying Existing Panels

Attributes of existing panels may be modified by double-clicking on the appropriate row of the table, making the required changes to the **Name**, **Product Code**, **Types** and/or **Capabilities**, and clicking **Save**.

Name	Product Code	Deprecated	Types	Capabilities	
CytoSure NGS Panel	502003	false	<input checked="" type="checkbox"/> CytoSure NGS (< 0.2.13) <input checked="" type="checkbox"/> CytoSure NGS (0.2.13+) <input type="checkbox"/> Somatic	<input type="checkbox"/> UMI Analysis <input checked="" type="checkbox"/> CNV Calling <input checked="" type="checkbox"/> PTD Calling <input checked="" type="checkbox"/> ITD Calling <input checked="" type="checkbox"/> Translocation Calling	Save Cancel

Figure 320: Editing attributes of a panel



Region Lists

Region lists are a set of defined genomic regions

Admin Controls

Overview

User Controls

Analysis

Manage Samples

Current Analyses

Protocols

Panels

Region Lists

Variant Lists

Classifications

Metric Sets

Manage Links

Filters

Preferred Transcripts

Reports

Guidelines

Software

Name	# Regions
------	-----------

+ Add Region List

Figure 321: The Region Lists start page

To add a region list users select the 'Add Region List' button which provides a form.

Add Region List

Name *

Add Regions

Import Regions from File

Choose file

No file chosen

Import

Add Genomic Feature

Type

Add Region Details

Genome Build

Chr

Start

End

Add

GRCh38

0

0

Added Regions *

Name	Chr	Start	End	Genome Build
------	-----	-------	-----	--------------

Save

Cancel

Figure 322: Form for adding a Region List

Regions can be defined using three parameters:

- 1. Regions can be imported from a file.
- 2. Regions can be defined as genomic features
- 3. Regions can be defined as chr:start-finish

The first step is defining the name of the region list

Figure 323: Setting the region list name

Importing Regions from a File

 Importing a file of defined regions requires the file be in BED or zipped BED formats.

Figure 324: The file chooser for importing regions from a file

Selecting 'Choose File' opens a file browser and when a file is selected it is shown in the window, as below.

Figure 325: Selection of a file to be imported

When the Import button is pressed the selected file is checked and if there is a problem with the file format the following is displayed.

Figure 326: File import error message

Clicking on the error will remove it.
If the file format is correct then the user needs to select the genome build from the drop-down menu.




Figure 327: Selecting the genome build

Finally, selecting Import completes the process




Figure 328: Select import to complete the upload

And the regions are now displayed in the table.

Import Regions from File

Choose file

Regions Test File.bed

Import

Add Genomic Feature

Type

Add Region Details

Genome Build

GRCh38

Chr

Start

0

End

0

Add

Added Regions *

	Name	Chr	Start	End	Genome Build
<input type="checkbox"/>	17:46094556-46094705	17	46,094,556	46,094,705	GRCh38
<input type="checkbox"/>	17:46170856-46225432	17	46,170,856	46,225,432	GRCh38
<input type="checkbox"/>	6:180970-335276	6	180,970	335,276	GRCh38
<input type="checkbox"/>	16:86568165-86568664	16	86,568,165	86,568,664	GRCh38

Figure 329: Imported regions displayed



The upload of regions can be edited by selecting the  icon on the side to remove the region on that row.

Defining a Region as a Genomic Feature

Regions can be added as genomic features selected from the drop-down menus. Currently available features are exons, genes, proteins and transcripts.

Add Genomic Feature

Type

Chr Start End Add

Chr 0 0

Chr	Start	End
-----	-------	-----

Figure 330: Types of genomic feature available

When a type has been chosen a second menu is provided containing the features of the chosen type that can be selected.

Users can scroll up and down the list to find the required feature

Add Genomic Feature

Type

Select Gene

Add Region Details

Genome Build

GRCh38

Added Regions *

Name

DDX11L1 - ENSG00000223972

WASH7P - ENSG00000227232

MIR6859-1 - ENSG00000278267

MIR1302-2HG - ENSG00000243485

MIR1302-2 - ENSG00000284332

FAM138A - ENSG00000237613

OR4G4P - ENSG00000268020

OR4G11P - ENSG00000240361

OR4F5 - ENSG00000186092

1-9/117501

End 0 Add

Name	Start	End
------	-------	-----

Figure 331: List of genes available to add as a region

Alternatively, typing text in the text field will automatically subset those features that contain the text.

In the figure below the user has typed DMT and the menu lists all features in which this text is found.

Add Genomic Feature

Type

Gene

Add Region Details

Genome Build

GRCh38

Added Regions *

	Name	Start	End
--	------	-------	-----

Select Gene

DMT

TRDMT1 - ENSG00000107614

DMTF1 - ENSG00000135164

DMTN - ENSG00000158856

TRDMT1 - ENSG00000107614

DMTF1 - ENSG00000135164

DMTN - ENSG00000158856

End

0

Add

Figure 332: Filtering genes available by gene symbol

Pressing return selects the feature and adds it to the region list.

Add Genomic Feature

Type

Gene

Add Region Details

Genome Build

GRCh38

Chr

Start

0

End

0

Add

Select Gene

DMTF1 - ENSG00000135164

Added Regions *

	Name	Chr	Start	End	Genome Build
-	DMTF1	7	87,152,361	87,196,337	GRCh38

Figure 333: Addition of a genomic feature to the region list

When selecting features, the drop-down lists are populated with the following information:

- Exon – the EnSEMBL ENSE number
- Gene – the gene name
- Protein – the EnSEMBL ENSP number
- Transcript – the EnSEMBL ENST number

Added Regions *					
	Name	Chr	Start	End	Genome Build
-	ENSE00002234944	1	11,869	12,227	GRCh38.94
-	DDX11L1	1	11,869	14,409	GRCh38.94
-	ENSP00000493376	1	65,565	65,573	GRCh38.94
-	ENST00000461467	1	35,245	36,073	GRCh38.94

Figure 334: A region list showing different genomic features

Regions can be defined as chromosome, start and end

The final mechanism by which to create a region is to manually specify the build, chromosome, start and end positions.

Add Region Details

Genome Build Chr Start End Add

 0 0

Figure 335: The add region detail section of the form

The drop-down menu provides a list if available genome builds.

Add Region Details

Genome Build Chr Start End Add

 0 0

Ad Chr Start

Figure 336: Specification of the required genome build

Subsequently, the chromosome can be selected from he populated list.

Add Region Details

Genome Build Chr Start End Add

 0 0

Added Regions * Name Chr Start

Figure 337: Selection of the chromosome from the drop-down menu

Lastly, the user can enter the start and end coordinates.

Add Region Details

Genome Build Chr Start End Add

 1,000,900 2,000,000

Figure 338: Selection of region details

Pressing Add appends the region to the table

Add Region Details

Genome Build

Chr

Start

0

End

0

Add

Added Regions *

	Name	Chr	Start	End	Genome Build
<div></div>	3:1000900-2000000	3	1,000,900	2,000,000	GRCh38

Figure 339: Addition of the region details to the new region list

Pressing save completes the process and you now see that there is a new Region List called Demo Region List that comprises 4 defined regions.

Admin Controls

Overview

User Controls

Analysis

Manage Samples

Current Analyses

Protocols

Panels

Region Lists

Variant Lists

Classifications

Metric Sets

Manage Links

Filters

Preferred Transcripts

Reports

Guidelines

Software

Name	# Regions
Demo Region List	4

+ Add Region List

Figure 340: Addition of the new region list

Using a Region List

When a Region List has been defined it can be used to filter variants.

Region Lists can be found in the 'Region/Variant List' category and then the 'Region List' type.

All available Region Lists will be presented in the drop down menu.

Variant Filters

CNV Filters

Create New Filters

Category

Prevalence

Population Frequencies

Classification/Pathogenicity

Genes, Exons and Proteins

Region/Variant Lists

Type

Region List(s)

Variant List(s)

Settings

Demo Region List

Figure 341: Availability of the new region list in the Filters

When it has been selected a Region List is added to the filter builder.




Figure 342: Addition of the Demo Region List to the filter

Alternatively, the newly created Region List can be used to dynamically filter variants.

In the figure below the Demo Region List has been deployed dynamically with the result that 3449 SNPs are filtered down to just 1 SNV.




Figure 343: Deployment of a region list in a dynamic filter

Variant Lists

Interpret allows users to generate variant lists that can be used as a filter in analysis of data generated for samples.




Figure 344: The Variant List start page

Selecting Add Variant List provides a popup where the user can enter the name of the new variant list.




Figure 345: Creating a new variant list

Once a variant list has been created, variants can be added from files. The format for uploading variants is either VCF or HGVS formats.




Figure 346: Uploading variants to a variant list

Once created and populated with variants the new variants list is available as filter.




Figure 347: Selecting the new variant list

The variant list can also be appended to within the results page. Right clicking on a variant provides an option to add a variant to an existing variant list.

	Sample	HGVSc (Gene Symbol)	Chromosome	Start	End	Ref	Alt	Minor Allele Frequency	Allele Frequency	Type	Ref Depth	Alt Depth	Total Depth
	5881	DVL1:c.*347T>C	1	1335795	1335795	A	G	0.05	100%	SNV	0	138	138
	5881	DVL1:c.366A>G	1	1342				0.00	99.57%	SNV	1	232	233
	5881	ATAD3A:c.671G>A	1	1520	DVL1:c.*347T>C			0.01	48.9%	SNV	185	177	362
	5881	ATAD3A:c.*154C>T	1	1534				0.42	43.24%	SNV	21	16	37
	5881	ATAD3A:c.*330G>C	1	1534					51.07%	SNV	183	191	374
	5881	SKI:c.*1802T>C	1	2308					100%	SNV	0	470	470
	5881	SKI:c.*3182dup	1	2309					86.16%	Insertion	22	137	166
	5881	CAMTA1:c.288C>T	1	6887					40.99%	SNV	190	132	322
	5881	CAMTA1:c.*419T>C	1	6888				0.3	47.55%	SNV	203	184	387
	5881	CAMTA1:c.1350G>A	1	7663				0.40	51.5%	SNV	194	206	400
	5881	CAMTA1:c.*890T>A	1	7767				0.19	51.69%	SNV	157	168	325
	5881	CAMTA1:c.*1364_*1365dup	1	7767843	7767843	C	CAA		94.44%	Insertion	9	153	224

Figure 348: Adding a variant to the new variant list



Classifications

It can be helpful to assign colours to particular classifications and in this page users are able to create new and modify existing assignments.

By default the software will ship with 5 classifications already created. These are listed below and displayed in Figure XX.

- 1. Benign
- 2. Uncertain significance (likely benign)
- 3. Uncertain significance
- 4. Uncertain significance (likely pathogenic)
- 5. Pathogenic




Figure 349: View of the default settings in Classifications

Adding a Classification



Figure 350: Setting a new classification name

By default the colour to be assigned is black but users can select a different colour by pressing on the Colour button.

This will produce a sub-window containing 3 tabs of different colour palettes called RGB, HSV and Swatches, each of which is displayed below.

Selecting on a name will display the tab or that palette.




Figure 351: Colour palettes available to setting a classification colour

Once the colour has been selected, you can click on OK to add the classification to the software



Figure 352: A new classification ready to add

And the new classification is displayed, as shown below:

Manage Classification Types

Name *

Colour

Save Clear

Name	Colour	
Benign	#008000	
Pathogenic	#9f0000	
Uncertain significance: likely benign	#a58f2f	
Extra Classification	#d324ff	
Uncertain significance: likely pathogenic	#d9642e	
Uncertain significance	#f79646	

Figure 353: New classification added

A classification can be removed by clicking on the waste bin icon on the row of the classification that needs to be removed.

	Sample	HGVSc (Gene Symbol)	Chromosome	Start	End	Ref	Alt	Allele Frequency	Type	Ref Depth	Alt Depth	Total Depth	Ref Quality	Alt Quality	Quality Score
<input checked="" type="checkbox"/>	5881	DVL1:c.*347T>C	1	1325795	1325795	A	G	100%	SNV	0	138	138	0.00	4771.00	4122.5
<input type="checkbox"/>	5881	DVL1:c.366A>G	1	1342153	1342153	T	C	99.57%	SNV	1	232	233	32.00	7624.00	6750.4
<input type="checkbox"/>	5881	ATAD3A:c.*154C>T	1	1520153	1520153	G	A	48.9%	SNV	185	177	362	6194.00	6050.00	4259.0
<input type="checkbox"/>	5881	ATAD3A:c.*330G>C	1	1534402	1534402	G	C			21	16	37	738.00	515.00	275.1
<input type="checkbox"/>	5881	SKI:c.*1802T>C	1	2308567	2308567	T	C			183	191	374	6338.00	6565.00	4665.6
<input type="checkbox"/>	5881	SKI:c.*3182dup	1	2309937	2309937	C	CT			0	470	470	0.00	15743.00	14005.0
<input type="checkbox"/>	5881	CAMTA1:c.288C>T	1	6887707	6887707	C	T			22	137	166	752.00	8035.00	6276.5
<input type="checkbox"/>	5881	CAMTA1:c.*419T>C	1	6888144	6888144	T	C						6523.00	4644.00	3142.3
<input type="checkbox"/>	5881	CAMTA1:c.1350G>A	1	7663897	7663897	G	A						7148.00	6375.00	4498.8
<input type="checkbox"/>	5881	CAMTA1:c.*890T>A	1	7767381	7767381	T	A						6530.00	7083.00	5107.8
<input type="checkbox"/>	5881	CAMTA1:c.*1364_*1365del	1	7767843	7767843	C	CAA	94.44%	Insertion				5559.00	5895.00	4270.7
<input type="checkbox"/>	5881	CAMTA1:c.*1640del	1	7768113	7768114	TA	T	70%	Deletion				291.00	8624.00	7109.4
<input type="checkbox"/>	5881	CAMTA1:c.*1738del	1	7768218	7768219	CA	C	92.83%	Deletion				856.00	3439.00	2613.2
<input type="checkbox"/>	5881	CAMTA1:c.*2522del	1	7769003	7769004	AT	A	55.33%	Deletion				537.00	12239.00	9947.0
<input type="checkbox"/>	5881	CAMTA1:c.*2556G>A	1	7769047	7769047	G	A	50.38%	SNV				5160.00	11149.00	8450.7
<input type="checkbox"/>	5881	CAMTA1:c.*2735A>G	1	7769226	7769226	A	G	45.71%	SNV				6783.00	7084.00	5113.0
<input type="checkbox"/>	5881	RERC:c.*2212G>A	1	8352875	8352875	C	T	51.22%	SNV	180	189	369	6045.00	4987.00	3481.8
													6259.00	6475.00	4658.6

Figure 354: The new classification can now be used to apply to a variant

The variants table is then updated to show the colour of the new classification.

	Sample	HGVSc (Gene Symbol)	Chromosome	Start	End	Ref	Alt	Allele Frequency	Type	Ref Depth	Alt Depth	Total Depth	Ref Quality	Alt Quality	Quality Score
<input checked="" type="checkbox"/>	5881	DVL1:c.*347T>C	1	1325795	1325795	A	G	100%	SNV	0	138	138	0.00	4771.00	4122.5
<input type="checkbox"/>	5881	DVL1:c.366A>G	1	1342153	1342153	T	C	99.57%	SNV	1	232	233	32.00	7624.00	6750.4
<input type="checkbox"/>	5881	ATAD3A:c.671G>A	1	1520153	1520153	G	A	48.9%	SNV	185	177	362	6194.00	6050.00	4259.0

Figure 355: A variant annotated with the newly generated classification type

Metric Sets

The metric sets define parameters to qualitatively assess an analysis. Selecting the Metric Sets page will display all the metric sets that have been created in Interpret.




Figure 356: The metric sets start page

The settings in the metric set are used in the batch reports for samples that have been analysed. Each metric reported is coloured red, blue or green depending on the values set as Excellent, Pass or Failed accordingly.

Completed Samples											
<input type="checkbox"/>	Sam...	View	# SNVs	# CNVs	# LOH	Mean Target Coverage	% Reads Aligned	% Duplication	Targets Not Cove...	% Usable On Target Reads	Report
<input type="checkbox"/>	5881	<input checked="" type="checkbox"/> SNVs <input type="checkbox"/> CNVs/LOH	2,754	8	16	344.21	99.55	9.08	138	71.95	
<input type="checkbox"/>	6937	<input checked="" type="checkbox"/> SNVs <input type="checkbox"/> CNVs/LOH	2,695	13	15	556.51	99.63	10.75	122	72.18	
<input type="checkbox"/>	7408	<input checked="" type="checkbox"/> SNVs <input type="checkbox"/> CNVs/LOH	2,740	7	12	381.59	99.64	8.6	142	72.41	
<input type="checkbox"/>	8210	<input checked="" type="checkbox"/> SNVs <input type="checkbox"/> CNVs/LOH	2,666	10	16	356.72	99.58	8.48	137	73.48	
<input type="checkbox"/>	10384	<input checked="" type="checkbox"/> SNVs <input type="checkbox"/> CNVs/LOH	2,650	4	17	313.08	99.63	8.05	150	73.06	
<input type="checkbox"/>	10847	<input checked="" type="checkbox"/> SNVs <input type="checkbox"/> CNVs/LOH	2,669	5	13	314.99	99.57	8.53	140	72.82	
<input type="checkbox"/>	11516	<input checked="" type="checkbox"/> SNVs <input type="checkbox"/> CNVs/LOH	2,571	7	16	407.31	99.63	9.29	136	72.07	
<input type="checkbox"/>	12878	<input checked="" type="checkbox"/> SNVs <input type="checkbox"/> CNVs/LOH	2,627	14	18	388.99	99.61	8.47	129	72.9	
<input type="checkbox"/>	14130	<input checked="" type="checkbox"/> SNVs <input type="checkbox"/> CNVs/LOH	2,614	18	14	396.32	99.61	8.41	150	72.49	

Figure 357: A table of completed samples showing metrics in different colours

For example the default metric set that is included in Interpret has the settings shown in the figure below.

Metric Set: Default

Created by admin

Created on 7 Jan 2020 12:13:19

Thresholds

Name	Excellent	Good	Poor
% Reads Aligned	>95.0	>80.0	<=80.0
% Duplication	<5.0	<20.0	>=20.0
Mean Target Coverage	>1500.0	>1000.0	<=1000.0
Targets Not Covered	<1.0	<10.0	>=10.0
Aligned Reads GC %	>70.0	>40.0	<=40.0
Aligned Reads Per Base Quality	>65.0	>40.0	<=40.0
% Usable On Target Reads	>55.0	>50.0	<=50.0
% Usable On Target Bases	>55.0	>50.0	<=50.0
% Usable On And Near Target Reads	>55.0	>50.0	<=50.0

OK

Figure 358: Settings in the default metric set

Selecting 'Add Metric Set' on the start page opens a window that allows users to set their own definitions.

Name	Excellent	Pass	Failed
% Reads Aligned	> 95	* > 80	* <=80.0
% Duplication	< 5	* < 20	* >=20.0
Mean Target Coverage	> 1,500	* > 1,000	* <=1000.0
Targets Not Covered	< 1	* < 10	* >=10.0
Aligned Reads GC %	> 70	* > 40	* <=40.0
Aligned Reads Per Base Quality	> 65	* > 40	* <=40.0
% Usable On Target Reads	> 55	* > 50	* <=50.0
% Usable On Target Bases	> 55	* > 50	* <=50.0
% Usable On And Near Target Reads	> 55	* > 50	* <=50.0
Uniformity of Coverage	< 0	* < 0	* >=0.0
Coverage Efficiency	< 0	* < 0	* >=0.0
Off Target Reads	< 0	* < 0	* >=0.0
% Reads Mapping Quality 0	< 0	* < 0	* >=0.0
Average Quality	< 0	* < 0	* >=0.0
Average Insert Size	< 0	* < 0	* >=0.0
Insert size std	< 0	* < 0	* >=0.0
Evenness	< 0	* < 0	* >=0.0
Uniformity	< 0	* < 0	* >=0.0
Sample Sex	< 0	* < 0	* >=0.0
# Exon Targets Not Covered	< 0	* < 0	* >=0.0
# SegDup Exon Targets Not Covered	< 0	* < 0	* >=0.0

Add Metric Set **Cancel**

Figure 359: Settings available for metric sets

Once the settings have been made, select 'Add metric Set' to complete and the new metric set will be displayed in the Metric Sets start page.

Name	Username	Created On
Default	admin	7 Jan 2020 12:13:19
New Metric Set	admin	27 Feb 2020 10:57:46

Figure 360: New metric set displayed in start page

Once a metric set has been created it is available to use in a protocol.

Figure 361: New metric set available for use in protocols

Manage links provides a means to define links to external resources such as ClinVar.

Figure 362: The Manage Links start page

New links can be generated by adding a description and URL in the table. Please contact OGT support to assist with this.

Manage Links to External Resources

Description *

Refseq (Gene)

URL *

https://www.ncbi.nlm.nih.gov/gene/?term=\$gene.getName

Type

All

Save

Clear

Figure 363: Defining a new link

A link type can be for only SNVs, only CNVs or both

Manage Links to External Resources

Description *

Refseq (Gene)

URL *

https://www.ncbi.nlm.nih.gov/gene/?term=\$gene.getName

Type

All

All

CNV

SNV

Save

Clear

Figure 364: Setting the link type

Once saved the link new link is displayed in the table on the start page.

Description	Link	Type
Ensembl (Position)	http://ensembl.org/Homo_sapiens/contigview?region=\$variant.getChr();start=\$variant.getStart();end=\$variant.getEnd()	All
Ensembl (Gene)	https://ensembl.org/Homo_sapiens/Gene/Summary?db=core;g=\$gene.getName()	SNV
Ensembl (Transcript)	https://ensembl.org/Homo_sapiens/Transcript/Summary?db=core;t=\$transcript.getName()	SNV
ClinVar (HGVSc)	https://www.ncbi.nlm.nih.gov/clinvar/variant/\$variant.getHgvscsGene()	SNV
ExAC (Gene)	http://exac.broadinstitute.org/variant/\$variant.getName()	SNV
ExAC (Transcript)	http://exac.broadinstitute.org/variant/\$variant.getName()	SNV
ExAC	http://exac.broadinstitute.org/variant/\$variant.getChr().getChrName()-\$variant.getStart()-\$variant.getRef()-\$variant.getAlt()	SNV
Refseq (Gene)	https://www.ncbi.nlm.nih.gov/gene/?term=\$gene.getName	SNV

Figure 365: New link displayed

The new link is now available to be used in the results page for variants.

Page 1 of 92 (1 - 30 of 2754)Page Size: 30Actions

Sample	HGVSc (Gene Symbol)	Chromosome	Start	End	Ref	Alt	Minor Allele Frequency	Allele Frequency	Type	Ref Depth	Alt Depth	Total Depth
5881	ATAD3A:c.*330G>C	1	1534402	1534402	G	C	0.05	51.07%	SNV	183	191	374
5881	SKI:c.*1802T>C	1	2308567	2308567	T	C	0.19	100%	SNV	0	470	470
5881	SKI:c.*3182dup	1	2309937	2309937				86.16%	Insertion	22	137	166
5881	CAMTA1:c.288C>T	1	6887707	6887707	C	T		40.99%	SNV	190	132	322
5881	CAMTA1:c.*419T>C	1	6888144	6888144	T	C	0.37	47.55%	SNV	203	184	387
5881	CAMTA1:c.1350G>A	1	7663897	7663897	G	A	0.40	51.5%	SNV	194	206	400
5881	CAMTA1:c.*890T>A	1	7767381	7767381	T	A	0.19	51.69%	SNV	157	168	325
5881	CAMTA1:c.*1364_*1365dup	1	7767843	7767843				94.44%	Insertion	9	153	224
5881	CAMTA1:c.*1640del	1	7768113	7768113				70%	Deletion	27	63	124
5881	CAMTA1:c.*1738del	1	7768218	7768218				92.83%	Deletion	14	363	239
5881	CAMTA1:c.*2522del	1	7769003	7769003				55.33%	Deletion			346
5881	CAMTA1:c.*2756G>A	1	7769047	7769047	G	A		50.38%	SNV			399
5881	CAMTA1:c.*2721G>A	1	7769226	7769226	A	G		45.71%	SNV			315
5881	REERC:c.*2212G>A	1	8352875	8352875	C	T		51.22%	SNV	189	193	369
5881	REERC:c.*2158T>G	1	8352929	8352929	A	G		100%	SNV	0	331	331
5881	REERC:c.*1830T>G	1	8353257	8353257	A	C		46.54%	SNV	232	202	434
5881	REERC:c.*1308T>C	1	8353779	8353779	A	G		46.26%	SNV	187	161	348
5881	REERC:c.*461del	1	8354625	8354625	AT	A		52.25%	Deletion	159	174	338

Figure 366: New link available for use in results tables



Filters

Filters provide the means to control which variants are displayed in the SNV and CNV variant tables.




Figure 367: The start window for the filters

The filters page is divided into 3 parts

1. Creation of new filters
2. The current filter
3. Existing filters




Figure 368: The different section of the filters page

New filters are created by selecting criteria and combining to generate a chain. Opening the Create New Filters tab shows category, type and settings.




Figure 369: Sections of the Create New Filters component

Within Category are a range of options and these are different for SNVs and CNVs.



Figure 370: SNV filter categories



Figure 371: CNV filter categories

Selecting one will then populate Type with the different values, for example below are some of the types of basic variant categories.



Figure 372: Example of displaying type available for a selected category

Selecting a type will then populate Settings with specific options for that type; for example below users can set a filter on selected chromosomes including mitochondrial (MT).

Figure 373: The setting available for the selected category type

To create a filter, users make the necessary selections and click on the on white plus in the green box to the side of the window.

Figure 374: Selection of settings to use in a filter

This will add the filter to the "Current Filter" window.

Figure 375: Display of the newly created filter

The filter options are numerous and very flexible and software allows users to chain filters together in order to tailor to a users needs.

Figure 376: Specifying a second filter

Adding another filter provides the user with the option how the second filter will be combined with the first filter.

The user must now decide to specify whether the options are to satisfy either filter with an OR statement...




Figure 377: Combining filters with an OR statement

...or to satisfy both filters with an AND statement.




Figure 378: Combining filters with an AND statement

In this case AND was selected and the current filter shows the chain of 2 filters.




Figure 379: A filter showing 2 settings combined with an AND statement

Additional filters can be added, for example adding a filter on Allele Frequency.




Figure 380: Adding an allele frequency filter




Figure 381: Adding the allele frequency filter with an AND statement




Figure 382: The updated filter is now displayed

As more filters are added they are grouped and this is reflected in the display by the drawing of dashed boxes around filter groups. Thus far there is a single filter contained within a single dashed box.

More complicated arrangements are possible whereby individual filters can selected and modified. In this chain of 3 filters the Allele Frequency > 80.0 has been selected and when selected the colour changes from blue to green.




Figure 383: Selection of a part of an existing filter

As such an additional filter, in this case for Allele Frequency < 20.0, will be added to the selected filter




Figure 384: Modifying part of an existing filter

When presented with the choice of AND or OR the new filter is only shown with the previously selected filter.




Figure 385: Updating part of an existing filter

The current filter will be updated to show the OR combination of the two filters.




Figure 386: Addition of the new filter

Clicking on the home or house icon returns the current filter to the full view with the hierarchy of filters defined by dashed boxes.




Figure 387: Updating of the display to show the new filter

As the filter builds the sections of dashed boxes increase.




Figure 388: An example of a filter with multiple sections

The user can modify individual parts or groups of filters by using the mouse. As the mouse moves over the filter the different parts will be highlighted in a different colour.




Figure 389: Highlighting different parts of a filter

If there is a modification to make to a part of a filter, then highlight the area to be updated.




Figure 390: Selection of part of the filter

And clicking on it will show the sub-section chosen.




Figure 391: A section of the filter

This sub-section can be modified




Figure 392: Modification of the selected part of a filter

Clicking on the home icon returns the view to the whole filter.

This can be continued until the required filter has been generated. The flexibility in Interpret allows a user to build highly complex queries.




Figure 393: An example of a complex filter generated in Interpret

Preferred Transcripts

It is usual for genes to have multiple transcripts. As users sometimes have specific transcripts of interest, Interpret provides the means to define a preferred transcript for each gene.




Figure 394: The Preferred Transcript start page

When no preferred transcript has been assigned the default behaviour is for the software to select the longest canonical transcript and report the annotation for that transcript alone.

To see a different annotation the user therefore needs to set their preferred transcript using either the admin controls prior to running an analysis, or from within the variants table post-analysis. The table will update with the different annotation once this is done.

To set a preferred transcript the user needs to enter a search term as in the figure below.




Figure 395: Adding a search term

The text entered in the field is matched across all the genes in the database. The search term must simply be present somewhere in the gene name. All matches will be displayed.



Preferred Transcript		
JAK		
Name	Gene Symbol	Preferred Transcript
ENSG00000162434	JAK1	<input type="text"/>
ENSG00000188385	JAKMIP3	<input type="text"/>
ENSG00000105639	JAK3	<input type="text"/>
ENSG00000152969	JAKMIP1	<input type="text"/>
ENSG00000280780	JAKMIP2-AS1	<input type="text"/>
ENSG00000176049	JAKMIP2	<input type="text"/>
ENSG00000096968	JAK2	<input type="text"/>
JAK1	JAK1	<input type="text"/>
JAKMIP1	JAKMIP1	<input type="text"/>
JAKMIP2-AS1	JAKMIP2-AS1	<input type="text"/>
JAKMIP2	JAKMIP2	<input type="text"/>
JAK2	JAK2	<input type="text"/>
JAKMIP3	JAKMIP3	<input type="text"/>
JAK3	JAK3	<input type="text"/>

Figure 396: Search results

For each gene returned by the search term there will be a drop down menu of transcripts available.

Preferred Transcript		
JAK		
Name	Gene Symbol	Preferred Transcript
ENSG00000162434	JAK1	<input type="text"/>
ENSG00000188385	JAKMIP3	<input type="text"/>
ENSG00000105639	JAK3	<input type="text"/>
ENSG00000152969	JAKMIP1	<input type="text"/>
ENSG00000280780	JAKMIP2-AS1	<input type="text"/>
ENSG00000176049	JAKMIP2	<input type="text"/>
ENSG00000096968	JAK2	<input type="text"/>
JAK1	JAK1	<input type="text"/>
JAKMIP1	JAKMIP1	<input type="text"/>
JAKMIP2-AS1	JAKMIP2-AS1	<input type="text"/>
JAKMIP2	JAKMIP2	<input type="text"/>
JAK2	JAK2	<input type="text"/>
JAKMIP3	JAKMIP3	<input type="text"/>
JAK3	JAK3	<input type="text"/>

- NM_001321856.1
- NM_001320923.1
- NM_001321852.2
- NM_001321854.2
- NM_001321853.2
- NM_001321857.2
- NM_002227.4
- NM_001321855.2

Figure 397: Selecting a transcript form the list available

When this is done the preferred transcripts page will update to reflect the selection.

Preferred Transcript		
JAK		
Name	Gene Symbol	Preferred Transcript
ENSG00000162434	JAK1	<div></div>
ENSG00000188385	JAKMIP3	<div></div>
ENSG00000105639	JAK3	<div></div>
ENSG00000152969	JAKMIP1	<div></div>
ENSG00000280780	JAKMIP2-AS1	<div></div>
ENSG00000176049	JAKMIP2	<div></div>
ENSG00000096968	JAK2	<div></div>
JAK1	JAK1	NM_001321856.1 <div>Change</div>
JAKMIP1	JAKMIP1	<div></div>
JAKMIP2-AS1	JAKMIP2-AS1	<div></div>
JAKMIP2	JAKMIP2	<div></div>
JAK2	JAK2	<div></div>
JAKMIP3	JAKMIP3	<div></div>
JAK3	JAK3	<div></div>

Figure 398: The selected transcript for a gene is displayed

Reports

The reporting section displays the reports and their templates currently loaded in the software; each type of report can have different templates associated with it, OGT is able to assist with the creation of customer specific reports and templates.

Admin Controls				
Overview	Available Reports			
User Controls	Name	Description	Type	
Analysis	Batch Report with template	A basic report for a batch that can be customised using templates	Report for a Batch	<div>Templates</div>
	PDF Report with Template	A basic pdf report that can be customised using templates	Report for Variants in Sample	<div>Templates</div>
	Sample Run Report with Template	A basic report for a sample run that can be customised using templates	Report for a Sample run	<div>Templates</div>
	VCF Report	A vcf report that can be customised using templates	Report for Variants in Sample	<div>Templates</div>
	Report with Template	A basic report that can be customised using templates	Report for Variants in Sample	<div>Templates</div>
	PDF Report with Template	A basic pdf report that can be customised using templates	Report for a Batch	<div>Templates</div>
Software	<div>+ Add Report</div>			

Figure 399: Initial view of the currently loaded reporting templates

New reports can be added by selecting the Add Report option and following the instructions in the popup window.




Figure 400: The Add Report popup window

Selecting the templates button for report shows the templates that are available to that particular report type




Figure 401: Templates associated with a report type

New templates can be added using the +Add button, selecting a template file and adding.




Figure 402: Add a template to a report

Guidelines

Guidelines provide a means to regulate classifications of variants by forcing the user to follow a question tree that will lead to a classification being generated. It is recommended to assess the template with reference to a user's own laboratory guidelines.

Using such a tree should remove any user variation with the result a classification will be generated consistently.

When first selected the guideline page is as follows:



Figure 403: The Guidelines start page

The guidelines page has two sections showing existing guidelines and the questions and transitions within them.




Figure 404: Existing guidelines and questions/transitions sections

Selecting ACMG guidelines will show a schematic of the questions and transitions. This is a complex decision tree as shown by the schematic below.

Selecting any of the questions, that are in the boxes, will show in the question in the box to the side as well as the rules associated with the question.





Figure 405: Schematic of the ACMG guidelines included with Interpret

Generating New Guidelines



Users can generate new guidelines to meet their own requirements. This is, currently, beyond the scope of this user manual. If this is required please contact OGT for further assistance.

Software

Advanced Settings

The advanced settings within Interpret are primarily for use by OGT support.

The default view shown in the figure; if any changes are required it is important to contact OGT in the first instance.

Advanced Settings

Unique Database Identifier

f0d899e2-7844-49ab-8fc2-b950b2d2373c

Processing Server Host Name

sureseq.local

Processing Server User Name

root

Processing Server User Password

.....

Processing Server File Storage Location

/Data/demo/filestore/

Previous Processing Server File Storage Locations

N/A

Processing Server Docker User

root

Web Application Server File Storage Location

/filestore

Previous Web Application Server File Storage Locations

N/A

Messaging Service Port

61618

Use Local FASTA and Cytoband Data for IGV

☒ Yes ☐ No

Install Default Panels

No

Install Default Tracks

No

Change default admin password

No

Processing Server Free Space Warning (bytes)

1,073,741,824

Temporary Storage Directory

/Data/demo/filestore/tmp

CytoSure Interpret Database Dialect

CytoSure Interpret Database Host

CytoSure Interpret Database Port

0

CytoSure Interpret Database Instance

CytoSure Interpret Database User

CytoSure Interpret Database Password

Interpret Base URL

Detect duplicate FASTQ uploads

☐ Yes ☒ No

Save Changes

Discard Changes

Logging

Root

WARN

Pipeline

INFO

Log Levels

Class	Log Level
velocity	WARN
com	
org	
java.sql.DatabaseMetaData	WARN
net.schmizz	
ro.fortsoft.pf4j	

Figure 406: Default advanced settings

ClinVar

Interpret is able to annotate detected variants with their pathogenicity classification according to ClinVar. ClinVar is updated weekly, and, in order to support up-to-date classification of variants, Interpret provides the option to manually update the ClinVar annotation used by the pipeline via **Admin Controls -> Software -> Advanced Settings**. To do this, simply click on the **Update** button next to the



appropriate **ClinVar Release Date** for the genome build in use. If the latest version of ClinVar is successfully downloaded and installed, the date associated with the ClinVar release will be displayed.

If a subsequently detected variant is found in the updated version of ClinVar, its pathogenicity will be displayed in the **ClinVar Significance** column in any variant tables.

Annotation

Annotation used within Interpret are managed on this page.
The initial table shows all the annotations loaded in the software.

Overview

User Controls

Analysis

Software

Software Overview

Annotation

Advanced Settings

Plug-ins

Admin Controls

Annotations:

Current Tracks:

Name	Last Updated	Track Type	Update
Genes	07-Jan-2020 12:32:45	Genes	
Aberration	03-Jan-2020 14:02:03	Aberration	
Gene	03-Jan-2020 14:57:59	Gene	
Exon	03-Jan-2020 14:02:48	Exon	
Copy Number Variation	03-Jan-2020 14:02:48	Copy Number Variation	
Segmental Duplication	03-Jan-2020 14:02:54	Segmental Duplication	
DECIPHER Patient	03-Jan-2020 14:02:54	DECIPHER Patient	
Syndrome	03-Jan-2020 15:00:44	Syndrome	
ClinGen Variation	03-Jan-2020 14:03:00	ClinGen Variation	
ClinGen Haploinsufficient Gene	03-Jan-2020 14:03:05	ClinGen Haploinsufficient Gene	
ClinGen Triplosensitive Gene	03-Jan-2020 14:03:06	ClinGen Triplosensitive Gene	
Cytocell FISH Probe	03-Jan-2020 14:03:11	Cytocell FISH Probe	

Create Track

Figure 407: The Annotation start page

Selecting **Create Track** provides a popup menu which allows a track to be created either from a file or by importing from an existing CytoSure Interpret installation.

Create Annotation Track:

Import From Cytosure Interpret

Upload File

Figure 408: Create track options

Plug-ins

The use of plug-ins is a mechanism whereby additional functionality can be provided.



In situations where a user requests some specific functionality it may be easier to generate a plug-in rather than a new version of the software.




Figure 409: The Plug-ins start page

Plug-ins can only be generated by OGT and if a plug-in provided then it can be loaded by selecting Add Plug-in and then using the file browser to select the file.




Figure 410: Load plug-in popup menu

Software Overview

The Software Overview provides version detail of the software installed.

Admin Controls

Overview

User Controls

Analysis

Software

Software Overview

Annotation

Advanced Settings

Plug-ins

Updates

Current SureSeq Version: 3.2.5

Contact Us

Name: *

Email: *

Reason: * Bug Report

Message: *

Send

Figure 411: The Software Overview start page

Reporting

Reporting has been discussed previously in the sections title Viewing Analysis Batches and Viewing Analysis Results where the details can be found.

Product-specific Guidance

Measurable Residual Disease

Overview

Detection and monitoring of Measurable Residual Disease (MRD) with the SureSeq Myeloid MRD Panel is made possible in Interpret through:

1. The ability to specify "Hotspots" (variants) which should be specifically interrogated by the pipeline for their presence at very low frequency.
2. The ability to visualise the change in allele frequency of these hotspots in multiple sequencing runs over time.

Discovery Mode

In order to identify candidate variants for use in "Monitoring Mode", where specific variants are interrogated by the pipeline in order to determine their allele frequency at very low depth, it may be necessary to process samples in "Discovery Mode".

Discovery Mode uses the standard SNV, Indel and ITD detection algorithms built into OGT's NGS analysis pipeline to report all variants present in a sample above a specific allele frequency and according to other quality-related criteria. To process samples in Discovery Mode:

1. Click on the **Batches** button in the toolbar and select **Run Batch**.
2. Enter a name for the batch in the **Batch Name** field.
3. Select the **SureSeq Myeloid MRD Panel** from the **Panel** drop-down list.
4. Select **Discovery Mode** from the **Protocol** drop-down list.
5. Select the samples to be processed from the list of available samples such that they are displayed in the **Selected Samples** table.
6. Click **Run Analysis**.
7. Click **OK**.

Once the batch has been started, the **Batch** page will be displayed showing the current status of the processing of the samples in the batch. The status of each sample can be updated by clicking the refresh button and, on completion, the **Completed Samples** table, displaying a summary of the results and relevant QC metrics, will appear.



Minimum Allele Frequency

By default, Discovery Mode is configured to detect variants at a minimum allele frequency of 1%. To reduce this value in order to increase the sensitivity, modify the Discovery Mode protocol as follows:

1. In the top-right corner of the screen, click on the user icon and select **Admin Controls**.
2. In the menu on the left-hand side, select **Analysis -> Protocols**.
3. In the **Protocols** list, select **Discovery Mode**.
4. Click the **Edit** button at the bottom of the screen.
5. Scroll down and select the **Advance Pipeline Configuration** tab.
6. In the **SNV Detection** section, modify the value of **Minimum Alt Fraction** as required.
7. Click **Save**.

Hotspot Monitoring

In order to visualise the results of hotspot monitoring:

1. Select **Tools -> Hotspot Monitoring Report**, and select the sample/source to be reported.




Figure 412: Selecting the sample(s) to be reported

- If necessary, enter the **Collection Date** of the sample(s). This only needs to be carried out once for each sample and will be remembered for future reports. Click **Next**.




Figure 413: Entering the collection date of the sample

- Select the hotspots to be reported using the same method described in [step 4e-f in the Selecting Hotspots](#) section and click **View**.




Figure 414: Selecting Hotspots to include in the report

- A graph containing the allele frequencies of all selected hotspots in all selected sample runs will be displayed, along with tabs allowing the user to view the results for individual hotspots. Graph images may be exported by right-clicking of the graph and selecting "Save Image". The table containing the data underlying the graph may be exported as a CSV via the **Export Data** button.




Figure 415: An example of a report



Monitoring Mode

To determine the allele frequency of hotspots in a batch of samples at very low depth, the samples should be processed using the "Monitoring Mode" protocol:

1. Upload the FASTQ files for the batch using the method described in the **Uploading FASTQ Files** section above.
2. Click on the **Batches** button in the toolbar and select **Run Batch**.
3. Enter a name for the batch in the **Batch Name** field.
4. Select the **SureSeq Myeloid MRD Panel** from the **Panel** drop-down list.
5. Select **Monitoring Mode** from the **Protocol** drop-down list.
6. Select the samples to be processed from the list of available samples such that they are displayed in the **Selected Samples** table.
7. Click **Run Analysis**.
8. Click **OK**.

Once the batch is complete, allele frequencies of hotspots selected in the Monitoring Mode protocol for a specific sample may be viewed by clicking on the **SNVs** button in the **Completed Samples**.

Selecting Hotspots

To select hotspots for use in Monitoring Mode, they must first be selected from variants identified in Discovery Mode:

1. In the **Batch** page for the Discovery Mode batch, click on the **SNVs** button in the **Completed Samples** table for a sample that may contain potential hotspots.
2. If necessary, filter the list of variants in the **Variants** page in order to identify potential hotspots more quickly.
3. For each variant to be monitored in Monitoring Mode:
 - a. Right-click on the variant in the table.
 - b. Select **Add to...**
 - i. If no Variant List has been created yet:
 1. Click **New List**
 2. Enter a **Name** for the list (e.g. "Hotspots")
 3. Click **Create**
 - ii. Otherwise, click on the name of the Variant List.
4. Once all required hotspots have been added to the Variant List, modify the "Monitoring Mode" protocol to use those hotspots:
 - a. In the top-right corner of the screen, click on the user icon and select **Admin Controls**.
 - b. In the menu on the left-hand side, select **Analysis -> Protocols**.
 - c. In the **Protocols** list, select **Monitoring Mode**.
 - d. Click the **Edit** button at the bottom of the screen.
 - e. Scroll down until the **Hotspots** table is displayed, and click on the name of the Variant List created in step 3 in the **List** column.

- f. In the **Variant** column, select all variants whose allele frequencies should be monitored, and click on the > button to add them to the **Selected Variants** table.
- g. When all variants have been added to the **Selected Variants** table, click the **Save** button.



Batch hotspot selection

The list of variants included in the Hotspots list for the Monitoring Mode protocol should cover all variants to be monitored in all samples in a batch. If different variants are relevant to different samples, (preferably) create the super-set of these variants in the protocol, or create separate protocols for each set of variants, and run the samples in different batches using the appropriate protocol.




Hotspots not detected in Discovery Mode

If the variant required for monitoring has not been detected in Discovery Mode, contact OGT for assistance to add the variant to a variant list.


Appendix

Attribute Definitions

Attribute	Variant Type	Description
#	SNV/Indel	Variant Database Identifier
# Alt Reads (-)	SNV/Indel	Number of alternative alleles on negative strand
# Alt Reads (+)	SNV/Indel	Number of alternative alleles on positive strand
# Markers	CNV/LOH	Number of bins (markers) used to identify CNVs
# Ref Reads (-)	SNV/Indel	Number of reference alleles on negative strand
# Ref Reads (+)	SNV/Indel	Number of reference alleles on positive strand
# Samples	SNV/Indel	Number of samples the variant is present in the database
% Samples	CNV/LOH	Percentage of samples the variant is present in the database
% Samples (Similar CNVs)	CNV/LOH	Number of samples in the database with an overlapping CNV
Allele Frequency	SNV/Indel	Allele Frequency
Alt	SNV/Indel	Alternate allele



Attribute	Variant Type	Description
Alt Depth	SNV/Indel	Number of reads supporting the alternative allele at the position
Alt Quality	SNV/Indel	Sum of alternative base qualities at the position
Alt Strand Bias	SNV/Indel	Sequencing bias in which one DNA strand is favoured over the other in the reads containing the alternative allele (Percentage)
Bands	CNV/LOH	Location of the variant on the chromosome
Batch Date	SNV/Indel	Date the batch was performed
Batch Name	SNV/Indel	Name of the batch containing the run
Canonical?	SNV/Indel	A flag indicating if the transcript is denoted as the canonical transcript for this gene
Chromosome	CNV/LOH	Chromosome of the CNV/LOH
Chromosome	SNV/Indel	Chromosome of the variant
Classification	CNV/LOH	User-assigned classification of the variant
ClinVar Significance	SNV/Indel	Clinical significance of variant according to ClinVar (e.g. benign, pathogenic, uncertain significance etc.)
Confidence	CNV/LOH	Confidence that the call is correct (e.g. High, Low) dependant on the Standard Error of Mean
Consequence Terms	SNV/Indel	Most severe outcome caused by the specific variant (e.g Frameshift variant, Stop gained, Synonymous variant etc.)
Context Length	SNV/Indel	Length of the genomic context overlapping the variant
Copy Number	CNV/LOH	Number of copies of the CNV event
Depth	CNV/LOH	Depth of coverage of the sample at position of the CNV
Description	Translocation	Donor and recipient gene symbol pair
Donor Breakpoint	Translocation	Position on the donor chromosome of the translocation
Donor Chromosome	Translocation	The chromosome number of the donor gene
Donor Gene	Translocation	Donor gene symbol where the translocation originated




Attribute	Variant Type	Description
Donor Locus Reads	Translocation	Read depth at the donor breakpoint position
Donor Orientation	Translocation	Orientation (strand) of the donor gene
Donor Reads Position	Translocation	Which side of the breakpoint the donor read lies (left or right)
End	CNV/LOH	Genomic position of end of CNV
End	SNV/Indel	Genomic position of end of variant
Estimated Tumour Content	CNV/LOH	Estimated tumour content (only used in cancer panels)
Exon ID	SNV/Indel	Unique ID for the exon
Exon Number	SNV/Indel	The number of the exon in the gene the variant is present
Frequency	CNV/LOH	Average Allele Frequency of common SNP's overlapping the CNV
Fusion Type	Translocation	The method of detection that highlighted the fusion, either: Over-expression of the gene (Expression), detection of a known fusion by sufficient supporting reads (Canonical), detection of an unknown fusion by sufficient supporting reads (Non-canonical)
Gene ID	SNV/Indel	Unique ID of the gene where the variant is
Gene Symbol	SNV/Indel	Gene symbol where the variant is
Genes	CNV/LOH	List of genes overlapping the CNV
Genome Build	CNV/LOH	Genome assembly version
Genomic Context	SNV/Indel	Genomic context the variant is overlapping (Low Complexity, Homopolymer, Simple Repeat)
Genotype	SNV/Indel	Genotype (Heterozygous/Homozygous)
gnomAD - African	SNV/Indel	The frequency variant appears on the Genome Aggregation Database from people of African decent
gnomAD - Ashkenazi Jewish	SNV/Indel	The frequency variant appears on the Genome Aggregation Database from people of Ashkenazi Jewish decent


Attribute	Variant Type	Description
gnomAD - East Asian	SNV/Indel	The frequency the variant appears on the Genome Aggregation Database from people of East Asian decent
gnomAD - European (Finnish)	SNV/Indel	The frequency the variant appears on the Genome Aggregation Database from people of Finnish decent
gnomAD - European (non-Finnish)	SNV/Indel	The frequency the variant appears on the Genome Aggregation Database from people of European (Non-Finnish) decent
gnomAD - Latino	SNV/Indel	The frequency variant appears on the Genome Aggregation Database from people of Latino decent
gnomAD - Other	SNV/Indel	The frequency the variant appears on the Genome Aggregation Database from people of another decent
gnomAD - South Asian	SNV/Indel	The frequency the variant appears on the Genome Aggregation Database from people of South Asian decent
gnomAD - Total	SNV/Indel	The frequency variant appears on the Genome Aggregation Database from all reference genomes
HGVSc	SNV/Indel	The HGVS coding sequence name
HGVSc (Gene Symbol)	SNV/Indel	The HGVS coding sequence name with the Transcript identifier replaced with its Gene Symbol
HGVSp	SNV/Indel	The HGVS protein sequence name
Homozygosity	SNV/Indel	Proportion of the genome covered by LOH regions larger than 5Mb
Impact	SNV/Indel	The Impact score according to Ensembl VEP of the genetic variation in the genetic sequence (e.g. LOW, MODERATE, HIGH etc.)
Inheritance	CNV/LOH	Estimated inheritance of the variant based on the presence of the variant in parental results, if available.
Inheritance	SNV/Indel	Estimated inheritance of the variant based on the presence of the variant in parental results, if available.
Inheritance	Translocation	Estimated inheritance of the variant based on the presence of the variant in parental results, if available.

Attribute	Variant Type	Description
ISCN	CNV/LOH	CNV/LOH variant encoded according to ISCN (International System for Human Cytogenomic Nomenclature)
ISCN	Translocation	Translocation variant encoded according to ISCN (International System for Human Cytogenomic Nomenclature)
Length	CNV/LOH	Length of CNV
Log Ratio	CNV/LOH	Mean log2 ratio of sample/reference of the CNV
Mean	CNV/LOH	Rescaled mean log2 of sample/reference of the CNV (only used in cancer panels)
Mean Standard Error	CNV/LOH	Standard Error of the Mean
Minor Allele	SNV/Indel	Base of the minor allele
Minor Allele Frequency	SNV/Indel	Rate at which the second most common allele occurs
Mosaicism	CNV/LOH	Estimate of the percentage of mosaicism observed in CNV region
Mosaicism Lower Bound	CNV/LOH	Estimate of the lower bound of mosaicism observed in the CNV region
Mosaicism Range	CNV/LOH	Estimate of the range of mosaicism observed in the sample
Mosaicism Upper Bound	CNV/LOH	Estimate of the upper bound of mosaicism observed in the CNV region
Most Severe Consequence	SNV/Indel	Most severe outcome caused by the specific variant (e.g Frameshift variant, Stop gained, Synonymous variant etc.)
Normalised Expression	Translocation	The expression of the baited gene relative to the housekeeping genes and normalised by total read count
Overlap	CNV/LOH	Genomic context of the CNV
P Value	Translocation	Probability of observing the translocation
Panel	SNV/Indel	Panel used for the analysis
PolyPhen Prediction	SNV/Indel	The prediction of how damaging a variant will be, based off the PolyPhen Score

Attribute	Variant Type	Description
PolyPhen Score	SNV/Indel	The probability that a substitution is damaging (e.g. 0.25 benign, 0.5 possibly damaging, 0.95 probably damaging)
Proportion	Translocation	Proportion of split reads over total reads at the donor breakpoint
Protein ID	SNV/Indel	Unique ID for the protein
Protocol	SNV/Indel	OGT Interpret software protocol used to analyse the run
Quality	CNV/LOH	(Not implemented)
Quality Score	SNV/Indel	Phred Quality score of the variant
Ratio	SNV/Indel	Ratio of depth observed in duplicated PTD exons compared to the exons in the rest of the gene
Read 1	SNV/Indel	File name of the FASTQ from R1 reads
Read 1 Size	SNV/Indel	Size of the FASTQ file from R1 reads
Read 2	SNV/Indel	File name of the FASTQ from R2 reads
Read 2 Size	SNV/Indel	Size of the FASTQ file from R2 reads
Reads Placed Left	SNV/Indel	Number of reads with supporting evidence to the left of the variant
Reads Placed Right	SNV/Indel	Number of reads with supporting evidence to the right of the variant
Recipient Breakpoint	Translocation	Position on the recipient chromosome of the translocation
Recipient Chromosome	Translocation	The chromosome of the recipient gene
Recipient Gene	Translocation	Recipient gene symbol where the translocation ended up
Recipient Locus Reads	Translocation	Read depth at the recipient breakpoint position
Recipient Orientation	Translocation	Orientation (strand) of the recipient gene
Recipient Reads Position	Translocation	Which side of the variant the donor read lies (left or right)
Ref	SNV/Indel	Reference nucleotide base
Ref Depth	SNV/Indel	Number of reads supporting the alternative allele at the position



Attribute	Variant Type	Description
Ref Quality	SNV/Indel	Sum of alternative reference qualities at the position
Ref Strand Bias	SNV/Indel	Sequencing bias in which one DNA strand is favoured over the other in the reads containing the reference allele (Percentage)
rsID	SNV/Indel	SNP id from NCBI dbSNP
Sample	CNV/LOH	ID of the sample containing the CNV
Sample	SNV/Indel	ID of the sample containing this variant
Sample	Translocation	ID of the sample containing this variant
Sample ID	SNV/Indel	ID of the sample containing this variant
Score	CNV/LOH	LOH score (Higher scores >30 indicate a higher confidence in the call)
Sex	SNV/Indel	Inferred chromosomal sex of the sample (Male, Female, Unknown)
SIFT Prediction	SNV/Indel	Prediction of how detrimental a variant will be to protein function (The opposite of Polyphen in terms of numbering)
SIFT Score	SNV/Indel	A score that predicts whether a variant will affect protein function (0 = deleterious , 1 = tolerated)
Source	CNV/LOH	Tool used for CNV identification
Start	CNV/LOH	Genomic position of start of CNV
Start	SNV/Indel	Genomic start position of variant
Supporting Reads	Translocation	The sum of split and discordant reads in support of the fusion call
Total Depth	SNV/Indel	Depth of coverage at the position
Transcript ID	SNV/Indel	Unique ID of the specific selected transcript
Transcript Resolution Method	SNV/Indel	Method used to determine which transcript to use
Type	CNV/LOH	Variant type (e.g. CNV, LOH)
Type	SNV/Indel	Variant type (e.g. SNV, ITD, PTD, etc.)
Type	Translocation	Variant type (e.g. Translocation)
User	SNV/Indel	Login name of user which ran the batch



Attribute	Variant Type	Description
VEP Version	SNV/Indel	Version of Ensembl Variant Effect Predictor
Zygosity	SNV/Indel	The degree at which both copies of the chromosome have the same genetic sequence (e.g. Homozygous or Heterozygous)

Table 9: Definitions of attributes displayed in various tables in Interpret

Legal Information

This handbook and its contents are © Oxford Gene Technology (Operations) Limited 2025. All rights reserved. Reproduction of all or any substantial part of its contents in any form is prohibited except that individual users may print or save portions of the protocol for their own personal use. This licence does not permit users to incorporate the material or any substantial part of it in any other work or publication, whether in hard copy or electronic or any other form. In particular (but without limitation), no substantial part of the handbook may be distributed or copied for any commercial purpose.

Interpret Software

Interpret software was developed by Oxford Gene Technology. The purchaser has the non-transferable right to use the product for RESEARCH USE ONLY AND NOT FOR DIAGNOSTICS PROCEDURES. It is not intended for use, and should not be used, for the diagnosis, prevention, monitoring, treatment or alleviation of any disease or condition, or for the investigation of any physiological process, in any identifiable human, or for any other medical purpose.

Customer's obligations

The Customer acknowledges that Oxford Gene Technology (Operations) Limited (or its group companies) owns all intellectual property rights in the design of the Product, including the choice and configuration of the oligonucleotide sequences used in the Product. The Product may only be reproduced or manufactured by Oxford Gene Technology (Operations) Limited or with its permission.



Contact information

UK +44 (0) 1865 856826

US +1 914 467 5285

Technical support: support@ogt.com

contact@ogt.com

ogt.com

Oxford Gene Technology Ltd.

Unit 5, Oxford Technology Park, 4A Technology Drive, Kidlington, Oxfordshire, OX5 1GN, UK



**What binds us,
makes us.**

SureSeq™ and CytoSure®: For Research Use Only; Not for Diagnostic Procedures. This document and its contents are © Oxford Gene Technology IP Limited – 2024. All rights reserved. OGT™, CytoSure® and SureSeq™ are trademarks of Oxford Gene Technology IP Limited.

Oxford Gene Technology (Operations) Ltd. Registered in England No: 03845432 Unit 5, Oxford Technology Park, 4A Technology Drive, Kidlington, Oxfordshire, OX5 1GN, UK